

COPY

EXHIBIT G-2  
THIS REPORT PREPARED PURSUANT TO DUTY IMPOSED BY A.A.C. R13-10-104 (A)

ARIZONA DEPARTMENT OF PUBLIC SAFETY

STANDARD QUALITY ASSURANCE PROCEDURES  
INTOXILYZER MODEL 8000

STANDARD CALIBRATION CHECK PROCEDURE

QA SPECIALIST J. Mural AGENCY Yavapai County Sheriff's Office  
DATE September 4, 2013 TIME 08:49:53  
INTOXILYZER SERIAL # 80-01734 LOCATION Seligman 103

- (✓) 1. Ensure that gas tank is attached to instrument and contains a standard alcohol concentration solution 100 AC.  
OR  
Pour a standard alcohol concentration solution ✓ AC, into a clean dry simulator and assemble the simulator. Ensure that a tight seal has been made. Turn on the simulator and allow temperature to reach  $34^{\circ}\text{C} \pm 0.2^{\circ}\text{C}$
- (✓) 2. Intoxilyzer 8000 display reads "PUSH BUTTON TO START"
- (✓) 3. Go to the "Control Testing Menu". Select "D" for dry control test or "W" for wet control test. After selection is made press ENTER.
- (✓) 4. Air blank completed.
- (✓) 5. Calibration check completed. Test results 0.096 AC.
- (✓) 6. Air blank completed.
- (✓) 7. Remove printed record. Attach the record to the completed checklist.

SIGNATURE J. Mural NSB/5597



INTOXILYZER 8000

Location: YCSO SELIGMAN

Serial Number: 80-001734

Core Version: 8105.48

09/01/2013

08:49:53

-----  
Standard Lot#: 09312100A1

Last Changed By: J. MURRAY #N56  
-----

QAS: J. MURRAY #N56

YAVAPAI YCSO  
-----

Operator: J. MURRAY #N56

YAVAPAI YCSO  
-----

Test	g/210L	Time
Air Blank	0.000	08:51:14
0.100 Cal Check	0.096	08:51:35
Air Blank	0.000	08:52:05

Preventative Maintenance  
Performed -- Timer Reset



EXHIBIT G-4  
THIS REPORT PREPARED PURSUANT TO DUTY IMPOSED BY A.A.C. R13-10-104 (A)

ARIZONA DEPARTMENT OF PUBLIC SAFETY

STANDARD QUALITY ASSURANCE PROCEDURES  
INTOXILYZER MODEL 8000

STANDARD QUALITY ASSURANCE PROCEDURE

QA SPECIALIST J. Murray NSB/5597 AGENCY Yavapai County Sheriff's Office  
DATE September 1, 2013 TIME 09:10:19  
INTOXILYZER SERIAL # 80-01734 LOCATION Seligman 103

1. Display Reads "PUSH BUTTON TO START"

DIAGNOSTIC TESTS

- 1. Clock time check.
- 2. Date check.

OPERATIONAL TESTS


- 1. Alcohol-free subject test result 0. 060 AC.
- 2. Error recognition logic system functioning.  
Not a Successfully Completed Test Sequence printed
- 3. Proper sample recognition system.  
Not a Successfully Completed Test Sequence printed  
Deficient sample printed.
- 4. Standard Calibration Check standard 0. 160 AC. Result 0. 096 AC.

Instrument is operating properly and accurately. YES  NO

COMMENTS

6/1/2014 800PST  
09312900A4

SIGNATURE J. Murray NSB/5597



INTOXILYZER 8000  
Location: YCSO SELIGMAN  
Serial Number: 80-001734  
Core Version: 8105.48  
09/01/2013

09:10:19

---

Standard Lot#: 09312100A1  
Last Changed By: J. MURRAY #N56

---

QAS: J. MURRAY #N56  
YAVAPAI YCSO

---

Operator: J. MURRAY #N56  
YAVAPAI YCSO

---

Subject: SQAP,TEST  
DOB: 10/20/1960  
Sex: F Weight: 115

---

15 Minute Deprivation Period? Yes

---

Test	g/210L	Time
Air Blank	0.000	09:11:16
Diagnostic Test	Pass	09:11:51
Air Blank	0.000	09:12:18
0.100 Cal Check	0.096	09:12:39
Air Blank	0.000	09:13:08
Reference	RFI*	09:13:13
Air Blank	0.000	09:13:41

\*RFI Detect

**Not a Successfully Completed Test Sequence**



INTOXILYZER 8000

Location: YCSO SELIGMAN

Serial Number: 80-001734

Core Version: 8105.48

09/01/2013

09:13:44

---

Standard Lot#: 09312100A1

Last Changed By: J. MURRAY

#N56

---

QAS: J. MURRAY

#N56

YAVAPAI YCSO

---

Operator: J. MURRAY

#N56

YAVAPAI YCSO

---

Subject: SQAP,TEST

DOB: 10/20/1960

Sex: F

Weight: 115

---

15 Minute Deprivation Period? Yes

---

Test	g/210L	Time
Air Blank	0.000	09:14:46
Diagnostic Test	Pass	09:15:21
Air Blank	0.000	09:15:49
0.100 Cal Check	0.097	09:16:09
Air Blank	0.000	09:16:39
Subject Test	0.000	09:17:08
Air Blank	0.000	09:17:36
Wait	ABT*	
Air Blank	RFI**	09:17:52

\*Sequence Aborted

\*\*RFI Detect

**Not a Successfully Completed Test Sequence**



INTOXILYZER 8000  
Location: YCSO SELIGMAN  
Serial Number: 80-001734  
Core Version: 8105.48  
09/01/2013

09:17:59

---

Standard Lot#: 09312100A1  
Last Changed By: J. MURRAY #N56

---

QAS: J. MURRAY #N56  
YAVAPAI YCSO

---

Operator: J. MURRAY #N56  
YAVAPAI YCSO

---

Subject: SQAP, TWST  
DOB: 10/20/1960  
Sex: F Weight: 115

---


15 Minute Deprivation Period? Yes

---

Test	g/210L	Time
Air Blank	IPS*	09:18:38
Air Blank	RFI**	09:18:46

\*Improper Sample  
\*\*RFI Detect

**Not a Successfully Completed Test Sequence**



INTOXILYZER 8000  
Location: YCSO SELIGMAN  
Serial Number: 80-001734  
Core Version: 8105.48  
09/01/2013

09:20:11

---

Standard Lot#: 09312100A1  
Last Changed By: J. MURRAY #N56

---

QAS: J. MURRAY #N56  
YAVAPAI YCSO

---

Operator: J. MURRAY #N56  
YAVAPAI YCSO

---

Subject: SQAP,TEST  
DOB: 10/20/1960  
Sex: F Weight: 115

---

15 Minute Deprivation Period? Yes

---

Test	g/210L	Time
Air Blank	0.000	09:21:12
Diagnostic Test	Pass	09:21:47
Air Blank	0.000	09:22:15
0.100 Cal Check	0.096	09:22:35
Air Blank	0.000	09:23:05
Subject Test	0.217*	09:26:12
Air Blank	0.000	09:26:46
Air Blank	0.000	09:27:13
Subject Test	XXX**	09:28:13
Air Blank	0.000	09:28:48

\*Deficient Sample

\*\*Invalid Sample

Begin new deprivation period

**Not a Successfully Completed Test Sequence**



EXHIBIT G-2  
THIS REPORT PREPARED PURSUANT TO DUTY IMPOSED BY A.A.C. R13-10-104 (A)

ARIZONA DEPARTMENT OF PUBLIC SAFETY  
STANDARD QUALITY ASSURANCE PROCEDURES  
INTOXILYZER MODEL 8000

STANDARD CALIBRATION CHECK PROCEDURE

QA SPECIALIST J. Murray NS6/5597 AGENCY Yavapai County Sheriff's Office  
DATE September 1, 2013 TIME 09:36:44  
INTOXILYZER SERIAL # 80-01734 LOCATION Selyman 103

- (x) 1. Ensure that gas tank is attached to instrument and contains a standard alcohol concentration solution .100 AC.  
OR  
Pour a standard alcohol concentration solution            AC, into a clean dry simulator and assemble the simulator. Ensure that a tight seal has been made. Turn on the simulator and allow temperature to reach  $34^{\circ}\text{C} \pm 0.2^{\circ}\text{C}$
- (x) 2. Intoxilyzer 8000 display reads "PUSH BUTTON TO START"
- (x) 3. Go to the "Control Testing Menu". Select "D" for dry control test or "W" for wet control test. After selection is made press ENTER.
- (x) 4. Air blank completed.
- (x) 5. Calibration check completed. Test results 0. 096 AC.
- (x) 6. Air blank completed.
- (x) 7. Remove printed record. Attach the record to the completed checklist.

SIGNATURE J. Murray NS6/5597



INTOXILYZER 8000

Location: YCSO SELIGMAN

Serial Number: 80-001734

Core Version: 8105.48

09/01/2013

09:36:44

-----  
Standard Lot#: 09312100A1

Last Changed By: J. MURRAY #N56  
-----

QAS: J. MURRAY #N56

YAVAPAI YCSO  
-----

Operator: J. MURRAY #N56

YAVAPAI YCSO  
-----

Test	g/210L	Time
Air Blank	0.000	09:37:26
0.100 Cal Check	0.096	09:37:47
Air Blank	0.000	09:38:16

Preventative Maintenance  
Performed -- Timer Reset

