

EXHIBIT G-2
THIS REPORT PREPARED PURSUANT TO DUTY IMPOSED BY A.A.C. R13-10-104 (A)

ARIZONA DEPARTMENT OF PUBLIC SAFETY
STANDARD QUALITY ASSURANCE PROCEDURES
INTOXILYZER MODEL 8000

STANDARD CALIBRATION CHECK PROCEDURE

QA SPECIALIST J. MURRAY NS6/5597 AGENCY Yavapai County Sheriff's Office
DATE 11/19/2014 TIME 14:31:10
INTOXILYZER SERIAL # 80-01734 LOCATION Seligman 103

- (X) 1. Ensure that gas tank is attached to instrument and contains a standard alcohol concentration solution, 100 AC.
OR
Pour a standard alcohol concentration solution AC, into a clean dry simulator and assemble the simulator. Ensure that a tight seal has been made. Turn on the simulator and allow temperature to reach 34°C ± 0.2°C
- (X) 2. Intoxilyzer 8000 display reads "PUSH BUTTON TO START"
- (X) 3. Go to the "Control Testing Menu". Select "D" for dry control test or "W" for wet control test. After selection is made press ENTER.
- (X) 4. Air blank completed.
- (X) 5. Calibration check completed. Test results 0. 099 AC.
- (X) 6. Air blank completed.
- (X) 7. Remove printed record. Attach the record to the completed checklist.

SIGNATURE J. Murray NS6/5597

DPS Form Exh G-2 (Rev 05-01)

QAS/SDAP Performed @ Prescott 103

11/19/2014

NS6

INTOXILYZER 8000
Location: YCSO SELIGMAN
Serial Number: 80-001734
Core Version: 8105.48
01/19/2014 14:31:10

Standard Lot#: 09312100A1
Last Changed By: J. MURRAY #N56

QAS: J. MURRAY5597 #N56
YAVAPAI YCSO

Operator: J. MURRAY5597 #N56
YAVAPAI YCSO

Test	g/210L	Time
Air Blank	0.000	14:31:52
0.100 Cal Check	0.099	14:32:13
Air Blank	0.000	14:32:42

Preventative Maintenance
Performed -- Timer Reset

EXHIBIT G-4
THIS REPORT PREPARED PURSUANT TO DUTY IMPOSED BY A.A.C. R13-10-104 (A)

ARIZONA DEPARTMENT OF PUBLIC SAFETY
STANDARD QUALITY ASSURANCE PROCEDURES
INTOXILYZER MODEL 8000

STANDARD QUALITY ASSURANCE PROCEDURE

QA SPECIALIST J. MULLAY NS6/5597 AGENCY Yavapai County Sheriff's Office
DATE 1/19/2014 TIME 14:35:08
INTOXILYZER SERIAL # 80-01234 LOCATION Seligman 103

1. Display Reads "PUSH BUTTON TO START"

DIAGNOSTIC TESTS

1. Clock time check.
 2. Date check.

OPERATIONAL TESTS

1. Alcohol-free subject test result 0. 000 AC.
 2. Error recognition logic system functioning.
Not a Successfully Completed Test Sequence printed
 3. Proper sample recognition system.
Not a Successfully Completed Test Sequence printed
Deficient sample printed.
 4. Standard Calibration Check standard 0. 100 AC. Result 0. 099 AC.

Instrument is operating properly and accurately. YES NO

COMMENTS

OAS/SDAP Performed @ Pursatt 103 1/19/2014

SIGNATURE J. MULLAY NS6/5597

INTOXILYZER 8000
Location: YCSO SELIGMAN
Serial Number: 80-001734
Core Version: 8105.48
01/19/2014 14:35:08

Standard Lot#: 09312100A1
Last Changed By: J. MURRAY #N56

QAS: J. MURRAY5597 #N56
YAVAPAI YCSO

Operator: J. MURRAY5597 #N56
YAVAPAI YCSO

Subject: SQAP,TEST
DOB: 10/20/1960
Sex: F Weight: 115

15 Minute Deprivation Period? Yes

Test	g/210L	Time
Air Blank	0.000	14:36:11
Diagnostic Test	Pass	14:36:46
Air Blank	0.000	14:37:13
0.100 Cal Check	0.099	14:37:34
Air Blank	0.000	14:38:03
Subject Test	0.000	14:38:27
Air Blank	0.000	14:38:55
Wait	ABT*	
Air Blank	IPS**	14:39:11

*Sequence Aborted

**Improper Sample

Not a Successfully Completed Test Sequence

INTOXILYZER 8000
Location: YCSO SELIGMAN
Serial Number: 80-001734
Core Version: 8105.48
01/19/2014

14:39:17

Standard Lot#: 09312100A1
Last Changed By: J. MURRAY #N56

QAS: J. MURRAY5597 #N56
YAVAPAI YCSO

Operator: J. MURRAY5597 #N56
YAVAPAI YCSO

Subject: SQAP,TEST
DOB: 10/20/1960
Sex: F Weight: 115

15 Minute Deprivation Period? Yes

Test	g/210L	Time
Air Blank	RFI*	14:39:55
Air Blank	IPS**	14:40:04

*RFI Detect
**Improper Sample

Not a Successfully Completed Test Sequence

INTOXILYZER 8000
Location: YCSO SELIGMAN
Serial Number: 80-001734
Core Version: 8105.48
01/19/2014

14:40:41

Standard Lot#: 09312100A1
Last Changed By: J. MURRAY #N56

QAS: J. MURRAY5597 #N56
YAVAPAI YCSO

Operator: J. MURRAY5597 #N56
YAVAPAI YCSO

Subject: SQAP,TEST
DOB: 10/20/1960
Sex: F Weight: 115

15 Minute Deprivation Period? Yes

Test	g/210L	Time
Air Blank	0.000	14:41:40
Diagnostic Test	Pass	14:42:15
Air Blank	0.000	14:42:42
0.100 Cal Check	0.099	14:43:03
Air Blank	0.000	14:43:33
Subject Test	0.071*	14:46:40
Air Blank	0.000	14:47:12
Air Blank	0.000	14:47:39
Subject Test	XXX**	14:48:16
Air Blank	0.000	14:48:50

*Deficient Sample

**Invalid Sample

Begin new deprivation period

Not a Successfully Completed Test Sequence

EXHIBIT G-2
THIS REPORT PREPARED PURSUANT TO DUTY IMPOSED BY A.A.C. R13-10-104 (A)

ARIZONA DEPARTMENT OF PUBLIC SAFETY
STANDARD QUALITY ASSURANCE PROCEDURES
INTOXILYZER MODEL 8000

STANDARD CALIBRATION CHECK PROCEDURE

QA SPECIALIST J. M. MUSTAY NS6/5597 AGENCY Yavapai County Sheriff's Office
DATE 1/19/2014 TIME 14:49:09
INTOXILYZER SERIAL # 80-01734 LOCATION Selma 103

- (X) 1. Ensure that gas tank is attached to instrument and contains a standard alcohol concentration solution 100 AC.
OR
Pour a standard alcohol concentration solution AC, into a clean dry simulator and assemble the simulator. Ensure that a tight seal has been made. Turn on the simulator and allow temperature to reach $34^{\circ}\text{C} \pm 0.2^{\circ}\text{C}$
- (+) 2. Intoxilyzer 8000 display reads "PUSH BUTTON TO START"
(+) 3. Go to the "Control Testing Menu". Select "D" for dry control test or "W" for wet control test. After selection is made press ENTER.
(+) 4. Air blank completed.
(+) 5. Calibration check completed. Test results 0.099 AC.
(+) 6. Air blank completed.
(X) 7. Remove printed record. Attach the record to the completed checklist.

SIGNATURE J. M. MUSTAY NS6/5597

DPS Form Exh G-2 (Rev 05-01)

QAS/SDAP performed @ Alcott 103
1/19/2014
(NS6)

INTOXILYZER 8000

Location: YCSO SELIGMAN

Serial Number: 80-001734

Core Version: 8105.48

01/19/2014

14:49:09

Standard Lot#: 09312100A1

Last Changed By: J. MURRAY #N56

QAS: J. MURRAY5597 #N56

YAVAPAI YCSO

Operator: J. MURRAY5597 #N56

YAVAPAI YCSO

Test	g/210L	Time
Air Blank	0.000	14:49:48
0.100 Cal Check	0.099	14:50:08
Air Blank	0.000	14:50:38

Preventative Maintenance
Performed -- Timer Reset