

COPY

EXHIBIT G-2
THIS REPORT PREPARED PURSUANT TO DUTY IMPOSED BY A.A.C. R13-10-104 (A)

ARIZONA DEPARTMENT OF PUBLIC SAFETY
STANDARD QUALITY ASSURANCE PROCEDURES
INTOXILYZER MODEL 8000
STANDARD CALIBRATION CHECK PROCEDURE

QUALITY SPECIALIST J. MURPHY AGENCY Yavapai County Sheriff's Office
DATE 6/28/15 TIME 14:23:16
INTOXILYZER SERIAL # 80-01734 LOCATION Selma

- 1. Ensure that gas tank is attached to instrument and contains a standard alcohol concentration solution 1.00 AC.
OR
Pour a standard alcohol concentration solution 1.00 AC into a clean dry simulator and assemble the simulator. Ensure that a tight seal has been made. Turn on the simulator and allow temperature to reach $34\text{ }^{\circ}\text{C} \pm 0.2\text{ }^{\circ}\text{C}$.
- 2. Intoxilyzer 8000 display reads "PUSH BUTTON TO START"
- 3. Go to the "Control Testing Menu". Select "D" for dry control test or "W" for wet control test. After selection is made press ENTER.
- 4. Air blank completed.
- 5. Calibration check completed. Test results 0.098 AC.
- 6. Air blank completed.
- 7. Remove printed record. Attach the record to the completed checklist.

SIGNATURE J. MURPHY NSB/SS97

EXHIBIT G-4
THIS REPORT PREPARED PURSUANT TO DUTY IMPOSED BY A.A.C. R13-10-104 (A)

ARIZONA DEPARTMENT OF PUBLIC SAFETY

STANDARD QUALITY ASSURANCE PROCEDURES
INTOXILYZER MODEL 8000

STANDARD QUALITY ASSURANCE PROCEDURE

NAME: J. MULLY

DATE: 6/28/15

TEST NUMBER: 80-01234

AGENCY: Yavapai County Sheriff's Office

TIME: 14:30:17

LOCATION: Selma

Intoxilyzer 8000 power ON TO START

Intoxilyzer 8000
Calibration Check

Intoxilyzer 8000 Model 0000 AC

Intoxilyzer 8000 recognition logic system functioning.

Intoxilyzer 8000 test sequence completed. Test sequence printed.

Intoxilyzer 8000 sample recognition system.

Intoxilyzer 8000 Not a Successfully completed. Test sequence printed.

Intoxilyzer 8000 Calibration Check standard 0. 100 AC. Result 0. 09.7 AC.

Intoxilyzer 8000 operating properly and accurately. YES NO

SN# 33913700A3 exp 1/1/2016 800 ASI

NAME: J. MULLY NSB/5597

TEST NUMBER: 80-01234

INTOXILYZER 8000
LOCATION: YCSO SELIGMAN
Serial Number: 80-001734
Cann Version: 8105.48
08/25/2015

14:30:17

Standard Lot#: 33913100A3
Last Changed By: J. MURRAY5597 #N56

OAS: J. MURRAY5597 #N56
YAVAPAI SHERIFFS

Operator: J. MURRAY5597 #N56
YAVAPAI SHERIFFS

Subject: SOAP, BEST
DOB: 10/20/1960

Sex: F Weight: 115

15 Minute Deprivation Period? Yes

Test	g/210L	Time
Air Blank	0.000	14:31:19
Magnetic Test	Pass	14:31:51
Air Blank	0.000	14:32:22
0.100 Air Check	0.097	14:32:42
Air Blank	0.000	14:33:12
Subject Test	0.000	14:33:36
Air Blank	0.000	14:34:04
Wait	ABT*	
Air Blank	IPS**	14:34:18

**Sequence Aborted
**Improper Sample

Not a Successfully Completed Test Sequence

INTOXILYZER 8000

Location: YCSO SELIGMAN

Serial Number: 80-001734

Core Version: 8105.48

06/28/2015

14:34:28

Standard Lot#: 33913100A3

Last Changed By: J. MURRAY5597 #N56

QAS: J. MURRAY5597 #N56

YAVAPAI SHERIFFS ..

Operator: J. MURRAY5597 #N56

YAVAPAI SHERIFFS

Subject: SQAP,TEST

DOB: 10/20/1960

Sex: F

Weight: 115

15 Minute Deprivation Period? Yes

Test	g/210L	Time
Air Blank	0.000	14:35:26
Diagnostic Test	Pass	14:36:01
Air Blank	RFI*	14:36:12
Air Blank	0.000	14:36:40

*RFI Detect

Not a Successfully Completed Test Sequence

INTOXILYZER 8000
Location: YCSO SELIGMAN
Serial Number: 80-001734
Core Version: 8105.48
06/28/2015

14:36:45

Standard Lot#: 33913100A3
Last Changed By: J. MURRAY5597 #N56

QAS: J. MURRAY5597 #N56
YAVAPAI SHERIFFS

Operator: J. MURRAY5597 #N56
YAVAPAI SHERIFFS

Subject: SQAP, TEST
DOB: 10/20/1960
Sex: F Weight: 115

15 Minute Deprivation Period? Yes

Test	g/210L	Time
Air Blank	0.000	14:37:44
Diagnostic Test	Pass	14:38:19
Air Blank	0.000	14:38:46
0.100 Cal Check	0.097	14:39:07
Air Blank	0.000	14:39:36
Subject Test	0.174*	14:42:44
Air Blank	0.000	14:43:16
Air Blank	0.000	14:43:43
Subject Test	XXX**	14:44:26
Air Blank	0.000	14:45:00

*Deficient Sample

**Invalid Sample

Begin new deprivation period

Not a Successfully Completed Test Sequence

EXHIBIT G-2
THIS REPORT PREPARED PURSUANT TO DUTY IMPOSED BY A.A.C. R13-10-104 (A)

ARIZONA DEPARTMENT OF PUBLIC SAFETY
STANDARD QUALITY ASSURANCE PROCEDURES
INTOXILYZER MODEL 8000

STANDARD CALIBRATION CHECK PROCEDURE

QA SPECIALIST J. Murray AGENCY Yavapai County Sheriff's Office
DATE: 6/28/15 TIME 14:45:14
INTOXILYZER SERIAL # 80-01734 LOCATION Seligman

- (X) 1. Ensure that gas tank is attached to instrument and contains a standard alcohol concentration solution 100 AC.
OR
Pour a standard alcohol concentration solution AC, into a clean dry simulator and assemble the simulator. Ensure that a tight seal has been made. Turn on the simulator and allow temperature to reach $34^{\circ}\text{C} \pm 0.2^{\circ}\text{C}$
- (X) 2. Intoxilyzer 8000 display reads "PUSH BUTTON TO START"
- (X) 3. Go to the "Control Testing Menu". Select "D" for dry control test or "W" for wet control test. After selection is made press ENTER.
- (X) 4. Air blank completed.
- (X) 5. Calibration check completed. Test results 0. 097 AC.
- (X) 6. Air blank completed.
- (X) 7. Remove printed record. Attach the record to the completed checklist.

SIGNATURE J. Murray N52/5597

INTOXILYZER 8000

Location: YCNO SELIGMAN

Serial Number: 80-00173a

Core Version: 8105.48

06/28/2015

14:45:14

Standard Lot#: 33913100A3

Last Changed By: J. MURRAY5597 #N56

QAS: J. MURRAY5597 #N56

YAVAPAI SHERIFFS

Operator: J. MURRAY5597 #N56

YAVAPAI SHERIFFS

Test	g/210L	Time
Air Blank	0.000	14:46:00
0.100 Cal Check	0.097	14:46:21
Air Blank	0.000	14:46:51

Preventative Maintenance
Performed -- Timer Reset