

DR 17-023711

RECEIVED

Copy

JUL 06 2017

YAVAPAI COUNTY ATTORNEY

EXHIBIT G-2  
THIS REPORT PREPARED PURSUANT TO DUTY IMPOSED BY A.A.C. R13-10-104 (A)

ARIZONA DEPARTMENT OF PUBLIC SAFETY  
STANDARD QUALITY ASSURANCE PROCEDURES  
INTOXILYZER MODEL 8000  
STANDARD CALIBRATION CHECK PROCEDURE

QA SPECIALIST Scott Blakley Sr AGENCY Yavapai County SO  
DATE 7-4-17 TIME 21:57:13  
INTOXILYZER SERIAL # 80-01734 LOCATION Seligman

- (X) 1. Ensure that gas tank is attached to instrument and contains a standard alcohol concentration solution 0.100 AC.  
OR  
Pour a standard alcohol concentration solution NA AC, into a clean dry simulator and assemble the simulator. Ensure that a tight seal has been made. Turn on the simulator and allow temperature to reach 34°C ± 0.2°C
- (X) 2. Intoxilyzer 8000 display reads "PUSH BUTTON TO START"
- (X) 3. Go to the "Control Testing Menu". Select "D" for dry control test or "W" for wet control test. After selection is made press ENTER.
- (X) 4. Air blank completed.
- (X) 5. Calibration check completed. Test results 0.098 AC.
- (X) 6. Air blank completed.
- (X) 7. Remove printed record. Attach the record to the completed checklist.

SIGNATURE Scott Blakley Sr

DPS Form Exh G-2 (Rev 05-01)

Tank Serial # 22315100A1  
Exp 09-05-2017  
Press. 350 psi

Inspection Beyond 31 day Due to Goodwin Fire  
Instrument Reads "Contact QAS"

DR 17-023711

INTOXILYZER 8000

Location: YCSO SELIGMAN

Serial Number: 80-001734

Core Version: 8105.48

07/04/2017

21:53:35

-----  
Standard Lot#: 22315100A1

Last Changed By: J. MURRAY5597 #N56  
-----

QAS: S. BLAKLEY SR #L6861

YAVAPAI CSO  
-----

Operator: S. BLAKLEY SR #L6861

YAVAPAI CSO  
-----

Test	g/210L	Time
Air Blank	0.000	21:54:00
0.100 Cal Check	0.099	21:54:27
Air Blank	0.000	21:54:57

Preventative Maintenance  
Performed -- Timer Reset

DR 17-023711

INTOXILYZER 8000

Location: ICSO SELIGMAN

Serial Number: 80-001734

Core Version: 8105.48

07/04/2017

21:57:13

-----  
Standard Lot#: 22315100A1

Last Changed By: J. MURRAY5597 #N56

-----  
QAS: S. BLAKLEY SR #L6861

YAVAPAI CSO

-----  
Operator: S. BLAKLEY SR #L6861

YAVAPAI CSO

-----  
Test                      g/210L                      Time  
-----  
Air Blank                      0.000                      21:57:43  
0.100 Cal Check              0.098                      21:58:04  
Air Blank                      0.000                      21:58:33

Preventative Maintenance  
Performed -- Timer Reset

EXHIBIT G-4  
THIS REPORT PREPARED PURSUANT TO DUTY IMPOSED BY A.A.C. R13-10-104 (A)

ARIZONA DEPARTMENT OF PUBLIC SAFETY

STANDARD QUALITY ASSURANCE PROCEDURES  
INTOXILYZER MODEL 8000

STANDARD QUALITY ASSURANCE PROCEDURE

QA SPECIALIST Scott E Blakley Sr. AGENCY Yavapai County SO  
DATE 7-4-17 TIME 22:01:02  
INTOXILYZER SERIAL # 80-01734 LOCATION Seligman

(X) 1. Display Reads "PUSH BUTTON TO START"

DIAGNOSTIC TESTS

- (X) 1. Clock time check.
- (X) 2. Date check.

OPERATIONAL TESTS

- (X) 1. Alcohol-free subject test result 0. 000 AC.
- (k) 2. Error recognition logic system functioning.  
Not a Successfully Completed Test Sequence printed
- (X) 3. Proper sample recognition system.  
Not a Successfully Completed Test Sequence printed  
Deficient sample printed.
- (X) 4. Standard Calibration Check standard 0. 100 AC. Result 0. 098 AC.

Instrument is operating properly and accurately. YES X NO     

COMMENTS

Inspection Past Due - Goodwin Fire  
\_\_\_\_\_  
\_\_\_\_\_  
\_\_\_\_\_

SIGNATURE Scott E Blakley Sr.

INTOXILYZER 8000  
Location: YCSO SELIGMAN  
Serial Number: 80-001734  
Core Version: 8105.48  
07/04/2017

22:01:02

Standard Lot#: 22315100A1  
Last Changed By: J. MURRAY5597 #N56

QAS: S. BLAKLEY SR #L6861  
YAVAPAI CSO

Operator: S. BLAKLEY SR #L6861  
YAVAPAI CSO

Subject: SQAP, TEST  
DOB: 10/20/1960  
Sex: F Weight: 115

15 Minute Deprivation Period? Yes

Test	g/210L	Time
Air Blank	0.000	22:02:22
Diagnostic Test	Pass	22:02:57
Air Blank	0.000	22:03:25
0.100 Cal Check	0.099	22:03:45
Air Blank	0.000	22:04:15
Subject Test	0.000	22:04:33
Air Blank	0.000	22:05:00
Wait	ABT*	
Air Blank	0.000	22:05:35

\*Sequence Aborted

Not a Successfully Completed Test Sequence

125 11-005 111

INTOXILYZER 8000

Location: YCSO SELIGMAN

Serial Number: 80-001734

Core Version: 8105.48

07/04/2017

22:07:08

Standard Lot#: 22315100A1

Last Changed By: J. MURRAY5597 #N56

QAS: S. BLAKLEY SR #L6861

YAVAPAI CSO

Operator: S. BLAKLEY SR #L6861

YAVAPAI CSO

Subject: SQAP, TEST

DOB: 10/20/1960

Sex: F Weight: 115

15 Minute Deprivation Period? Yes

Test	g/210L	Time
Air Blank	IPS*	22:07:55
Air Blank	0.000	22:08:24

\*Improper Sample

Not a Successfully Completed Test Sequence

INTOXILYZER 8000  
Location: YCSO SELIGMAN  
Serial Number: 80-001734  
Core Version: 8105.48  
07/04/2017

22:08:34

Standard Lot#: 22315100A1  
Last Changed By: J. MURRAY5597 #N56

QAS: S. BLAKLEY SR #L6861  
YAVAPAI CSO

Operator: S. BLAKLEY SR #L6861  
YAVAPAI CSO

Subject: SOAP, TEST  
DOB: 10/20/1960  
Sex: F Weight: 115

15 Minute Deprivation Period? Yes

Test	g/210L	Time
Air Blank	RFI*	22:09:41
Air Blank	0.000	22:10:10

\*RFI Detect

**Not a Successfully Completed Test Sequence**

DR 17-023711

INTOXILYZER 8000

Location: YCSO SELIGMAN

Serial Number: 80-001734

Core Version: 8105.48

07/04/2017

22:11:01

---

Standard Lot#: 22315100A1

Last Changed By: J. MURRAY5597 #N56

---

OAS: S. BLAKLEY SR #L6861  
YAVAPAI CSO

---

Operator: S. BLAKLEY SR #L6861  
YAVAPAI CSO

---

Subject: SQAP,TEST

DOB: 10/20/1960

Sex: F Weight: 115

---

15 Minute Deprivation Period? Yes

---

Test	g/210L	Time
Air Blank	0.000	22:11:59
Diagnostic Test	Pass	22:12:34
Air Blank	0.000	22:13:01
0.100 Cal Check	0.098	22:13:22
Air Blank	0.000	22:13:51
Subject Test	0.116*	22:16:59
Air Blank	0.000	22:17:31
Air Blank	0.000	22:17:58
Subject Test	RGE**	22:18:20
Air Blank	0.000	22:18:59

\*Deficient Sample

\*\*Range Exceeded

**Not a Successfully Completed Test Sequence**

DXK 17-02371)

INTOXILYZER 8000

Location: YCSO SELIGMAN

Serial Number: 80-001734

Core Version: 8105.48

07/04/2017

22:19:58

-----  
Standard Lot#: 22315100A1

Last Changed By: J. MURRAY5597 #D56

-----  
QAS: S. BLAKLEY SR #L6861

YAVAPAI CSO

-----  
Operator: S. BLAKLEY SR #L6861

YAVAPAI CSO

-----  
Test                      g/210L                      Time  
-----  
Air Blank                      0.000                      22:20:29  
0.100 Cal Check                0.098                      22:20:49  
Air Blank                      0.000                      22:21:19

Preventative Maintenance  
Performed -- Timer Reset

EXHIBIT G-2  
THIS REPORT PREPARED PURSUANT TO DUTY IMPOSED BY A.A.C. R13-10-104 (A)

ARIZONA DEPARTMENT OF PUBLIC SAFETY  
STANDARD QUALITY ASSURANCE PROCEDURES  
INTOXILYZER MODEL 8000  
STANDARD CALIBRATION CHECK PROCEDURE

QA SPECIALIST Scott Blakley Sr AGENCY Yavapai County SO  
DATE 7-4-17 TIME 22:29:35  
INTOXILYZER SERIAL # 80-01734 LOCATION Seligman

- (X) 1. Ensure that gas tank is attached to instrument and contains a standard alcohol concentration solution 0.100 AC.  
OR  
Pour a standard alcohol concentration solution N/A AC, into a clean dry simulator and assemble the simulator. Ensure that a tight seal has been made. Turn on the simulator and allow temperature to reach 34°C ± 0.2°C
- (X) 2. Intoxilyzer 8000 display reads "PUSH BUTTON TO START"
- (X) 3. Go to the "Control Testing Menu". Select "D" for dry control test or "W" for wet control test. After selection is made press ENTER.
- (X) 4. Air blank completed.
- (X) 5. Calibration check completed. Test results 0. 098 AC.
- (X) 6. Air blank completed.
- (X) 7. Remove printed record. Attach the record to the completed checklist.

SIGNATURE [Signature]

DR 17-023711

INTOXILYZER 8000

Location: YCSO SELIGMAN

Serial Number: 80-001734

Core Version: 8105.48

07/04/2017

22:29:35

-----  
Standard Lot#: 22315100A1

Last Changed By: J. MURRAY5597 #N56  
-----

QAS: S. BLAKLEY SR #L6861

YAVAPAI CSO  
-----

Operator: S. BLAKLEY SR #L6861

YAVAPAI CSO  
-----

Test	g/210L	Time
Air Blank	0.000	22:30:05
0.100 Cal Check	0.098	22:30:26
Air Blank	0.000	22:30:55

Preventative Maintenance  
Performed -- Timer Reset