

EXHIBIT G-2  
THIS REPORT PREPARED PURSUANT TO DUTY IMPOSED BY A.A.C. R13-10-104 (A)

ARIZONA DEPARTMENT OF PUBLIC SAFETY  
STANDARD QUALITY ASSURANCE PROCEDURES  
INTOXILYZER MODEL 8000

STANDARD CALIBRATION CHECK PROCEDURE

QA SPECIALIST Debbie Hughes AGENCY Camp Verde Marshal's Ofc  
DATE 12-18-15 TIME 1033  
INTOXILYZER SERIAL # 80-1736 LOCATION CVMO-Booking Ofc

- (✓) 1. Ensure that gas tank is attached to instrument and contains a standard alcohol concentration solution 0.100 AC.  
OR  
Pour a standard alcohol concentration solution \_\_\_\_\_ AC, into a clean dry simulator and assemble the simulator. Ensure that a tight seal has been made. Turn on the simulator and allow temperature to reach  $34^{\circ}\text{C} \pm 0.2^{\circ}\text{C}$
- (✓) 2. Intoxilyzer 8000 display reads "PUSH BUTTON TO START"
- (✓) 3. Go to the "Control Testing Menu". Select "D" for dry control test or "W" for wet control test. After selection is made press ENTER.
- (✓) 4. Air blank completed.
- (✓) 5. Calibration check completed. Test results 0. 098 AC.
- (✓) 6. Air blank completed.
- (✓) 7. Remove printed record. Attach the record to the completed checklist.

SIGNATURE  #53

INTOXILYZER 8000

Location: CVMO

Serial Number: 80-001736

Core Version: 8105.48

12/18/2015

10:35:58

-----  
Standard Lot#: 14014100A1

Last Changed By: S. SOTO #6432  
-----

QAS: D. HUGHES #53

CAMP VERDE MO  
-----

Operator: D. HUGHES #53

CAMP VERDE MO  
-----

Test	g/210L	Time
Air Blank	0.000	10:36:32
0.100 Cal Check	0.098	10:36:53
Air Blank	0.000	10:37:22

Preventative Maintenance  
Performed -- Timer Reset