

EXHIBIT G-2
THIS REPORT PREPARED PURSUANT TO DUTY IMPOSED BY A.A.C. R13-10-104 (A)

ARIZONA DEPARTMENT OF PUBLIC SAFETY
STANDARD QUALITY ASSURANCE PROCEDURES
INTOXILYZER MODEL 8000

STANDARD CALIBRATION CHECK PROCEDURE

QA SPECIALIST Sgt. D. Jacobs AGENCY Camp Verde Marshals
DATE 3/6/17 TIME 1159
INTOXILYZER SERIAL # 80-001736 LOCATION CVMO Evidence Room

- (/) 1. Ensure that gas tank is attached to instrument and contains a standard alcohol concentration solution 0.100 AC.
OR
Pour a standard alcohol concentration solution _____ AC, into a clean dry simulator and assemble the simulator. Ensure that a tight seal has been made. Turn on the simulator and allow temperature to reach $34^{\circ}\text{C} \pm 0.2^{\circ}\text{C}$
- (/) 2. Intoxilyzer 8000 display reads "PUSH BUTTON TO START"
(/) 3. Go to the "Control Testing Menu". Select "D" for dry control test or "W" for wet control test. After selection is made press ENTER.
(/) 4. Air blank completed.
(/) 5. Calibration check completed. Test results 0.098 AC.
(X) 6. Air blank completed.
(X) 7. Remove printed record. Attach the record to the completed checklist.

SIGNATURE  CV08

INTOXILYZER 8000

Location: CVMO

Serial Number: 80-001736

Core Version: 8105.48

03/06/2017

11:59:13

Standard Lot#: 09216100A1

Last Changed By: D. HUGHES #53

QAS: D. JACOBS #CV08

CAMP VERDE MO

Operator: D. JACOBS #CV08

CAMP VERDE MO

Test	g/210L	Time
Air Blank	0.000	11:59:47
0.100 Cal Check	0.098	12:00:08
Air Blank	0.000	12:00:37

Preventative Maintenance
Performed -- Timer Reset

EXHIBIT G-4
THIS REPORT PREPARED PURSUANT TO DUTY IMPOSED BY A.A.C. R13-10-104 (A)

ARIZONA DEPARTMENT OF PUBLIC SAFETY

STANDARD QUALITY ASSURANCE PROCEDURES
INTOXILYZER MODEL 8000

STANDARD QUALITY ASSURANCE PROCEDURE

QA SPECIALIST Sgt. D. Jacobs AGENCY Camp Verde Marshals
DATE 3/6/17 TIME 1201
INTOXILYZER SERIAL # 80-001736 LOCATION Cvma Evidence Room

(/) 1. Display Reads "PUSH BUTTON TO START"

DIAGNOSTIC TESTS

- (/) 1. Clock time check.
(/) 2. Date check.

OPERATIONAL TESTS

- (X) 1. Alcohol-free subject test result 0. 000 AC.
(X) 2. Error recognition logic system functioning.
Not a Successfully Completed Test Sequence printed
(/) 3. Proper sample recognition system.
Not a Successfully Completed Test Sequence printed
Deficient sample printed.
(X) 4. Standard Calibration Check standard 0. 100 AC. Result 0. 098 AC.

Instrument is operating properly and accurately. YES X NO

COMMENTS

SIGNATURE  CV08

INTOXILYZER 8000
Location: CVMO
Serial Number: 80-001736
Core Version: 8105.48
03/06/2017 12:01:07

Standard Lot#: 09216100A1
Last Changed By: D. HUGHES #53

QAS: D. JACOBS #CV08
CAMP VERDE MO

Operator: D. JACOBS #CV08
CAMP VERDE MO

Subject: TEST,TEST
DOB: 01/02/1989
Sex: M Weight: 200

15 Minute Deprivation Period? Yes

Test	g/210L	Time
Air Blank	0.000	12:02:09
Diagnostic Test	Pass	12:02:44
Air Blank	0.000	12:03:11
0.100 Cal Check	0.098	12:03:32
Air Blank	0.000	12:04:01
Subject Test	0.000	12:04:19
Air Blank	ABT*	12:04:30
Air Blank	IPS**	12:04:38

*Sequence Aborted
**Improper Sample

Not a Successfully Completed Test Sequence

INTOXILYZER 8000
Location: CVMO
Serial Number: 80-001736
Core Version: 8105.48
03/06/2017 12:05:28

Standard Lot#: 09216100A1
Last Changed By: D. HUGHES #53

QAS: D. JACOBS #CV08
CAMP VERDE MO

Operator: D. JACOBS #CV08
CAMP VERDE MO

Subject: TEST,TEST
DOB: 01/02/1989
Sex: M Weight: 200

15 Minute Deprivation Period? Yes

Test	g/210L	Time
Air Blank	0.000	12:06:28
Diagnostic Test	Pass	12:07:03
Air Blank	0.000	12:07:30
0.100 Cal Check	0.098	12:07:51
Air Blank	0.000	12:08:21
Subject Test	0.026*	12:11:28
Air Blank	0.000	12:11:59
Air Blank	0.000	12:12:27
Subject Test	XXX**	12:13:06
Air Blank	RFI***	12:13:16

*Deficient Sample
**Invalid Sample
 Begin new deprivation period
***RFI Detect

Not a Successfully Completed Test Sequence

EXHIBIT G-2
THIS REPORT PREPARED PURSUANT TO DUTY IMPOSED BY A.A.C. R13-10-104 (A)

ARIZONA DEPARTMENT OF PUBLIC SAFETY
STANDARD QUALITY ASSURANCE PROCEDURES
INTOXILYZER MODEL 8000

STANDARD CALIBRATION CHECK PROCEDURE

QA SPECIALIST Sgt. D. Jacobs AGENCY Camp Verde Marshals
DATE 3/6/17 TIME 1214
INTOXILYZER SERIAL # 80-001736 LOCATION CVMA Evidence Room

- (X) 1. Ensure that gas tank is attached to instrument and contains a standard alcohol concentration solution 0.100 AC.
OR
Pour a standard alcohol concentration solution 0.100 AC, into a clean dry simulator and assemble the simulator. Ensure that a tight seal has been made. Turn on the simulator and allow temperature to reach $34^{\circ}\text{C} \pm 0.2^{\circ}\text{C}$
- (X) 2. Intoxilyzer 8000 display reads "PUSH BUTTON TO START"
- (X) 3. Go to the "Control Testing Menu". Select "D" for dry control test or "W" for wet control test. After selection is made press ENTER.
- (X) 4. Air blank completed.
- (X) 5. Calibration check completed. Test results 0. 098 AC.
- (X) 6. Air blank completed.
- (X) 7. Remove printed record. Attach the record to the completed checklist.

SIGNATURE [Signature] CV08

INTOXILYZER 8000

Location: CVMO

Serial Number: 80-001736

Core Version: 8105.48

03/06/2017

12:14:02

Standard Lot#: 09216100A1

Last Changed By: D. HUGHES #53

QAS: D. JACOBS #CV08

CAMP VERDE MO

Operator: D. JACOBS #CV08

CAMP VERDE MO

Test	g/210L	Time
Air Blank	0.000	12:14:36
0.100 Cal Check	0.098	12:14:57
Air Blank	0.000	12:15:26

Preventative Maintenance
Performed -- Timer Reset