

COPY

EXHIBIT G-2  
THIS REPORT PREPARED PURSUANT TO DUTY IMPOSED BY A.A.C. R13-10-104 (A)

ARIZONA DEPARTMENT OF PUBLIC SAFETY  
STANDARD QUALITY ASSURANCE PROCEDURES  
INTOXILYZER MODEL 8000

STANDARD CALIBRATION CHECK PROCEDURE

QA SPECIALIST J. Murray N56 AGENCY Yavapai County Sheriff's  
DATE 9/25/11 TIME 13:37:04  
INTOXILYZER SERIAL # 80-000289 LOCATION Pulscott

- (X) 1. Ensure that gas tank is attached to instrument and contains a standard alcohol concentration solution .100 AC.  
OR  
Pour a standard alcohol concentration solution        AC, into a clean dry simulator and assemble the simulator. Ensure that a tight seal has been made. Turn on the simulator and allow temperature to reach 34°C ± 0.2°C
- (+) 2. Intoxilyzer 8000 display reads "PUSH BUTTON TO START"
- (X) 3. Go to the "Control Testing Menu". Select "D" for dry control test or "W" for wet control test. After selection is made press ENTER.
- (+) 4. Air blank completed.
- (+) 5. Calibration check completed. Test results 0. 094 AC.
- (+) 6. Air blank completed.
- (+) 7. Remove printed record. Attach the record to the completed checklist.

SIGNATURE J. Murray N56

INTOXILYZER 8000

Location: YAVAPAI CO. SHERIFF

Serial Number: 80-000289

Core Version: 8105.48

09/25/2011

13:37:04

-----  
Standard Lot#: 01011100A1

Last Changed By: J. MURRAY #N56  
-----

QAS: J. MURRAY #N56

YAVAPAI YCSO  
-----

Operator: J. MURRAY #N56

YAVAPAI YCSO  
-----

Test	g/210L	Time
Air Blank	0.000	13:37:55
0.100 Cal Check	0.094	13:38:15
Air Blank	0.000	13:38:45

Preventative Maintenance  
Performed -- Timer Reset

EXHIBIT G-4  
THIS REPORT PREPARED PURSUANT TO DUTY IMPOSED BY A.A.C. R13-10-104 (A)

ARIZONA DEPARTMENT OF PUBLIC SAFETY  
STANDARD QUALITY ASSURANCE PROCEDURES  
INTOXILYZER MODEL 8000

STANDARD QUALITY ASSURANCE PROCEDURE

QA SPECIALIST JMurray NSB AGENCY Yavapai County Sheriff's  
DATE 9/25/11 TIME 13:39:29  
INTOXILYZER SERIAL # 80-000289 LOCATION Prescott

1. Display Reads "PUSH BUTTON TO START"

DIAGNOSTIC TESTS

1. Clock time check.  
 2. Date check.

OPERATIONAL TESTS

1. Alcohol-free subject test result 0. 000 AC.  
 2. Error recognition logic system functioning.  
Not a Successfully Completed Test Sequence printed  
 3. Proper sample recognition system.  
Not a Successfully Completed Test Sequence printed  
Deficient sample printed.  
 4. Standard Calibration Check standard 0. 100 AC. Result 0. 094 AC.

Instrument is operating properly and accurately. YES  NO

COMMENTS

Trk# 01011100A1 1000 PSL  
\_\_\_\_\_  
\_\_\_\_\_  
\_\_\_\_\_

SIGNATURE JMurray NSB

INTOXILYZER 8000  
Location: YAVAPAI CO. SHERIFF  
Serial Number: 80-000289  
Core Version: 8105.48  
09/25/2011 13:39:29

---

Standard Lot#: 01011100A1  
Last Changed By: J. MURRAY #N56

---

QAS: J. MURRAY #N56  
YAVAPAI YCSO

---

Operator: J. MURRAY #N56  
YAVAPAI YCSO

---

Subject: SQAP,TEST  
DOB: 10/20/1960  
Sex: F Weight: 115

---

15 Minute Deprivation Period? Yes

---

Test	g/210L	Time
Air Blank	0.000	13:40:29
Diagnostic Test	Pass	13:41:04
Air Blank	0.000	13:41:32
0.100 Cal Check	0.094	13:41:53
Air Blank	0.000	13:42:22
Subject Test	0.000	13:43:36
Air Blank	0.000	13:44:03
Wait	ABT*	
Air Blank	RFI**	13:44:37

\*Sequence Aborted

\*\*RFI Detect

Not a Successfully Completed Test Sequence

INTOXILYZER 8000  
Location: YAVAPAI CO. SHERIFF  
Serial Number: 80-000289  
Core Version: 8105.48  
09/25/2011 13:44:41

---

Standard Lot#: 01011100A1  
Last Changed By: J. MURRAY #N56

---

QAS: J. MURRAY #N56  
YAVAPAI YCSO

---

Operator: J. MURRAY #N56  
YAVAPAI YCSO

---

Subject: SQAP,TEST  
DOB: 10/20/1960  
Sex: F Weight: 115

---

15 Minute Deprivation Period? Yes

---

Test	g/210L	Time
Air Blank	IPS*	13:45:27
Air Blank	RFI**	13:45:39

\*Improper Sample  
\*\*RFI Detect

Not a Successfully Completed Test Sequence

INTOXILYZER 8000  
Location: YAVAPAI CO. SHERIFF  
Serial Number: 80-000289  
Core Version: 8105.48  
09/25/2011 13:46:02

---

Standard Lot#: 01011100A1  
Last Changed By: J. MURRAY #N56

---

QAS: J. MURRAY #N56  
YAVAPAI YCSO

---

Operator: J. MURRAY #N56  
YAVAPAI YCSO

---

Subject: SQAP,TEST  
DOB: 10/20/1960  
Sex: F Weight: 115

---

15 Minute Deprivation Period? Yes

---

Test	g/210L	Time
Air Blank	0.000	13:46:58
Diagnostic Test	Pass	13:47:33
Air Blank	0.000	13:48:00
0.100 Cal Check	0.094	13:48:21
Air Blank	0.000	13:48:50
Subject Test	0.136*	13:51:58
Air Blank	0.000	13:52:30
Air Blank	0.000	13:52:57
Subject Test	XXX**	13:53:23
Air Blank	0.000	13:53:56

\*Deficient Sample

\*\*Invalid Sample

Begin new deprivation period

**Not a Successfully Completed Test Sequence**

EXHIBIT G-2  
THIS REPORT PREPARED PURSUANT TO DUTY IMPOSED BY A.A.C. R13-10-104 (A)

ARIZONA DEPARTMENT OF PUBLIC SAFETY  
STANDARD QUALITY ASSURANCE PROCEDURES  
INTOXILYZER MODEL 8000

STANDARD CALIBRATION CHECK PROCEDURE

QA SPECIALIST J. Murray NSB AGENCY Yavapai County Sheriff's  
DATE 9/25/11 TIME 13:54:47  
INTOXILYZER SERIAL # 80-000287 LOCATION Prescott

- (X) 1. Ensure that gas tank is attached to instrument and contains a standard alcohol concentration solution 100 AC.  
OR  
Pour a standard alcohol concentration solution        AC, into a clean dry simulator and assemble the simulator. Ensure that a tight seal has been made. Turn on the simulator and allow temperature to reach  $34^{\circ}\text{C} \pm 0.2^{\circ}\text{C}$
- (X) 2. Intoxilyzer 8000 display reads "PUSH BUTTON TO START"
- (X) 3. Go to the "Control Testing Menu". Select "D" for dry control test or "W" for wet control test. After selection is made press ENTER.
- (X) 4. Air blank completed.
- (X) 5. Calibration check completed. Test results 0.094 AC.
- (X) 6. Air blank completed.
- (X) 7. Remove printed record. Attach the record to the completed checklist.

SIGNATURE J. Murray NSB

INTOXILYZER 8000  
Location: YAVAPAI CO. SHERIFF  
Serial Number: 80-000289  
Core Version: 8105.48  
09/25/2011 13:54:47

-----  
Standard Lot#: 01011100A1  
Last Changed By: J. MURRAY #N56  
-----

QAS: J. MURRAY #N56  
YAVAPAI YCSO  
-----

Operator: J. MURRAY #N56  
YAVAPAI YCSO  
-----

Test	g/210L	Time
Air Blank	0.000	13:55:23
0.100 Cal Check	0.094	13:55:44
Air Blank	0.000	13:56:13

Preventative Maintenance  
Performed -- Timer Reset