

EXHIBIT G-2
THIS REPORT PREPARED PURSUANT TO DUTY IMPOSED BY A.A.C. R13-10-104(A)

ARIZONA DEPARTMENT OF PUBLIC SAFETY
STANDARD QUALITY ASSURANCE PROCEDURES
INTOXILYZER MODEL 8000
STANDARD CALIBRATION CHECK PROCEDURE

QA SPECIALIST S. A. Soto #6432 AGENCY Arizona DPS

DATE 07/31/2015 2223
INTOXILYZER SERIAL # 80-000289 LOCATION YCSO Black Canyon City TIME

- (SAP) 1. Ensure that gas tank is attached to instrument and contains a standard alcohol concentration solution 0.100 AC.
OR
Pour a standard alcohol concentration solution _____ AC, into a clean dry simulator and assemble the simulator. Ensure that a tight seal has been made. Turn on the simulator and allow temperature to reach 34° C ± 0.2° C
- (SAP) 2. Intoxilyzer 8000 display reads "PUSH BUTTON TO START"
- (SAP) 3. Go to the "Control Testing Menu". Select "D" for dry control test or "W" for wet control test. After selection is made press ENTER.
- (SAP) 4. Air blank completed.
- (SAP) 5. Calibration check completed. Test results 0. 0.099 AC.
- (SAP) 6. Air blank completed.
- (SAP) 7. Remove printed record. Attach the record to the completed checklist.

SIGNATURE S. A. Soto #6432

INTOXILYZER 6030
Location: YC50 - BLACK CANYON
Serial Number: 86-000289
Core Version: 8105.48
07/31/2015 22:23:12

Standard Lot#: 25414100A4
Last Changed By: L. MCCOY #7486

DAS: S. SOTO #6432
AZ DPS

Operator: S. SOTO #6432
AZ DPS

Test	g/210	Time
Air Blank	0.000	22:23:45
0.100 Cal Check	0.099	22:24:05
Air Blank	0.000	22:24:35

Preventative Maintenance
Performed -- Timer Reset