

EXHIBIT G-2  
THIS REPORT PREPARED PURSUANT TO DUTY IMPOSED BY A.A.C. R13-10-104 (A)  
ARIZONA DEPARTMENT OF PUBLIC SAFETY  
STANDARD QUALITY ASSURANCE PROCEDURES  
INTOXILYZER MODEL 8000  
STANDARD CALIBRATION CHECK PROCEDURE

QA SPECIALIST Logan Moody #7486 AGENCY Yavapai County Sheriff's Office  
DATE 06/22/2017 TIME 2132  
INTOXILYZER SERIAL # 80-000289 LOCATION YCSO - Black Canyon

- 1. Ensure that gas tank is attached to instrument and contains a standard alcohol concentration solution .100 AC.  
OR  
~~Pour a standard alcohol concentration solution \_\_\_\_\_ AC, into a clean dry simulator and assemble the simulator. Ensure that a tight seal has been made. Turn on the simulator and allow temperature to reach 34°C ± 0.2°C~~
- 2. Intoxilyzer 8000 display reads "PUSH BUTTON TO START"
- 3. Go to the "Control Testing Menu". Select "D" for dry control test or "W" for wet control test. After selection is made press ENTER.
- 4. Air blank completed.
- 5. Calibration check completed. Test results 0. 099 AC.
- 6. Air blank completed.
- 7. Remove printed record. Attach the record to the completed checklist.

SIGNATURE Logan Moody



EXHIBIT G-4  
THIS REPORT PREPARED PURSUANT TO DUTY IMPOSED BY A.A.C. R13-10-104 (A)

ARIZONA DEPARTMENT OF PUBLIC SAFETY  
STANDARD QUALITY ASSURANCE PROCEDURES  
INTOXILYZER MODEL 8000

STANDARD QUALITY ASSURANCE PROCEDURE

QA SPECIALIST Logan Moody #7486 AGENCY Yavapai County Sheriff's Office

DATE 06/22/2017 TIME 2134

INTOXILYZER SERIAL # 80-000289 LOCATION YCSO - Black Canyon

1. Display Reads "PUSH BUTTON TO START"

DIAGNOSTIC TESTS

- 1. Clock time check.
- 2. Date check.

OPERATIONAL TESTS

- 1. Alcohol-free subject test result 0. 000 AC.
- 2. Error recognition logic system functioning.  
Not a Successfully Completed Test Sequence printed
- 3. Proper sample recognition system.  
Not a Successfully Completed Test Sequence printed  
Deficient sample printed.
- 4. Standard Calibration Check standard 0. 100 AC. Result 0. 099 AC.

Instrument is operating properly and accurately. YES  NO

COMMENTS

---



---



---



---

SIGNATURE Logan Moody #7486

17-022216

INTOXILYZER 8000  
Location: YCSO - BLACK CANYON  
Serial Number: 80-000289  
Core Version: 8105.48  
06/22/2017 21:34:04

---

Standard Lot#: 23016100A1  
Last Changed By: L. MOODY #7486

---

QAS: L. MOODY #7486  
YAVAPAI CSO

---

Operator: L. MOODY #7486  
YAVAPAI CSO

---

Subject: SQAP,TEST  
DOB: 01/01/1990  
Sex: M Weight: 200

---

15 Minute Deprivation Period? Yes

---

Test	g/210L	Time
Air Blank	0.000	21:34:54
Diagnostic Test	Pass	21:35:29
Air Blank	0.000	21:35:56
0.100 Cal Check	0.100	21:36:17
Air Blank	0.000	21:36:46
Subject Test	0.000	21:37:09
Air Blank	ABT*	21:37:19
Air Blank	RFI**	21:37:26

\*Sequence Aborted  
\*\*RFI Detect

**Not a Successfully Completed Test Sequence**

INTOXILYZER 8000  
 Location: YCSO - BLACK CANYON  
 Serial Number: 80-000289  
 Core Version: 8105.48  
 06/22/2017 21:37:34

---

Standard Lot#: 23016100A1  
 Last Changed By: L. MOODY #7486

---

QAS: L. MOODY #7486  
 YAVAPAI CSO

---

Operator: L. MOODY #7486  
 YAVAPAI CSO

---

Subject: SQAP, TEST  
 DOB: 01/01/1990  
 Sex: M Weight: 200

---

15 Minute Deprivation Period? Yes

---

Test	g/210L	Time
Air Blank	0.000	21:38:24
Diagnostic Test	Pass	21:38:59
Air Blank	0.000	21:39:26
0.100 Cal Check	0.100	21:39:47
Air Blank	0.000	21:40:17
Subject Test	0.117*	21:43:24
Air Blank	0.000	21:43:56
Air Blank	0.000	21:44:23
Subject Test	XXX**	21:44:57
Air Blank	IPS***	21:45:04

\*Deficient Sample

\*\*Invalid Sample

Begin new deprivation period

\*\*\*Improper Sample

**Not a Successfully Completed Test Sequence**

**EXHIBIT G-2  
THIS REPORT PREPARED PURSUANT TO DUTY IMPOSED BY A.A.C. R13-10-104 (A)**

**ARIZONA DEPARTMENT OF PUBLIC SAFETY**

**STANDARD QUALITY ASSURANCE PROCEDURES  
INTOXILYZER MODEL 8000**

**STANDARD CALIBRATION CHECK PROCEDURE**

QA SPECIALIST Logan Moody #7486 AGENCY Yavapai County Sheriff's Office  
 DATE 06/22/2017 TIME 2145  
 INTOXILYZER SERIAL # 80-000289 LOCATION YCSO - Black Canyon

- (C) 1. Ensure that gas tank is attached to instrument and contains a standard alcohol concentration solution 100 AC.  
 OR  
 Pour a standard alcohol concentration solution \_\_\_\_\_ AC, into a clean dry simulator and assemble the simulator. Ensure that a tight seal has been made. Turn on the simulator and allow temperature to reach  $34^{\circ}\text{C} \pm 0.2^{\circ}\text{C}$
- (C) 2. Intoxilyzer 8000 display reads "PUSH BUTTON TO START"
- (C) 3. Go to the "Control Testing Menu". Select "D" for dry control test or "W" for wet control test. After selection is made press ENTER.
- (C) 4. Air blank completed.
- (C) 5. Calibration check completed. Test results 0. 099 AC.
- (C) 6. Air blank completed.
- (C) 7. Remove printed record. Attach the record to the completed checklist.

SIGNATURE Logan Moody

