

EXHIBIT G-2
THIS REPORT PREPARED PURSUANT TO DUTY IMPOSED BY A.A.C. R13-10-104 (A)

ARIZONA DEPARTMENT OF PUBLIC SAFETY
STANDARD QUALITY ASSURANCE PROCEDURES
INTOXILYZER MODEL 8000
STANDARD CALIBRATION CHECK PROCEDURE

QA SPECIALIST Logan Moody #7486 AGENCY Yavapai County Sheriff's Office

DATE 01/05/2018 TIME 1417

INTOXILYZER SERIAL # 80-000289 LOCATION YCSO - Black Canyon

- (4) 1. Ensure that gas tank is attached to instrument and contains a standard alcohol concentration solution, 100 AC.
OR
Pour a standard alcohol concentration solution _____ AC, into a clean dry simulator and assemble the simulator. ~~Ensure that a tight seal has been made.~~ Turn on the simulator and allow temperature to reach 34°C ± 0.2°C
- (4) 2. Intoxilyzer 8000 display reads "PUSH BUTTON TO START"
- (4) 3. Go to the "Control Testing Menu". Select "D" for dry control test or "W" for wet control test. After selection is made press ENTER.
- (4) 4. Air blank completed.
- (4) 5. Calibration check completed. Test results 0. 097 AC.
- (4) 6. Air blank completed.
- (4) 7. Remove printed record. Attach the record to the completed checklist.

SIGNATURE Logan Moody #7486

18-000515

INTOXILYZER 8000
Location: YCSO - BLACK CANYON
Serial Number: 80-000289
Core Version: 8105.48
01/05/2018 14:17:11

Standard Lot#: 23016100A1
Last Changed By: L. MOODY #7486

QAS: L. MOODY #7486
YAVAPAI CSO

Operator: L. MOODY #7486
YAVAPAI CSO

Test	g/210L	Time
Air Blank	0.000	14:17:39
0.100 Cal Check	0.097	14:18:00
Air Blank	0.000	14:18:30

Preventative Maintenance
Performed -- Timer Reset

18-000515

EXHIBIT G-2
THIS REPORT PREPARED PURSUANT TO DUTY IMPOSED BY A.A.C. R13-10-104 (A)

ARIZONA DEPARTMENT OF PUBLIC SAFETY
STANDARD QUALITY ASSURANCE PROCEDURES
INTOXILYZER MODEL 8000

STANDARD CALIBRATION CHECK PROCEDURE

QA SPECIALIST Logan Moody #7486 AGENCY Yavapai County Sheriff's Office

DATE 01/05/2018 TIME 1403

INTOXILYZER SERIAL # 80-000289 LOCATION YCSO - Black Canyon

- (18) 1. Ensure that gas tank is attached to instrument and contains a standard alcohol concentration solution .100 AC.
OR
Pour a standard alcohol concentration solution _____ AC, into a clean dry simulator and assemble the simulator. ~~Ensure that a tight seal has been made.~~ Turn on the simulator and allow temperature to reach $34^{\circ}\text{C} \pm 0.2^{\circ}\text{C}$
- (18) 2. Intoxilyzer 8000 display reads "PUSH BUTTON TO START"
- (18) 3. Go to the "Control Testing Menu". Select "D" for dry control test or "W" for wet control test. After selection is made press ENTER.
- (18) 4. Air blank completed.
- (18) 5. Calibration check completed. Test results 0. 097 AC.
- (18) 6. Air blank completed.
- (18) 7. Remove printed record. Attach the record to the completed checklist.

SIGNATURE Logan Moody #7486

18-000515

INTOXILYZER 8000

Location: YCSO - BLACK CANYON

Serial Number: 80-000289

Core Version: 8105.48

01/05/2018

14:03:34

Standard Lot#: 23016100A1

Last Changed By: L. MOODY #7486

QAS: L. MOODY #7486

YAVAPAI CSO

Operator: L. MOODY #7486

YAVAPAI CSO

Test	g/210L	Time
Air Blank	0.000	14:04:02
0.100 Cal Check	0.097	14:04:23
Air Blank	0.000	14:04:52

Preventative Maintenance
Performed -- Timer Reset

EXHIBIT G-4
THIS REPORT PREPARED PURSUANT TO DUTY IMPOSED BY A.A.C. R13-10-104 (A)

ARIZONA DEPARTMENT OF PUBLIC SAFETY

STANDARD QUALITY ASSURANCE PROCEDURES
INTOXILYZER MODEL 8000

STANDARD QUALITY ASSURANCE PROCEDURE

QA SPECIALIST Logan Moody #7486 AGENCY Yavapai County Sheriff's Office

DATE 01/05/2018 TIME 1405

INTOXILYZER SERIAL # 80-000289 LOCATION YCSO - Black Canyon

1. Display Reads "PUSH BUTTON TO START"

DIAGNOSTIC TESTS

- 1. Clock time check.
- 2. Date check.

OPERATIONAL TESTS

- 1. Alcohol-free subject test result 0. 000 AC.
- 2. Error recognition logic system functioning.
Not a Successfully Completed Test Sequence printed
- 3. Proper sample recognition system.
Not a Successfully Completed Test Sequence printed
Deficient sample printed.
- 4. Standard Calibration Check standard 0. 100 AC. Result 0. 097 AC.

Instrument is operating properly and accurately. YES NO

COMMENTS

950 psi

SIGNATURE Logan Moody 7486

INTOXILYZER 8000
Location: YCSO - BLACK CANYON
Serial Number: 80-000289
Core Version: 8105.48
01/05/2018 14:05:10

Standard Lot#: 23016100A1
Last Changed By: L. MOODY #7486

QAS: L. MOODY #7486
YAVAPAI CSO

Operator: L. MOODY #7486
YAVAPAI CSO

Subject: SQAP,TEST
DOB: 01/01/1990
Sex: M Weight: 200

15 Minute Deprivation Period? Yes

Test	g/210L	Time
Air Blank	0.000	14:06:07
Diagnostic Test	Pass	14:06:42
Air Blank	0.000	14:07:09
0.100 Cal Check	0.097	14:07:30
Air Blank	0.000	14:07:59
Subject Test	0.000	14:08:21
Air Blank	ABT*	14:08:28
Air Blank	RFI**	14:08:36

*Sequence Aborted
**RFI Detect

Not a Successfully Completed Test Sequence

INTOXILYZER 8000

Location: YCSO - BLACK CANYON

Serial Number: 80-000289

Core Version: 8105.48

01/05/2018

14:08:42

Standard Lot#: 23016100A1

Last Changed By: L. MOODY

#7486

QAS: L. MOODY

#7486

YAVAPAI CSO

Operator: L. MOODY

#7486

YAVAPAI CSO

Subject: SQAP,TEST

DOB: 01/01/1990

Sex: M

Weight: 200

15 Minute Deprivation Period? Yes

Test	g/210L	Time
Air Blank	0.000	14:09:26
Diagnostic Test	Pass	14:10:01
Air Blank	0.000	14:10:29
0.100 Cal Check	0.097	14:10:50
Air Blank	0.000	14:11:19
Subject Test	0.049*	14:14:26
Air Blank	0.000	14:14:57
Air Blank	0.000	14:15:25
Subject Test	XXX**	14:15:44
Air Blank	IPS***	14:15:51

*Deficient Sample

**Invalid Sample

Begin new deprivation period

***Improper Sample

Not a Successfully Completed Test Sequence

INTOXILYZER 8000

Location: YCSO - BLACK CANYON

Serial Number: 80-000289

Core Version: 8105.48

01/05/2018

14:16:04

Standard Lot#: 23016100A1

Last Changed By: L. MOODY #7486

QAS: L. MOODY #7486

YAVAPAI CSO

Operator: L. MOODY #7486YAVAPAI CSO

Test	g/210L	Time
Air Blank	0.000	14:16:32
Reference	UNS*	14:16:40
Air Blank	0.000	14:17:08

*Unstable Signal