

EXHIBIT G-2
THIS REPORT PREPARED PURSUANT TO DUTY IMPOSED BY A.A.C. R13-10-104 (A)

ARIZONA DEPARTMENT OF PUBLIC SAFETY
STANDARD QUALITY ASSURANCE PROCEDURES
INTOXILYZER MODEL 8000

STANDARD CALIBRATION CHECK PROCEDURE

QA SPECIALIST Tobin #73 AGENCY PVPD
DATE 01-13-2008 TIME 0625
INTOXILYZER SERIAL # 80-000293 LOCATION PVPD

- JPT 1. Ensure that gas tank is attached to instrument and contains a standard alcohol concentration solution .100 AC.
OR
Pour a standard alcohol concentration solution N/A AC, into a clean dry simulator and assemble the simulator. Ensure that a tight seal has been made. Turn on the simulator and allow temperature to reach $34^{\circ}\text{C} \pm 0.2^{\circ}\text{C}$
- JPT 2. Intoxilyzer 8000 display reads "PUSH BUTTON TO START"
- JPT 3. Go to the "Control Testing Menu". Select "D" for dry control test or "W" for wet control test. After selection is made press ENTER.
- JPT 4. Air blank completed.
- JPT 5. Calibration check completed. Test results 0. 098 AC.
- JPT 6. Air blank completed.
- JPT 7. Remove printed record. Attach the record to the completed checklist.

SIGNATURE J. Smith #73

INTOXILYZER 8000

Location: PRESCOTT VALLEY PD

Serial Number: 80-000293

01/13/2008

06:27:59

Standard Lot#: 7065021

Last Changed By: J. TOBIN #73

QAS: J. TOBIN #73

PRESCOTT VALLEY PD

Test g/210L Time

Air Blank 0.000 06:28:27
Cal Check 0.098 06:28:46
Air Blank 0.000 06:29:15

J. Tobin #73

Operator's Signature