

EXHIBIT G-2
THIS REPORT PREPARED PURSUANT TO DUTY IMPOSED BY A.A.C. R13-10-104 (A)

ARIZONA DEPARTMENT OF PUBLIC SAFETY
STANDARD QUALITY ASSURANCE PROCEDURES
INTOXILYZER MODEL 8000

STANDARD CALIBRATION CHECK PROCEDURE

QA SPECIALIST C. Krogedal AGENCY PVPD
DATE 5/11/10 TIME 0403
INTOXILYZER SERIAL # 80-000293 LOCATION Prescott Valley

- (CK) 1. Ensure that gas tank is attached to instrument and contains a standard alcohol concentration solution .100 AC.
OR
Pour a standard alcohol concentration solution _____ AC, into a clean dry simulator and assemble the simulator. Ensure that a tight seal has been made. Turn on the simulator and allow temperature to reach $34^{\circ}\text{C} \pm 0.2^{\circ}\text{C}$
- (CK) 2. Intoxilyzer 8000 display reads "PUSH BUTTON TO START"
- (CK) 3. Go to the "Control Testing Menu". Select "D" for dry control test or "W" for wet control test. After selection is made press ENTER.
- (CK) 4. Air blank completed.
- (CK) 5. Calibration check completed. Test results 0. 096 AC.
- (CK) 6. Air blank completed.
- (CK) 7. Remove printed record. Attach the record to the completed checklist.

SIGNATURE C. Krogedal #90 / QAS #33200

INTOXILYZER 8000

Location: PRESCOTT VALLEY PD

Serial Number: 80-000293

Core Version: 8105.48

05/01/2010

04:05:34

Standard Lot#: 585220

Last Changed By: A. LAWRENCE #61

QAS: C. KROGEDAL #90

PRESCOTT VALLEY PD

Operator: C. KROGEDAL #90

PRESCOTT VALLEY PD

Test	g/210L	Time
Air Blank	0.000	04:06:32
0.100 Cal Check	0.096	04:06:53
Air Blank	0.000	04:07:22

Preventative Maintenance
Performed -- Timer Reset

C. Krogedal #90

QAS #33200