

EXHIBIT G-2
THIS REPORT PREPARED PURSUANT TO DUTY IMPOSED BY A.A.C. R13-10-104 (A)

ARIZONA DEPARTMENT OF PUBLIC SAFETY
STANDARD QUALITY ASSURANCE PROCEDURES
INTOXILYZER MODEL 8000

STANDARD CALIBRATION CHECK PROCEDURE

QA SPECIALIST C. Krogedal #90 AGENCY PVPD
DATE 8/25/10 TIME 2330
INTOXILYZER SERIAL # 80-000293 LOCATION PVPD Booking Room

- () 1. Ensure that gas tank is attached to instrument and contains a standard alcohol concentration solution .100 AC.
OR
Pour a standard alcohol concentration solution _____ AC, into a clean dry simulator and assemble the simulator. Ensure that a tight seal has been made. Turn on the simulator and allow temperature to reach $34^{\circ}\text{C} \pm 0.2^{\circ}\text{C}$
- () 2. Intoxilyzer 8000 display reads "PUSH BUTTON TO START"
- () 3. Go to the "Control Testing Menu". Select "D" for dry control test or "W" for wet control test. After selection is made press ENTER.
- () 4. Air blank completed.
- () 5. Calibration check completed. Test results 0. 096 AC.
- () 6. Air blank completed.
- () 7. Remove printed record. Attach the record to the completed checklist.

SIGNATURE

C. Krogedal #90 / QAS # 33200

INTOXILYZER 8000
Location: PRESCOTT VALLEY PD
Serial Number: 80-000293
Core Version: 8105.48
08/25/2010 23:28:08

Standard Lot#: 685215
Last Changed By: C. KROGEDAL #90

QAS: C. KROGEDAL #90
PRESCOTT VALLEY PD

Operator: C. KROGEDAL #90
PRESCOTT VALLEY PD

Test	g/210L	Time
Air Blank	0.000	23:28:44
0.100 Cal Check	0.096	23:29:04
Air Blank	0.000	23:29:34

Preventative Maintenance
Performed -- Timer Reset

C. Krogedal #90
PVPD QAS # 33200