

EXHIBIT G-2  
THIS REPORT PREPARED PURSUANT TO DUTY IMPOSED BY A.A.C. R13-10-104 (A)

ARIZONA DEPARTMENT OF PUBLIC SAFETY  
STANDARD QUALITY ASSURANCE PROCEDURES  
INTOXILYZER MODEL 8000

STANDARD CALIBRATION CHECK PROCEDURE

QA SPECIALIST C. Krogedal AGENCY PVPD  
DATE 5/13/11 TIME 1400 2155  
INTOXILYZER SERIAL # 80-000293 LOCATION PVPD

- (~~CA~~) 1. Ensure that gas tank is attached to instrument and contains a standard alcohol concentration solution 100 AC.  
OR  
Pour a standard alcohol concentration solution \_\_\_\_\_ AC, into a clean dry simulator and assemble the simulator. Ensure that a tight seal has been made. Turn on the simulator and allow temperature to reach  $34^{\circ}\text{C} \pm 0.2^{\circ}\text{C}$
- (~~CA~~) 2. Intoxilyzer 8000 display reads "PUSH BUTTON TO START"
- (~~CA~~) 3. Go to the "Control Testing Menu". Select "D" for dry control test or "W" for wet control test. After selection is made press ENTER.
- (~~CA~~) 4. Air blank completed.
- (~~CA~~) 5. Calibration check completed. Test results 0. 095 AC.
- (~~CA~~) 6. Air blank completed.
- (~~CA~~) 7. Remove printed record. Attach the record to the completed checklist.

SIGNATURE C. Krogedal #90 / QAS # 33200

INTOXILYZER 8000

Location: PRESCOTT VALLEY PD

Serial Number: 80-000293

Core Version: 8105.48

09/13/2011

21:57:32

-----  
Standard Lot#: 779386

Last Changed By: C. KROGEDAL #90  
-----

QAS: C. KROGEDAL #90

PRESCOTT VALLEY PD  
-----

Operator: C. KROGEDAL #90

PRESCOTT VALLEY PD  
-----

Test	g/210L	Time
Air Blank	0.000	21:58:04
0.100 Cal Check	0.095	21:58:25
Air Blank	0.000	21:58:54

Preventative Maintenance  
Performed -- Timer Reset