

EXHIBIT G-2
THIS REPORT PREPARED PURSUANT TO DUTY IMPOSED BY A.A.C. R13-10-104 (A)

ARIZONA DEPARTMENT OF PUBLIC SAFETY
STANDARD QUALITY ASSURANCE PROCEDURES
INTOXILYZER MODEL 8000
STANDARD CALIBRATION CHECK PROCEDURE

QA SPECIALIST C. Krogel AGENCY PVPD

DATE 10/17/11 TIME 1800

INTOXILYZER SERIAL # 20-000293 LOCATION PVPD

1. Ensure that gas tank is attached to instrument and contains a standard alcohol concentration solution _____ AC.
OR

Pour a standard alcohol concentration solution _____ AC, into a clean dry simulator and assemble the simulator. Ensure that a tight seal has been made. Turn on the simulator and allow temperature to reach $34^{\circ}\text{C} \pm 0.2^{\circ}\text{C}$

2. Intoxilyzer 8000 display reads "PUSH BUTTON TO START"

3. Go to the "Control Testing Menu". Select "D" for dry control test or "W" for wet control test. After selection is made press ENTER.

4. Air blank completed.

5. Calibration check completed. Test results 0. 095 AC.

6. Air blank completed.

7. Remove printed record. Attach the record to the completed checklist.

SIGNATURE C. Krogel #90 / OAS # 33200

INTOXILYZER 8000

Location: PRESCOTT VALLEY PD

Serial Number: 80-000293

Core Version: 8105.48

10/17/2011

18:00:19

Standard Lot#: 779386

Last Changed By: C. KROGEDAL #90

QAS: C. KROGEDAL #90
PRESCOTT VALLEY PD

Operator: C. KROGEDAL #90
PRESCOTT VALLEY PD

Test	g/210L	Time
Air Blank	0.000	18:00:54
0.100 Cal Check	0.095	18:01:15
Air Blank	0.000	18:01:44

Preventative Maintenance
Performed -- Timer Reset