


EXHIBIT G-2  
THIS REPORT PREPARED PURSUANT TO DUTY IMPOSED BY A.A.C. R13-10-104 (A)

ARIZONA DEPARTMENT OF PUBLIC SAFETY  
STANDARD QUALITY ASSURANCE PROCEDURES  
INTOXILYZER MODEL 8000

STANDARD CALIBRATION CHECK PROCEDURE

QA SPECIALIST G. KOESTER AGENCY PRESCOTT PD.  
DATE 7-3-08 TIME 1348  
INTOXILYZER SERIAL # 80-000337 LOCATION PRESCOTT PD

- (GK) 1. Ensure that gas tank is attached to instrument and contains a standard alcohol concentration solution 0.100 AC.  
OR  
Pour a standard alcohol concentration solution \_\_\_\_\_ AC, into a clean dry simulator and assemble the simulator. Ensure that a tight seal has been made. Turn on the simulator and allow temperature to reach  $34^{\circ}\text{C} \pm 0.2^{\circ}\text{C}$
- (GK) 2. Intoxilyzer 8000 display reads "PUSH BUTTON TO START"
- (GK) 3. Go to the "Control Testing Menu". Select "D" for dry control test or "W" for wet control test. After selection is made press ENTER.
- (GK) 4. Air blank completed.
- (GK) 5. Calibration check completed. Test results 0.098 AC.
- (GK) 6. Air blank completed.
- (GK) 7. Remove printed record. Attach the record to the completed checklist.

SIGNATURE G. KOESTER 

INTOXILYZER 8000

Location: PRESCOTT P.D

Serial Number: 80-000337

07/03/2008

13:48:29


-----  
Standard Lot#: 727501I

Last Changed By: T. ATTENBERGER #252  
-----

QAS: G. KOESTER #212

PRESCOTT PD  
-----

Test	g/210L	Time
Air Blank	0.000	13:48:57
Cal Check	0.098	13:49:16
Air Blank	0.000	13:49:45

G. KOESTER 

Operator's Signature

EXHIBIT G-4  
THIS REPORT PREPARED PURSUANT TO DUTY IMPOSED BY A.A.C. R13-10-104 (A)

ARIZONA DEPARTMENT OF PUBLIC SAFETY  
STANDARD QUALITY ASSURANCE PROCEDURES  
INTOXILYZER MODEL 8000

STANDARD QUALITY ASSURANCE PROCEDURE

QA SPECIALIST G. KOESTER AGENCY PRESOTT PD  
DATE 7-3-08 TIME 1350  
INTOXILYZER SERIAL # 80-000337 LOCATION PRESOTT PD

1. Display Reads "PUSH BUTTON TO START"

DIAGNOSTIC TESTS

1. Clock time check.  
 2. Date check.

OPERATIONAL TESTS

1. Alcohol-free subject test result 0. 0.00 AC.  
 2. Error recognition logic system functioning.  
Not a Successfully Completed Test Sequence printed  
 3. Proper sample recognition system.  
Not a Successfully Completed Test Sequence printed  
Deficient sample printed.  
( ) 4. Standard Calibration Check standard 0. 0.100 AC. Result 0. 0.097 AC.

Instrument is operating properly and accurately. YES  NO

COMMENTS

ALL OPERATING PROPERLY

SIGNATURE G. KOESTER 

INTOXILYZER 8000  
Location: PRESCOTT P.D  
Serial Number : 80-000337  
07/03/2008 13:50:59

---

Standard Lot#: 727501I  
Last Changed By: T. ATTENBERGER #252

---

QAS: G. KOESTER #212  
PRESCOTT PD

---

Operator: G. KOESTER #212  
PRESCOTT PD

---

Subject: SQAP,TEST  
DOB: 01/01/1901  
Sex: M Weight: 200

---

15 Minute Deprivation Period? Yes

---

Test	g/210L	Time
Air Blank	0.000	13:52:28
Diagnostic Test	Pass	13:52:55
Air Blank	0.000	13:53:22
0.100 Cal Check	0.098	13:53:41
Air Blank	0.000	13:54:11
Subject Test	0.000	13:54:28
Air Blank	0.000	13:54:56
Five Minute Wait	ABT*	
Air Blank	0.000	13:55:31

\*Sequence Aborted

**Not a Successfully Completed Test Sequence**

INTOXILYZER 8000  
Location: PRESCOTT P.D  
Serial Number : 80-000337  
07/03/2008 13:55:39

---

Standard Lot#: 727501I  
Last Changed By: T. ATTENBERGER #252

---

QAS: G. KOESTER #212  
PRESCOTT PD

---

Operator: G. KOESTER #212  
PRESCOTT PD

---

Subject: SQAP,TEST  
DOB: 01/01/1901  
Sex: M Weight: 200

---

15 Minute Deprivation Period? Yes

---

Test	g/210L	Time
Air Blank	IPS*	13:56:16
Air Blank	0.000	13:56:44

\*Improper Sample

**Not a Successfully Completed Test Sequence**

INTOXILYZER 8000  
Location: PRESCOTT P.D  
Serial Number : 80-000337  
07/03/2008 13:56:54

---

Standard Lot#: 727501I  
Last Changed By: T. ATTENBERGER #252

---

QAS: G. KOESTER #212  
PRESCOTT PD

---

Operator: G. KOESTER #212  
PRESCOTT PD

---

Subject: SQAP,TEST  
DOB: 01/01/1901  
Sex: M Weight: 200

---

15 Minute Deprivation Period? Yes

---

Test	g/210L	Time
Air Blank	RFI*	13:57:28
Air Blank	0.000	13:57:56

\*RFI Detect

**Not a Successfully Completed Test Sequence**

INTOXILYZER 8000  
Location: PRESCOTT P.D  
Serial Number : 80-000337  
07/03/2008

13:58:10

---

Standard Lot#: 727501I  
Last Changed By: T. ATTENBERGER #252

---

QAS: G. KOESTER #212  
PRESCOTT PD

---

Operator: G. KOESTER #212  
PRESCOTT PD

---

Subject: SQAP,TEST  
DOB: 01/01/1901  
Sex: M Weight: 200

---

15 Minute Deprivation Period? Yes

---

Test	g/210L	Time
Air Blank	0.000	13:59:10
Diagnostic Test	Pass	13:59:37
Air Blank	0.000	14:00:04
0.100 Cal Check	0.098	14:00:23
Air Blank	0.000	14:00:53
Subject Test	0.039*	14:04:01
Air Blank	0.000	14:04:32
Air Blank	0.000	14:04:59
Subject Test	XXX**	14:05:28
Air Blank	0.000	14:06:00

\*Deficient Sample

\*\*Invalid Sample

Begin new deprivation period


**Not a Successfully Completed Test Sequence**

EXHIBIT G-2  
THIS REPORT PREPARED PURSUANT TO DUTY IMPOSED BY A.A.C. R13-10-104 (A)

ARIZONA DEPARTMENT OF PUBLIC SAFETY  
STANDARD QUALITY ASSURANCE PROCEDURES  
INTOXILYZER MODEL 8000  
STANDARD CALIBRATION CHECK PROCEDURE

QA SPECIALIST G. KOESTER AGENCY PRESCOTT PD.  
DATE 7-3-08 TIME 1408  
INTOXILYZER SERIAL # 80-000337 LOCATION PRESCOTT PD

- (GK) 1. Ensure that gas tank is attached to instrument and contains a standard alcohol concentration solution .100 AC.  
OR  
Pour a standard alcohol concentration solution \_\_\_\_\_ AC, into a clean dry simulator and assemble the simulator. Ensure that a tight seal has been made. Turn on the simulator and allow temperature to reach  $34^{\circ}\text{C} \pm 0.2^{\circ}\text{C}$
- (GK) 2. Intoxilyzer 8000 display reads "PUSH BUTTON TO START"
- (GK) 3. Go to the "Control Testing Menu". Select "D" for dry control test or "W" for wet control test. After selection is made press ENTER.
- (GK) 4. Air blank completed.
- (GK) 5. Calibration check completed. Test results 0. 097 AC.
- (GK) 6. Air blank completed.
- (GK) 7. Remove printed record. Attach the record to the completed checklist.

SIGNATURE G. KOESTER 

INTOXILYZER 8000

Location: PRESCOTT P.D

Serial Number: 80-000337

07/03/2008

14:07:38

-----  
Standard Lot#: 727501I

Last Changed By: T. ATTENBERGER #252  
-----

QAS: G. KOESTER #212

PRESCOTT PD  
-----

Test	g/210L	Time
Air Blank	0.000	14:08:06
Cal Check	0.097	14:08:24
Air Blank	0.000	14:08:54

G. KOESTER [Signature]  
Operator's Signature