

EXHIBIT G-2  
THIS REPORT PREPARED PURSUANT TO DUTY IMPOSED BY A.A.C. R13-10-104 (A)

ARIZONA DEPARTMENT OF PUBLIC SAFETY  
STANDARD QUALITY ASSURANCE PROCEDURES  
INTOXILYZER MODEL 8000

STANDARD CALIBRATION CHECK PROCEDURE

QA SPECIALIST T. ARZENBERG SC #250 AGENCY Prescott P.D.

DATE 092508 TIME 1625

INTOXILYZER SERIAL # 80-000337 LOCATION P.P.D.

- (✓) 1. Ensure that gas tank is attached to instrument and contains a standard alcohol concentration solution 100 AC.  
OR  
Pour a standard alcohol concentration solution \_\_\_\_\_ AC, into a clean dry simulator and assemble the simulator. Ensure that a tight seal has been made. Turn on the simulator and allow temperature to reach  $34^{\circ}\text{C} \pm 0.2^{\circ}\text{C}$
- (✓) 2. Intoxilyzer 8000 display reads "PUSH BUTTON TO START"
- (✓) 3. Go to the "Control Testing Menu". Select "D" for dry control test or "W" for wet control test. After selection is made press ENTER.
- (✓) 4. Air blank completed.
- (✓) 5. Calibration check completed. Test results 0. 097 AC.
- (✓) 6. Air blank completed.
- (✓) 7. Remove printed record. Attach the record to the completed checklist.

SIGNATURE [Signature]

Prescott PD CHANGE 5-12-06

INTOXILYZER 8000

Location: PRESCOTT P.D

Serial Number: 80-000337

09/25/2008

16:25:56

-----  
Standard Lot#: 727501I

Last Changed By: T. ATTENBERGER #252  
-----

QAS: T. ATTENBERGER #252

PRESCOTT PD  
-----

Test	g/210L	Time
Air Blank	0.000	16:26:23
Cal Check	0.097	16:26:42
Air Blank	0.000	16:27:12


  
-----  
Operator's Signature

EXHIBIT G-4  
THIS REPORT PREPARED PURSUANT TO DUTY IMPOSED BY A.A.C. R13-10-104 (A)

ARIZONA DEPARTMENT OF PUBLIC SAFETY

STANDARD QUALITY ASSURANCE PROCEDURES  
INTOXILYZER MODEL 8000

STANDARD QUALITY ASSURANCE PROCEDURE

QA SPECIALIST T. ATTENBERGER #202 AGENCY Prescott P.D.  
DATE 092508 TIME 1628  
INTOXILYZER SERIAL # 60-002337 LOCATION P.P.D.

() 1. Display Reads "PUSH BUTTON TO START"

DIAGNOSTIC TESTS

- () 1. Clock time check.  
() 2. Date check.

OPERATIONAL TESTS

- () 1. Alcohol-free subject test result 0. 000 AC.  
() 2. Error recognition logic system functioning.  
Not a Successfully Completed Test Sequence printed  
() 3. Proper sample recognition system.  
Not a Successfully Completed Test Sequence printed  
Deficient sample printed.  
() 4. Standard Calibration Check standard 0. 100 AC. Result 0. .097 AC.

Instrument is operating properly and accurately. YES  NO

COMMENTS

ALL FUNCTIONING PROPERLY.

SIGNATURE T. Attenger #202

INTOXILYZER 8000  
Location: PRESCOTT P.D  
Serial Number : 80-000337  
09/25/2008 16:28:20

---

Standard Lot#: 727501I  
Last Changed By: T. ATTENBERGER #252

---

QAS: T. ATTENBERGER #252  
PRESCOTT PD

---

Operator: T. ATTENBERGER #252  
PRESCOTT PD

---

Subject: SQAP, TEST  
DOB: 01/01/1900  
Sex: M Weight: 300

---

15 Minute Deprivation Period? Yes

---

Test	g/210L	Time
Air Blank	0.000	16:30:06
Diagnostic Test	Pass	16:30:33
Air Blank	0.000	16:31:00
0.100 Cal Check	0.097	16:31:19
Air Blank	0.000	16:31:49
Subject Test	0.000	16:32:07
Air Blank	0.000	16:32:34
Five Minute Wait	ABT*	
Air Blank	0.000	16:33:09

\*Sequence Aborted

A handwritten signature in black ink, appearing to be 'T. Attenberger', with a large flourish at the end.

**Not a Successfully Completed Test Sequence**

INTOXILYZER 8000  
Location: PRESCOTT P.D  
Serial Number : 80-000337  
09/25/2008

16:33:57

---

Standard Lot#: 727501I  
Last Changed By: T. ATTENBERGER #252

---

QAS: T. ATTENBERGER #252  
PRESCOTT PD

---

Operator: T. ATTENBERGER #252  
PRESCOTT PD

---

Subject: SQAP, TEST  
DOB: 01/01/1900  
Sex: M Weight: 300

---

15 Minute Deprivation Period? Yes

---

Test	g/210L	Time
Air Blank	IPS*	16:35:02
Air Blank	0.000	16:35:31

\*Improper Sample

A handwritten signature in black ink, appearing to be "T. Attenberger" with a date "9/25/08" written to the right of the signature.

**Not a Successfully Completed Test Sequence**

INTOXILYZER 8000  
Location: PRESCOTT P.D  
Serial Number : 80-000337  
09/25/2008

16:35:56

---

Standard Lot#: 727501I  
Last Changed By: T. ATTENBERGER #252

---

QAS: T. ATTENBERGER #252  
PRESCOTT PD

---

Operator: T. ATTENBERGER #252  
PRESCOTT PD

---

Subject: SQAP,TEST  
DOB: 01/01/1900  
Sex: M Weight: 300


---

15 Minute Deprivation Period? Yes

---

Test	g/210L	Time
Air Blank	RFI*	16:36:39
Air Blank	0.000	16:37:08

\*RFI Detect



**Not a Successfully Completed Test Sequence**

INTOXILYZER 8000  
Location: PRESCOTT P.D  
Serial Number : 80-000337  
09/25/2008

16:37:33

---

Standard Lot#: 727501I  
Last Changed By: T. ATTENBERGER #252

---

QAS: T. ATTENBERGER #252  
PRESCOTT PD

---

Operator: T. ATTENBERGER #252  
PRESCOTT PD

---

Subject: SQAP, TEST  
DOB: 01/01/1900  
Sex: M Weight: 300

---

15 Minute Deprivation Period? Yes

---

Test	g/210L	Time
Air Blank	0.000	16:38:37
Diagnostic Test	Pass	16:39:04
Air Blank	0.000	16:39:32
0.100 Cal Check	0.096	16:39:50
Air Blank	0.000	16:40:20
Subject Test	0.023*	16:43:28
Air Blank	0.000	16:43:58
Air Blank	0.000	16:44:25
Subject Test	XXX**	16:44:57
Air Blank	0.000	16:45:31

\*Deficient Sample

\*\*Invalid Sample

Begin new deprivation period



**Not a Successfully Completed Test Sequence**

EXHIBIT G-2  
THIS REPORT PREPARED PURSUANT TO DUTY IMPOSED BY A.A.C. R13-10-104 (A)

ARIZONA DEPARTMENT OF PUBLIC SAFETY  
STANDARD QUALITY ASSURANCE PROCEDURES  
INTOXILYZER MODEL 8000  
STANDARD CALIBRATION CHECK PROCEDURE

QA SPECIALIST T. ATENBERGER #202 AGENCY Prescott P.D.  
DATE 092508 TIME 1647  
INTOXILYZER SERIAL # 80-000337 LOCATION P.P.D.

- () 1. Ensure that gas tank is attached to instrument and contains a standard alcohol concentration solution 100 AC.  
OR  
Pour a standard alcohol concentration solution \_\_\_\_\_ AC, into a clean dry simulator and assemble the simulator. Ensure that a tight seal has been made. Turn on the simulator and allow temperature to reach 34°C ± 0.2°C
- () 2. Intoxilyzer 8000 display reads "PUSH BUTTON TO START"
- () 3. Go to the "Control Testing Menu". Select "D" for dry control test or "W" for wet control test. After selection is made press ENTER.
- () 4. Air blank completed.
- () 5. Calibration check completed. Test results 0. 097 AC.
- () 6. Air blank completed.
- () 7. Remove printed record. Attach the record to the completed checklist.

SIGNATURE  #202

DPS Form Exh G-2 (Rev 05-01)


Prescott PD CHANGE 5-12-06

INTOXILYZER 8000  
Location: PRESCOTT P.D  
Serial Number: 80-000337  
09/25/2008 16:47:47

-----  
Standard Lot#: 727501I  
Last Changed By: T. ATTENBERGER #252

-----  
QAS: T. ATTENBERGER #252  
PRESCOTT PD  
-----

Test	g/210L	Time
Air Blank	0.000	16:48:14
Cal Check	0.097	16:48:33
Air Blank	0.000	16:49:02

  
-----  
Operator's Signature