

EXHIBIT G-2
THIS REPORT PREPARED PURSUANT TO DUTY IMPOSED BY A.A.C. R13-10-104 (A)

ARIZONA DEPARTMENT OF PUBLIC SAFETY
STANDARD QUALITY ASSURANCE PROCEDURES
INTOXILYZER MODEL 8000


STANDARD CALIBRATION CHECK PROCEDURE

QA SPECIALIST G. KOESTER AGENCY PRESCOTT PD

DATE 2-18-10 TIME 1607

INTOXILYZER SERIAL # 80-000337 LOCATION PRESCOTT PD

- (X) 1. Ensure that gas tank is attached to instrument and contains a standard alcohol concentration solution .100 AC.
OR
Pour a standard alcohol concentration solution _____ AC, into a clean dry simulator and assemble the simulator. Ensure that a tight seal has been made. Turn on the simulator and allow temperature to reach $34^{\circ}\text{C} \pm 0.2^{\circ}\text{C}$
- (X) 2. Intoxilyzer 8000 display reads "PUSH BUTTON TO START"
- (X) 3. Go to the "Control Testing Menu". Select "D" for dry control test or "W" for wet control test. After selection is made press ENTER.
- (X) 4. Air blank completed.
- (X) 5. Calibration check completed. Test results 0. .099 AC.
- (X) 6. Air blank completed.
- (X) 7. Remove printed record. Attach the record to the completed checklist.

SIGNATURE G. KOESTER 

INTOXILYZER 8000

Location: PRESCOTT P.D

Serial Number: 80-000337

Core Version: 8105.48

02/18/2010

16:06:29

Standard Lot#: 818201I

Last Changed By: G. KOESTER #212

QAS: G. KOESTER #212

PRESCOTT PD

Operator: G. KOESTER #212

PRESCOTT PD

Test g/210L Time

Air Blank 0.000 16:07:00
0.100 Cal Check 0.099 16:07:21
Air Blank 0.000 16:07:51

Preventative Maintenance
Performed -- Timer Reset

EXHIBIT G-4
THIS REPORT PREPARED PURSUANT TO DUTY IMPOSED BY A.A.C. R13-10-104 (A)

ARIZONA DEPARTMENT OF PUBLIC SAFETY

STANDARD QUALITY ASSURANCE PROCEDURES
INTOXILYZER MODEL 8000

STANDARD QUALITY ASSURANCE PROCEDURE

QA SPECIALIST G. KOESTER AGENCY PRESCOTT PD

DATE 2-18-10 TIME 1610

INTOXILYZER SERIAL # 80-000337 LOCATION PRESCOTT PD

(GK) 1. Display Reads "PUSH BUTTON TO START"

DIAGNOSTIC TESTS

(GK) 1. Clock time check.

(GK) 2. Date check.

OPERATIONAL TESTS

(GK) 1. Alcohol-free subject test result 0. 0.000 AC.

(GK) 2. Error recognition logic system functioning.
Not a Successfully Completed Test Sequence printed

(GK) 3. Proper sample recognition system.
Not a Successfully Completed Test Sequence printed
Deficient sample printed.

(GK) 4. Standard Calibration Check standard 0. .100 AC. Result 0. .099 AC.

Instrument is operating properly and accurately. YES X NO

COMMENTS

QAS COMPLETED AFTER ADMINISTRATIVE SOFTWARE UPDATE.

INSTRUMENT BEING PLACED BACK INTO SERVICE

SIGNATURE G. KOESTER

INTOXILYZER 8000
Location: PRESCOTT P.D
Serial Number: 80-000337
Core Version: 8105.48
02/18/2010

16:09:55

Standard Lot#: 818201I
Last Changed By: G. KOESTER #212

QAS: G. KOESTER #212
PRESCOTT PD

Operator: G. KOESTER #212
PRESCOTT PD

Subject: SQAP, TEST
DOB: 01/01/1901
Sex: M Weight: 200

15 Minute Deprivation Period? Yes

Test	g/210L	Time
Air Blank	0.000	16:10:43
Diagnostic Test	Pass	16:11:17
Air Blank	0.000	16:11:45
0.100 Cal Check	0.099	16:12:05
Air Blank	0.000	16:12:35
Subject Test	0.000	16:12:52
Air Blank	ABT*	16:12:58
Air Blank	0.000	16:13:27

*Sequence Aborted

Not a Successfully Completed Test Sequence

INTOXILYZER 8000
Location: PRESCOTT P.D
Serial Number: 80-000337
Core Version: 8105.48
02/18/2010

16:13:30

Standard Lot#: 818201I
Last Changed By: G. KOESTER #212

QAS: G. KOESTER #212
PRESCOTT PD

Operator: G. KOESTER #212
PRESCOTT PD

Subject: SQAP, TEST
DOB: 01/01/1901
Sex: M Weight: 200

15 Minute Deprivation Period? Yes

Test	g/210L	Time
Air Blank	IPS*	16:13:59
Air Blank	0.000	16:14:27

*Improper Sample

Not a Successfully Completed Test Sequence

INTOXILYZER 8000
Location: PRESCOTT P.D
Serial Number: 80-000337
Core Version: 8105.48
02/18/2010

16:14:31

Standard Lot#: 818201I
Last Changed By: G. KOESTER #212

QAS: G. KOESTER #212
PRESCOTT PD

Operator: G. KOESTER #212
PRESCOTT PD

Subject: SQAP, TEST
DOB: 01/01/1901
Sex: M Weight: 200

15 Minute Deprivation Period? Yes

Test	g/210L	Time
Air Blank	RFI*	16:15:01
Air Blank	0.000	16:15:30

*RFI Detect

Not a Successfully Completed Test Sequence

INTOXILYZER 8000
Location: PRESCOTT P.D
Serial Number: 80-000337
Core Version: 8105.48
02/18/2010

16:15:34

Standard Lot#: 818201I
Last Changed By: G. KOESTER #212

QAS: G. KOESTER #212
PRESCOTT PD

Operator: G. KOESTER #212
PRESCOTT PD

Subject: SQAP,TEST
DOB: 01/01/1901
Sex: M Weight: 200

15 Minute Deprivation Period? Yes

Test	g/210L	Time
Air Blank	0.000	16:16:27
Diagnostic Test	Pass	16:17:01
Air Blank	0.000	16:17:29
0.100 Cal Check	0.100	16:17:50
Air Blank	0.000	16:18:19
Subject Test	0.097*	16:21:27
Air Blank	0.000	16:21:59
Air Blank	0.000	16:22:26
Subject Test	XXX**	16:23:04
Air Blank	0.000	16:23:43

*Deficient Sample

**Invalid Sample

Begin new deprivation period

Not a Successfully Completed Test Sequence

EXHIBIT G-2
THIS REPORT PREPARED PURSUANT TO DUTY IMPOSED BY A.A.C. R13-10-104 (A)

ARIZONA DEPARTMENT OF PUBLIC SAFETY
STANDARD QUALITY ASSURANCE PROCEDURES
INTOXILYZER MODEL 8000

STANDARD CALIBRATION CHECK PROCEDURE

QA SPECIALIST G. KOESTER AGENCY PRESOTT PD

DATE 2-18-10 TIME 1624

INTOXILYZER SERIAL # 80-020337 LOCATION PRESOTT PD

- (GK) 1. Ensure that gas tank is attached to instrument and contains a standard alcohol concentration solution .100 AC.
OR
Pour a standard alcohol concentration solution _____ AC, into a clean dry simulator and assemble the simulator. Ensure that a tight seal has been made. Turn on the simulator and allow temperature to reach $34^{\circ}\text{C} \pm 0.2^{\circ}\text{C}$
- (GK) 2. Intoxilyzer 8000 display reads "PUSH BUTTON TO START"
- (GK) 3. Go to the "Control Testing Menu". Select "D" for dry control test or "W" for wet control test. After selection is made press ENTER.
- (GK) 4. Air blank completed.
- (GK) 5. Calibration check completed. Test results 0. .099 AC.
- (GK) 6. Air blank completed.
- (GK) 7. Remove printed record. Attach the record to the completed checklist.

SIGNATURE G. KOESTER 

INTOXILYZER 8000

Location: PRESCOTT P.D

Serial Number: 80-000337

Core Version: 8105.48

02/18/2010

16:24:03

Standard Lot#: 818201I

Last Changed By: G. KOESTER #212

QAS: G. KOESTER #212

PRESCOTT PD

Operator: G. KOESTER #212

PRESCOTT PD

Test g/210L Time

Air Blank 0.000 16:24:33
0.100 Cal Check 0.099 16:24:54
Air Blank 0.000 16:25:23

Preventative Maintenance
Performed -- Timer Reset

EXHIBIT G-2
THIS REPORT PREPARED PURSUANT TO DUTY IMPOSED BY A.A.C. R13-10-104 (A)

ARIZONA DEPARTMENT OF PUBLIC SAFETY
STANDARD QUALITY ASSURANCE PROCEDURES
INTOXILYZER MODEL 8000

STANDARD CALIBRATION CHECK PROCEDURE

QA SPECIALIST G. ROESTER AGENCY PRESCOTT PD

DATE 2-18-10 TIME 1255

INTOXILYZER SERIAL # 80-000337 LOCATION FLAGSTAFF PD

- (GK) 1. Ensure that gas tank is attached to instrument and contains a standard alcohol concentration solution .100 AC.
OR
Pour a standard alcohol concentration solution _____ AC, into a clean dry simulator and assemble the simulator. Ensure that a tight seal has been made. Turn on the simulator and allow temperature to reach $34^{\circ}\text{C} \pm 0.2^{\circ}\text{C}$
- (GK) 2. Intoxilyzer 8000 display reads "PUSH BUTTON TO START"
- (GK) 3. Go to the "Control Testing Menu". Select "D" for dry control test or "W" for wet control test. After selection is made press ENTER.
- (GK) 4. Air blank completed.
- (GK) 5. Calibration check completed. Test results 0. .099 AC.
- (GK) 6. Air blank completed.
- (GK) 7. Remove printed record. Attach the record to the completed checklist.

SIGNATURE G. ROESTER  2/18

INTOXILYZER 8000
Location: PRESCOTT P.D
Serial Number: 80-000337
Core Version: 8105.48
02/18/2010

12:55:05

Standard Lot#: 818201I
Last Changed By: G. KOESTER #212

QAS: G. KOESTER #212
PRESCOTT PD

Operator: G. KOESTER #212
PRESCOTT PD

Test	g/210L	Time
Air Blank	0.000	12:55:38
0.100 Cal Check	0.099	12:55:59
Air Blank	0.000	12:56:29

Preventative Maintenance
Performed -- Timer Reset

EXHIBIT G-2
THIS REPORT PREPARED PURSUANT TO DUTY IMPOSED BY A.A.C. R13-10-104 (A)

ARIZONA DEPARTMENT OF PUBLIC SAFETY
STANDARD QUALITY ASSURANCE PROCEDURES
INTOXILYZER MODEL 8000

STANDARD CALIBRATION CHECK PROCEDURE

QA SPECIALIST G. KOESTER AGENCY PRESTON PD

DATE 2-18-16 TIME 0915

INTOXILYZER SERIAL # 80-000337 LOCATION PRESTON PD

- (GK) 1. Ensure that gas tank is attached to instrument and contains a standard alcohol concentration solution 0.100 AC.
OR
Pour a standard alcohol concentration solution _____ AC, into a clean dry simulator and assemble the simulator. Ensure that a tight seal has been made. Turn on the simulator and allow temperature to reach $34^{\circ}\text{C} \pm 0.2^{\circ}\text{C}$
- (GK) 2. Intoxilyzer 8000 display reads "PUSH BUTTON TO START"
- (GK) 3. Go to the "Control Testing Menu". Select "D" for dry control test or "W" for wet control test. After selection is made press ENTER.
- (GK) 4. Air blank completed.
- (GK) 5. Calibration check completed. Test results 0. 0.100 AC.
- (GK) 6. Air blank completed.
- (GK) 7. Remove printed record. Attach the record to the completed checklist.

SIGNATURE G. KOESTER

EXHIBIT G-4
THIS REPORT PREPARED PURSUANT TO DUTY IMPOSED BY A.A.C. R13-10-104 (A)

ARIZONA DEPARTMENT OF PUBLIC SAFETY
STANDARD QUALITY ASSURANCE PROCEDURES
INTOXILYZER MODEL 8000

STANDARD QUALITY ASSURANCE PROCEDURE

QA SPECIALIST G. KOESTER AGENCY PRESOTT PD

DATE 2-18-10 TIME 0917

INTOXILYZER SERIAL # 80-000337 LOCATION PRESOTT PD

1. Display Reads "PUSH BUTTON TO START"

DIAGNOSTIC TESTS

1. Clock time check.

2. Date check.

OPERATIONAL TESTS

1. Alcohol-free subject test result 0. 0.000 AC.

2. Error recognition logic system functioning.
Not a Successfully Completed Test Sequence printed

3. Proper sample recognition system.
Not a Successfully Completed Test Sequence printed
Deficient sample printed.

4. Standard Calibration Check standard 0. 1.00 AC. Result 0. 0.99 AC.

Instrument is operating properly and accurately. YES NO


COMMENTS

QAS COMPLETED PRIOR TO INSTRUMENT BEING TAKEN OUT OF

SERVICE FOR ADMINISTRATIVE SOFTWARE UPDATE.

ALL SYSTEMS FUNCTIONING PROPERLY. PROPER SAMPLE TEST EXCEEDED

RANGE DUE TO TOO MUCH MOUTH ALCOHOL.

SIGNATURE G. KOESTER 

Prescott PD 02-18-10

INTOXILYZER 8000
Location: PRESCOTT P.D
Serial Number : 80-000337
02/18/2010

09:16:50

Standard Lot#: 818201I
Last Changed By: G. KOESTER #212

QAS: G. KOESTER #212
PRESCOTT PD

Operator: G. KOESTER #212
PRESCOTT PD

Subject: SQAP,TEST
DOB: 01/01/1901
Sex: M Weight: 200

15 Minute Deprivation Period? Yes

Test	g/210L	Time
Air Blank	0.000	09:17:38
Diagnostic Test	Pass	09:18:06
Air Blank	0.000	09:18:33
0.100 Cal Check	0.100	09:18:52
Air Blank	0.000	09:19:22
Subject Test	0.000	09:19:39
Air Blank	0.000	09:20:06
Five Minute Wait	ABT*	
Air Blank	0.000	09:20:41

*Sequence Aborted

Not a Successfully Completed Test Sequence

INTOXILYZER 8000
Location: PRESCOTT P.D
Serial Number : 80-000337
02/18/2010

09:21:28

Standard Lot#: 818201I
Last Changed By: G. KOESTER #212

QAS: G. KOESTER #212
PRESCOTT PD

Operator: G. KOESTER #212
PRESCOTT PD

Subject: SQAP, TEST
DOB: 01/01/1901
Sex: M Weight: 200

15 Minute Deprivation Period? Yes

Test	g/210L	Time
Air Blank	IPS*	09:21:55
Air Blank	0.000	09:22:24

*Improper Sample

Not a Successfully Completed Test Sequence

INTOXILYZER 8000
Location: PRESCOTT P.D
Serial Number : 80-000337
02/18/2010

09:22:30

Standard Lot#: 818201I
Last Changed By: G. KOESTER #212

QAS: G. KOESTER #212
PRESCOTT PD

Operator: G. KOESTER #212
PRESCOTT PD

Subject: SQAP,TEST
DOB: 01/01/1901
Sex: M Weight: 200

15 Minute Deprivation Period? Yes

Test	g/210L	Time
Air Blank	RFI*	09:23:13
Air Blank	0.000	09:23:41

*RFI Detect

Not a Successfully Completed Test Sequence

INTOXILYZER 8000
Location: PRESCOTT P.D
Serial Number : 80-000337
02/18/2010

09:24:09

Standard Lot#: 818201I
Last Changed By: G. KOESTER #212

QAS: G. KOESTER #212
PRESCOTT PD

Operator: G. KOESTER #212
PRESCOTT PD

Subject: SQAP,TEST
DOB: 01/01/1901
Sex: M Weight: 200

15 Minute Deprivation Period? Yes

Test	g/210L	Time
Air Blank	0.000	09:24:59
Diagnostic Test	Pass	09:25:27
Air Blank	0.000	09:25:54
0.100 Cal Check	0.100	09:26:13
Air Blank	0.000	09:26:42
Subject Test	0.097*	09:29:51
Air Blank	0.000	09:30:22
Air Blank	0.000	09:30:50
Subject Test	RGE**	09:31:08
Air Blank	0.000	09:31:48

*Deficient Sample

**Range Exceeded

Not a Successfully Completed Test Sequence

INTOXILYZER 8000
Location: PRESCOTT P.D
Serial Number : 80-000337
02/18/2010

09:31:55

Standard Lot#: 818201I
Last Changed By: G. KOESTER #212

QAS: G. KOESTER #212
PRESCOTT PD

Operator: G. KOESTER #212
PRESCOTT PD

Subject: SQAP, TEST
DOB: 01/01/1901
Sex: M Weight: 200

15 Minute Deprivation Period? Yes

Test	g/210L	Time
Air Blank	0.000	09:32:40
Diagnostic Test	Pass	09:33:07
Air Blank	0.000	09:33:35
0.100 Cal Check	0.099	09:33:53
Air Blank	0.000	09:34:23
Subject Test	RGE*	09:34:41
Air Blank	0.000	09:35:16

*Range Exceeded

Not a Successfully Completed Test Sequence

INTOXILYZER 8000
Location: PRESCOTT P.D
Serial Number : 80-000337
02/18/2010

09:35:22

Standard Lot#: 818201I
Last Changed By: G. KOESTER #212

QAS: G. KOESTER #212
PRESCOTT PD

Operator: G. KOESTER #212
PRESCOTT PD

Subject: SQAP,TEST
DOB: 01/01/1901
Sex: M Weight: 200

15 Minute Deprivation Period? Yes

Test	g/210L	Time
Air Blank	0.000	09:36:09
Diagnostic Test	Pass	09:36:36
Air Blank	0.000	09:37:03
0.100 Cal Check	0.099	09:37:22
Air Blank	0.000	09:37:52
Subject Test	XXX*	09:38:13
Air Blank	0.000	09:38:45

*Invalid Sample
Begin new deprivation period

Not a Successfully Completed Test Sequence

EXHIBIT G-2
THIS REPORT PREPARED PURSUANT TO DUTY IMPOSED BY A.A.C. R13-10-104 (A)

ARIZONA DEPARTMENT OF PUBLIC SAFETY

STANDARD QUALITY ASSURANCE PROCEDURES
INTOXILYZER MODEL 8000

STANDARD CALIBRATION CHECK PROCEDURE

QA SPECIALIST G. KOESTER AGENCY PRESOTT PD

DATE 2-18-10 TIME 0939

INTOXILYZER SERIAL # 80-080337 LOCATION PRESOTT PD

(K) 1. Ensure that gas tank is attached to instrument and contains a standard alcohol concentration solution .100 AC.
OR

Pour a standard alcohol concentration solution _____ AC, into a clean dry simulator and assemble the simulator. Ensure that a tight seal has been made. Turn on the simulator and allow temperature to reach $34^{\circ}\text{C} \pm 0.2^{\circ}\text{C}$

(K) 2. Intoxilyzer 8000 display reads "PUSH BUTTON TO START"

(GK) 3. Go to the "Control Testing Menu". Select "D" for dry control test or "W" for wet control test. After selection is made press ENTER.

(GK) 4. Air blank completed.

(GK) 5. Calibration check completed. Test results 0.099 AC.

(GK) 6. Air blank completed.

(GK) 7. Remove printed record. Attach the record to the completed checklist.

SIGNATURE G. KOESTER

INTOXILYZER 8000

Location: PRESCOTT P.D

Serial Number: 80-000337

02/18/2010

09:39:09

Standard Lot#: 818201I

Last Changed By: G. KOESTER #212

QAS: G. KOESTER #212

PRESCOTT PD

Test	g/210L	Time
Air Blank	0.000	09:39:36
Cal Check	0.099	09:39:55
Air Blank	0.000	09:40:24



Operator's Signature