

EXHIBIT G-2
THIS REPORT PREPARED PURSUANT TO DUTY IMPOSED BY A.A.C. R13-10-104 (A)

ARIZONA DEPARTMENT OF PUBLIC SAFETY
STANDARD QUALITY ASSURANCE PROCEDURES
INTOXILYZER MODEL 8000

STANDARD CALIBRATION CHECK PROCEDURE

QA SPECIALIST M. MEDINA AGENCY PRESCOTT P.D.
DATE 2/9/12 TIME 1920
INTOXILYZER SERIAL # 80-000337 LOCATION PRESCOTT P.D.

- () 1. Ensure that gas tank is attached to instrument and contains a standard alcohol concentration solution .100 AC.
OR
Pour a standard alcohol concentration solution _____ AC, into a clean dry simulator and assemble the simulator. Ensure that a tight seal has been made. Turn on the simulator and allow temperature to reach $34^{\circ}\text{C} \pm 0.2^{\circ}\text{C}$
- () 2. Intoxilyzer 8000 display reads "PUSH BUTTON TO START"
- () 3. Go to the "Control Testing Menu". Select "D" for dry control test or "W" for wet control test. After selection is made press ENTER.
- () 4. Air blank completed.
- () 5. Calibration check completed. Test results 0. 101 AC.
- () 6. Air blank completed.
- () 7. Remove printed record. Attach the record to the completed checklist.

SIGNATURE 

M. MEDINA 356
2/9/12

INTOXILYZER 8000
Location: PRESCOTT P.D
Serial Number: 80-000337
Core Version: 8105.48
02/09/2012 19:20:51

Standard Lot#: 03811100A2
Last Changed By: M. MEDINA #356

QAS: M. MEDINA #356
PRESCOTT PD

Operator: M. MEDINA #356
PRESCOTT PD

Test	g/210L	Time
Air Blank	0.000	19:21:32
0.100 Cal Check	0.101	19:21:53
Air Blank	0.000	19:22:23

Preventative Maintenance
Performed -- Timer Reset

INTOXILYZER 8000
Location: PRESCOTT P.D
Serial Number: 80-000337
Core Version: 8105.48
02/09/2012

19:26:42

M. MEDINA
356
2/9/12

Standard Lot#: 17911100A1
Last Changed By: M. MEDINA #356

QAS: M. MEDINA #356
PRESCOTT PD

Operator: M. MEDINA #356
PRESCOTT PD

Test	g/210L	Time
Air Blank	0.000	19:27:14
0.100 Cal Check	0.100	19:27:35
Air Blank	0.000	19:28:05

Preventative Maintenance
Performed -- Timer Reset

EXHIBIT G-4
THIS REPORT PREPARED PURSUANT TO DUTY IMPOSED BY A.A.C. R13-10-104 (A)

ARIZONA DEPARTMENT OF PUBLIC SAFETY
STANDARD QUALITY ASSURANCE PROCEDURES
INTOXILYZER MODEL 8000

STANDARD QUALITY ASSURANCE PROCEDURE

QA SPECIALIST M. MEDINA AGENCY PRESOTT
DATE 2/9/12 TIME 1948
INTOXILYZER SERIAL # 80-000337 LOCATION PRESOTT P.D.

() 1. Display Reads "PUSH BUTTON TO START"

DIAGNOSTIC TESTS

- () 1. Clock time check.
- () 2. Date check.

OPERATIONAL TESTS

- () 1. Alcohol-free subject test result 0. 000 AC.
- () 2. Error recognition logic system functioning.
Not a Successfully Completed Test Sequence printed
- () 3. Proper sample recognition system.
Not a Successfully Completed Test Sequence printed
Deficient sample printed.
- () 4. Standard Calibration Check standard 0. 100 AC. Result 0. 099 AC.

Instrument is operating properly and accurately. YES NO

COMMENTS

SIGNATURE 

M. MEDINA
2/9/12

INTOXILYZER 8000
Location: PRESCOTT P.D
Serial Number: 80-000337
Core Version: 8105.48
02/09/2012

19:45:54

Standard Lot#: 17911100A1
Last Changed By: M. MEDINA #356

QAS: M. MEDINA #356
PRESCOTT PD

Operator: M. MEDINA #356
PRESCOTT PD

Subject: SQAP, TEST
DOB: 02/02/2000
Sex: M Weight: 180

15 Minute Deprivation Period? Yes

Test	g/210L	Time
Air Blank	0.000	19:47:15
Diagnostic Test	Pass	19:47:49
Air Blank	0.000	19:48:17
0.100 Cal Check	0.099	19:48:38
Air Blank	0.000	19:49:07
Subject Test	0.000	19:49:30
Air Blank	0.000	19:49:57
Wait	ABT*	
Air Blank	RFI**	19:50:20

*Sequence Aborted
**RFI Detect

Not a Successfully Completed Test Sequence

M. MEDINA 356
2/9/12

INTOXILYZER 8000
Location: PRESCOTT P.D
Serial Number: 80-000337
Core Version: 8105.48
02/09/2012

19:51:22

Standard Lot#: 17911100A1
Last Changed By: M. MEDINA #356

QAS: M. MEDINA #356
PRESCOTT PD

Operator: M. MEDINA #356
PRESCOTT PD

Subject: SQAP, TEST
DOB: 02/02/2000
Sex: M Weight: 180

15 Minute Deprivation Period? Yes

Test	g/210L	Time
Air Blank	0.000	19:52:26
Diagnostic Test	Pass	19:53:00
Air Blank	0.000	19:53:28
0.100 Cal Check	0.099	19:53:49
Air Blank	0.000	19:54:18
Subject Test	0.063*	19:57:26
Air Blank	0.000	19:57:58
Air Blank	0.000	19:58:25
Subject Test	XXX**	19:58:49
Air Blank	IPS***	19:59:03

*Deficient Sample

**Invalid Sample

Begin new deprivation period

***Improper Sample

Not a Successfully Completed Test Sequence

M. MEDINA
356
2/9/12

INTOXILYZER 8000
Location: PRESCOTT P.D
Serial Number: 80-000337
Core Version: 8105.48
02/09/2012 20:04:38

Standard Lot#: 17911100A1
Last Changed By: M. MEDINA #356

QAS: M. MEDINA #356
PRESCOTT PD

Operator: M. MEDINA #356
PRESCOTT PD

Test	g/210L	Time
Air Blank	0.000	20:05:09
0.100 Cal Check	0.099	20:05:30
Air Blank	0.000	20:06:00

Preventative Maintenance
Performed -- Timer Reset

EXHIBIT G-2
THIS REPORT PREPARED PURSUANT TO DUTY IMPOSED BY A.A.C. R13-10-104 (A)

ARIZONA DEPARTMENT OF PUBLIC SAFETY
STANDARD QUALITY ASSURANCE PROCEDURES
INTOXILYZER MODEL 8000

STANDARD CALIBRATION CHECK PROCEDURE

QA SPECIALIST M. MEDINA AGENCY PRESCOTT P.D.
DATE 2/9/12 TIME 2000
INTOXILYZER SERIAL # 80-000337 LOCATION PRESCOTT P.D.

- () 1. Ensure that gas tank is attached to instrument and contains a standard alcohol concentration solution 0.100 AC.
OR
Pour a standard alcohol concentration solution _____ AC, into a clean dry simulator and assemble the simulator. Ensure that a tight seal has been made. Turn on the simulator and allow temperature to reach $34^{\circ}\text{C} \pm 0.2^{\circ}\text{C}$
- () 2. Intoxilyzer 8000 display reads "PUSH BUTTON TO START"
- () 3. Go to the "Control Testing Menu". Select "D" for dry control test or "W" for wet control test. After selection is made press ENTER.
- () 4. Air blank completed.
- () 5. Calibration check completed. Test results 0. 099 AC.
- () 6. Air blank completed.
- () 7. Remove printed record. Attach the record to the completed checklist.

SIGNATURE 

DPS Form Exh G-2 (Rev 05-01)

1ST AIR BLANK "AMBIENT FAIL"

M. MEDINA
356
2/9/12

INTOXILYZER 8000
Location: PRESCOTT P.D
Serial Number: 80-000337
Core Version: 8105.48
02/09/2012 20:03:23

Standard Lot#: 17911100A1
Last Changed By: M. MEDINA #356

QAS: M. MEDINA #356
PRESCOTT PD

Operator: M. MEDINA #356
PRESCOTT PD

Test	g/210L	Time
Air Blank	AMB*	20:03:55
Air Blank	0.000	20:04:23

*Ambient Fail

PRESCOTT PD
2/9/12