

EXHIBIT G-2  
THIS REPORT PREPARED PURSUANT TO DUTY IMPOSED BY A.A.C. R13-10-104 (A)

ARIZONA DEPARTMENT OF PUBLIC SAFETY  
STANDARD QUALITY ASSURANCE PROCEDURES  
INTOXILYZER MODEL 8000

STANDARD CALIBRATION CHECK PROCEDURE

QA SPECIALIST D. [Signature] AGENCY Puerco P.D.

DATE 05/31/12 TIME 0749

INTOXILYZER SERIAL # 80-000337 LOCATION P.P.D.

- () 1. Ensure that gas tank is attached to instrument and contains a standard alcohol concentration solution .100 AC.  
OR  
Pour a standard alcohol concentration solution \_\_\_\_\_ AC, into a clean dry simulator and assemble the simulator. Ensure that a tight seal has been made. Turn on the simulator and allow temperature to reach  $34^{\circ}\text{C} \pm 0.2^{\circ}\text{C}$
- () 2. Intoxilyzer 8000 display reads "PUSH BUTTON TO START"
- () 3. Go to the "Control Testing Menu". Select "D" for dry control test or "W" for wet control test. After selection is made press ENTER.
- () 4. Air blank completed.
- () 5. Calibration check completed. Test results 0.099 AC.
- () 6. Air blank completed.
- () 7. Remove printed record. Attach the record to the completed checklist.

SIGNATURE [Signature]

INTOXILYZER 8000

Location: PRESCOTT P.D

Serial Number: 80-000337

Core Version: 8105.48

05/31/2012

07:49:50

-----  
Standard Lot#: 17911100A1

Last Changed By: M. MEDINA #356  
-----

QAS: T. ATTENBERGER #252

PRESCOTT PD  
-----

Operator: T. ATTENBERGER #252

PRESCOTT PD  
-----

Test	g/210L	Time
Air Blank	0.000	07:50:22
0.100 Cal Check	0.099	07:50:43
Air Blank	0.000	07:51:12

Preventative Maintenance  
Performed -- Timer Reset



EXHIBIT G-4  
THIS REPORT PREPARED PURSUANT TO DUTY IMPOSED BY A.A.C. R13-10-104 (A)

ARIZONA DEPARTMENT OF PUBLIC SAFETY

STANDARD QUALITY ASSURANCE PROCEDURES  
INTOXILYZER MODEL 8000

STANDARD QUALITY ASSURANCE PROCEDURE

QA SPECIALIST D. Agnew AGENCY Prescott PD  
DATE 05/31/12 TIME 0753  
INTOXILYZER SERIAL # 80-005337 LOCATION P.P.D.

1. Display Reads "PUSH BUTTON TO START"

DIAGNOSTIC TESTS

1. Clock time check.  
 2. Date check.

OPERATIONAL TESTS

1. Alcohol-free subject test result 0. 000 AC.  
 2. Error recognition logic system functioning.  
Not a Successfully Completed Test Sequence printed  
 3. Proper sample recognition system.  
Not a Successfully Completed Test Sequence printed  
Deficient sample printed.  
 4. Standard Calibration Check standard 0. 100 AC. Result 0. 099 AC.

Instrument is operating properly and accurately. YES  NO

COMMENTS

ALL FUNCTIONING CORRECTLY.

SIGNATURE D. Agnew

INTOXILYZER 8000

Location: PRESCOTT P.D

Serial Number: 80-000337

Core Version: 8105.48

05/31/2012

07:53:12

---

Standard Lot#: 17911100A1

Last Changed By: M. MEDINA

#356

---

QAS: T. ATTENBERGER

#252

PRESCOTT PD

---

Operator: T. ATTENBERGER

#252

PRESCOTT PD

---

Subject: SQAP, TEST

DOB: 01/01/1900

Sex: M

Weight: 300

---

15 Minute Deprivation Period? Yes

---

Test	g/210L	Time
Air Blank	0.000	07:54:34
Diagnostic Test	Pass	07:55:08
Air Blank	0.000	07:55:36
0.100 Cal Check	0.099	07:55:56
Air Blank	0.000	07:56:26
Subject Test	0.000	07:56:48
Air Blank	0.000	07:57:16
Wait	ABT*	
Air Blank	0.000	07:57:51

\*Sequence Aborted



Not a Successfully Completed Test Sequence

INTOXILYZER 8000

Location: PRESCOTT P.D

Serial Number: 80-000337

Core Version: 8105.48

05/31/2012

07:58:39

---

Standard Lot#: 17911100A1

Last Changed By: M. MEDINA

#356

---

QAS: T. ATTENBERGER

#252

PRESCOTT PD

---

Operator: T. ATTENBERGER

#252

PRESCOTT PD

---

Subject: SQAP, TEST

DOB: 01/01/1900

Sex: M

Weight: 300

---

15 Minute Deprivation Period? Yes

---

Test	g/210L	Time
Air Blank	IPS*	07:59:25
Air Blank	0.000	07:59:53

\*Improper Sample

A handwritten signature in black ink, appearing to be "T. Attenberger" with a date "05/31/2012" written next to it.

Not a Successfully Completed Test Sequence

INTOXILYZER 8000  
Location: PRESCOTT P.D  
Serial Number: 80-000337  
Core Version: 8105.48  
05/31/2012

08:00:18

---

Standard Lot#: 17911100A1  
Last Changed By: M. MEDINA #356

---

QAS: T. ATTENBERGER #252  
PRESCOTT PD

---

Operator: T. ATTENBERGER #252  
PRESCOTT PD

---

Subject: SQAP, TEST  
DOB: 01/01/1900  
Sex: M Weight: 300

---

15 Minute Deprivation Period? Yes

---

Test	g/210L	Time
Air Blank	RFI*	08:01:12
Air Blank	0.000	08:01:41

\*RFI Detect



Not a Successfully Completed Test Sequence

INTOXILYZER 8000  
Location: PRESCOTT P.D  
Serial Number: 80-000337  
Core Version: 8105.48  
05/31/2012

08:02:13

---

Standard Lot#: 17911100A1  
Last Changed By: M. MEDINA #356

---

QAS: T. ATTENBERGER #252  
PRESCOTT PD

---

Operator: T. ATTENBERGER #252  
PRESCOTT PD

---

Subject: SQAP, TEST  
DOB: 01/01/1900  
Sex: M Weight: 300

---

15 Minute Deprivation Period? Yes

---

Test	g/210L	Time
Air Blank	0.000	08:03:10
Diagnostic Test	Pass	08:03:44
Air Blank	0.000	08:04:12
0.100 Cal Check	0.099	08:04:32
Air Blank	0.000	08:05:02
Subject Test	0.034	08:05:37
Air Blank	0.000	08:06:10
Wait	ABT*	
Air Blank	0.000	08:06:52

\*Sequence Aborted ←

*BLEW TOO LONG ON DEPR SAMPLE*



Not a Successfully Completed Test Sequence

INTOXILYZER 8000  
Location: PRESCOTT P.D  
Serial Number: 80-000337  
Core Version: 8105.48  
05/31/2012

08:07:39

---

Standard Lot#: 17911100A1  
Last Changed By: M. MEDINA #356

---

QAS: T. ATTENBERGER #252  
PRESCOTT PD

---

Operator: T. ATTENBERGER #252  
PRESCOTT PD

---

Subject: SQAP,TEST  
DOB: 01/01/1900  
Sex: M Weight: 300

---

15 Minute Deprivation Period? Yes


---

Test	g/210L	Time
Air Blank	0.000	08:08:31
Diagnostic Test	Pass	08:09:05
Air Blank	0.000	08:09:33
0.100 Cal Check	0.099	08:09:53
Air Blank	0.000	08:10:23
Subject Test	0.040*	08:13:30
Air Blank	0.000	08:14:01
Air Blank	0.000	08:14:29
Subject Test	XXX**	08:15:11
Air Blank	0.000	08:15:44

\*Deficient Sample

\*\*Invalid Sample

Begin new deprivation period



Not a Successfully Completed Test Sequence

EXHIBIT G-2  
THIS REPORT PREPARED PURSUANT TO DUTY IMPOSED BY A.A.C. R13-10-104 (A)

ARIZONA DEPARTMENT OF PUBLIC SAFETY  
STANDARD QUALITY ASSURANCE PROCEDURES  
INTOXILYZER MODEL 8000

STANDARD CALIBRATION CHECK PROCEDURE

QA SPECIALIST T. Aram... 1234 AGENCY P. P. D.

DATE 03/31/12 TIME 0817

INTOXILYZER SERIAL # 80-000337 LOCATION P. P. D.

- () 1. Ensure that gas tank is attached to instrument and contains a standard alcohol concentration solution .100 AC.  
OR  
Pour a standard alcohol concentration solution \_\_\_\_\_ AC, into a clean dry simulator and assemble the simulator. Ensure that a tight seal has been made. Turn on the simulator and allow temperature to reach 34°C ± 0.2°C
- () 2. Intoxilyzer 8000 display reads "PUSH BUTTON TO START"
- () 3. Go to the "Control Testing Menu". Select "D" for dry control test or "W" for wet control test. After selection is made press ENTER.
- () 4. Air blank completed.
- () 5. Calibration check completed. Test results 0.099 AC.
- () 6. Air blank completed.
- () 7. Remove printed record. Attach the record to the completed checklist.

SIGNATURE [Signature]

INTOXILYZER 8000  
Location: PRESCOTT P.D  
Serial Number: 80-000337  
Core Version: 8105.48  
05/31/2012

08:17:10


-----  
Standard Lot#: 17911100A1  
Last Changed By: M. MEDINA #356  
-----

QAS: T. ATTENBERGER #252  
PRESCOTT PD  
-----

Operator: T. ATTENBERGER #252  
PRESCOTT PD  
-----

Test	g/210L	Time
Air Blank	0.000	08:17:42
0.100 Cal Check	0.099	08:18:03
Air Blank	0.000	08:18:32

Preventative Maintenance  
Performed -- Timer Reset



RECEIVED  
SERIES 1 P. YAM  
MAY 31 2012