


EXHIBIT G-2
THIS REPORT PREPARED PURSUANT TO DUTY IMPOSED BY A.A.C. R13-10-104 (A)

ARIZONA DEPARTMENT OF PUBLIC SAFETY
STANDARD QUALITY ASSURANCE PROCEDURES
INTOXILYZER MODEL 8000

STANDARD CALIBRATION CHECK PROCEDURE

QA SPECIALIST N. MEDINA AGENCY PRESCOTT P.D.
DATE 5/30/13 TIME 0929
INTOXILYZER SERIAL # 80-000337 LOCATION PRESCOTT P.D.

- () 1. Ensure that gas tank is attached to instrument and contains a standard alcohol concentration solution .100 AC.
OR
Pour a standard alcohol concentration solution _____ AC, into a clean dry simulator and assemble the simulator. Ensure that a tight seal has been made. Turn on the simulator and allow temperature to reach $34^{\circ}\text{C} \pm 0.2^{\circ}\text{C}$.
- () 2. Intoxilyzer 8000 display reads "PUSH BUTTON TO START"
- () 3. Go to the "Control Testing Menu". Select "D" for dry control test or "W" for wet control test. After selection is made press ENTER.
- () 4. Air blank completed.
- () 5. Calibration check completed. Test results 0. 097 AC.
- () 6. Air blank completed.
- () 7. Remove printed record. Attach the record to the completed checklist.

SIGNATURE 

DPS Form Exh G-2 (Rev 05-01)

M. MEDINA
5/30/13

INTOXILYZER 8000
Location: PRESCOTT P.D
Serial Number: 80-000337
Core Version: 8105.48
05/30/2013 09:29:13

Standard Lot#: 34612100A1
Last Changed By: M. MEDINA #356

QAS: M. MEDINA #356
PRESCOTT PD

Operator: M. MEDINA #356
PRESCOTT PD

Test	g/210L	Time
Air Blank	0.000	09:30:01
0.100 Cal Check	0.097	09:30:22
Air Blank	0.000	09:30:52

Preventative Maintenance
Performed -- Timer Reset

EXHIBIT G-4
THIS REPORT PREPARED PURSUANT TO DUTY IMPOSED BY A.A.C. R13-10-104 (A)

ARIZONA DEPARTMENT OF PUBLIC SAFETY
STANDARD QUALITY ASSURANCE PROCEDURES
INTOXILYZER MODEL 8000
STANDARD QUALITY ASSURANCE PROCEDURE

QA SPECIALIST M. MEDINA AGENCY PRESOTT P.D.
DATE 5/30/13 TIME 0933
INTOXILYZER SERIAL # 80-000337 LOCATION PRESOTT P.D.

1. Display Reads "PUSH BUTTON TO START"

DIAGNOSTIC TESTS

1. Clock time check.
 2. Date check.

OPERATIONAL TESTS

1. Alcohol-free subject test result 0. 000 AC.
 2. Error recognition logic system functioning.
Not a Successfully Completed Test Sequence printed
 3. Proper sample recognition system.
Not a Successfully Completed Test Sequence printed
Deficient sample printed.
 4. Standard Calibration Check standard 0. 100 AC. Result 0. 096 AC.

Instrument is operating properly and accurately. YES NO

COMMENTS

SIGNATURE [Signature]

DPS Form Exh G-4 (Rev 05-01)

PRESOTT PD CALL# 5-18-06
SD NATION #27

M. MEDINA
5/30/13

INTOXILYZER 8000
Location: PRESCOTT P.D
Serial Number: 80-000337
Core Version: 8105.48
05/30/2013

09:33:14

Standard Lot#: 34612100A1
Last Changed By: M. MEDINA #356

QAS: M. MEDINA #356
PRESCOTT PD

Operator: M. MEDINA #356
PRESCOTT PD

Subject: SQAP, TEST
DOB: 02/02/2000
Sex: M Weight: 180

15 Minute Deprivation Period? Yes

Test	g/210L	Time
Air Blank	0.000	09:34:30
Diagnostic Test	Pass	09:35:04
Air Blank	0.000	09:35:32
0.100 Cal Check	0.097	09:35:52
Air Blank	0.000	09:36:22
Subject Test	0.000	09:36:44
Air Blank	0.000	09:37:12
Wait	ABT*	
Air Blank	RFI**	09:37:30

*Sequence Aborted

**RFI Detect

Not a Successfully Completed Test Sequence

M. MEDINA
5/30/13

INTOXILYZER 8000
Location: PRESCOTT P.D
Serial Number: 80-000337
Core Version: 8105.48
05/30/2013 09:38:06

Standard Lot#: 34612100A1
Last Changed By: M. MEDINA #356

QAS: M. MEDINA #356
PRESCOTT PD

Operator: M. MEDINA #356
PRESCOTT PD

Subject: SQAP, TEST
DOB: 02/02/2000
Sex: M Weight: 180

15 Minute Deprivation Period? Yes

Test	g/210L	Time
Air Blank	IPS*	09:39:26
Air Blank	0.000	09:39:54

*Improper Sample

Not a Successfully Completed Test Sequence

M. MEDINA 354
5/30/13

INTOXILYZER 8000
Location: PRESCOTT P.D
Serial Number: 80-000337
Core Version: 8105.48
05/30/2013

09:40:22

Standard Lot#: 34612100A1
Last Changed By: M. MEDINA #356

QAS: M. MEDINA #356
PRESCOTT PD

Operator: M. MEDINA #356
PRESCOTT PD

Subject: SQAP, TEST
DOB: 02/02/2000
Sex: M Weight: 180

15 Minute Deprivation Period? Yes

Test	g/210L	Time
Air Blank	0.000	09:41:19
Diagnostic Test	Pass	09:41:53
Air Blank	0.000	09:42:20
0.100 Cal Check	0.097	09:42:41
Air Blank	0.000	09:43:11
Subject Test	0.043*	09:46:18
Air Blank	0.000	09:46:49
Air Blank	0.000	09:47:17
Subject Test	XXX**	09:47:43
Air Blank	0.000	09:48:16

*Deficient Sample

**Invalid Sample

Begin new deprivation period

Not a Successfully Completed Test Sequence

EXHIBIT G-2
THIS REPORT PREPARED PURSUANT TO DUTY IMPOSED BY A.A.C. R13-10-104 (A)

ARIZONA DEPARTMENT OF PUBLIC SAFETY
STANDARD QUALITY ASSURANCE PROCEDURES
INTOXILYZER MODEL 8000
STANDARD CALIBRATION CHECK PROCEDURE

QA SPECIALIST M. MEDINA AGENCY PRESOTT
DATE 5/30/13 TIME 0949
INTOXILYZER SERIAL # 80 000337 LOCATION PRESOTT P.D.

- () 1. Ensure that gas tank is attached to instrument and contains a standard alcohol concentration solution, 100 AC.
OR
Pour a standard alcohol concentration solution _____ AC, into a clean dry simulator and assemble the simulator. Ensure that a tight seal has been made. Turn on the simulator and allow temperature to reach $34^{\circ}\text{C} \pm 0.2^{\circ}\text{C}$
- () 2. Intoxilyzer 8000 display reads "PUSH BUTTON TO START"
- () 3. Go to the "Control Testing Menu". Select "D" for dry control test or "W" for wet control test. After selection is made press ENTER.
- () 4. Air blank completed.
- () 5. Calibration check completed. Test results 0. 096 AC.
- () 6. Air blank completed.
- () 7. Remove printed record. Attach the record to the completed checklist.

SIGNATURE 

DPS Form Exh G-2 (Rev 05-01)

PRESOTT PD CHANGE 5-18-06
7-10-000337

M. MEDINA
5/30/13

INTOXILYZER 8000
Location: PRESCOTT P.D
Serial Number: 80-000337
Core Version: 8105.48
05/30/2013 09:49:31

Standard Lot#: 34612100A1
Last Changed By: M. MEDINA #356

QAS: M. MEDINA #356
PRESCOTT PD

Operator: M. MEDINA #356
PRESCOTT PD

Test	g/210L	Time
Air Blank	0.000	09:50:02
0.100 Cal Check	0.096	09:50:22
Air Blank	0.000	09:50:52

Preventative Maintenance
Performed -- Timer Reset