

Routing: YCAO

EXHIBIT G-2
THIS REPORT PREPARED PURSUANT TO DUTY IMPOSED BY A.A.C. R13-10-104 (A)

ARIZONA DEPARTMENT OF PUBLIC SAFETY
STANDARD QUALITY ASSURANCE PROCEDURES
INTOXILYZER MODEL 8000

STANDARD CALIBRATION CHECK PROCEDURE

QA SPECIALIST S. Hutton AGENCY Prescott P.D.

DATE 4-1-06 TIME 2125

INTOXILYZER SERIAL # 80-010337 LOCATION 222 S. MARICOPA ST.

BT 1. Ensure that gas tank is attached to instrument and contains a standard alcohol concentration solution 0.100 AC.

OR

Pour a standard alcohol concentration solution 0.100 AC, into a clean dry simulator and assemble the simulator. Ensure that a tight seal has been made. Turn on the simulator and allow temperature to reach $34^{\circ}\text{C} \pm 0.2^{\circ}\text{C}$

BT
BT

2. Intoxilyzer 8000 display reads "PUSH BUTTON TO START"

3. Go to the "Control Testing Menu". Select "D" for dry control test or "W" for wet control test. After selection is made press ENTER.

BT

4. Air blank completed.

BT

5. Calibration check completed. Test results 0. 097 AC.

BT

6. Air blank completed.

BT

7. Remove printed record. Attach the record to the completed checklist.

SIGNATURE S. Hutton 372

Prescott PD CHANGE 5-18-06
JD NAUTER #127

INTOXILYZER 8000
Location: PRESCOTT P.D
Serial Number: 80-000337
Core Version: 8105.48
01/10/2016 21:27:04

Standard Lot#: 25514100A2
Last Changed By: M. MEDINA #356

QAS: S. HILTON #372
PRESCOTT PD

Operator: S. HILTON #372
PRESCOTT PD

Test	g/210L	Time
Air Blank	0.000	21:27:37
0.100 Cal Check	0.097	21:27:58
Air Blank	0.000	21:28:28

Preventative Maintenance
Performed -- Timer Reset

SHT 372

EXHIBIT G-4
THIS REPORT PREPARED PURSUANT TO DUTY IMPOSED BY A.A.C. R13-10-104 (A)

ARIZONA DEPARTMENT OF PUBLIC SAFETY
STANDARD QUALITY ASSURANCE PROCEDURES
INTOXILYZER MODEL 8000

STANDARD QUALITY ASSURANCE PROCEDURE

QA SPECIALIST S. BROWN AGENCY PRESCOTT P.D.
DATE 1/10/16 TIME 2:29
INTOXILYZER SERIAL # 80-00337 LOCATION 222 S. MAIN ST.

1. Display Reads "PUSH BUTTON TO START"

DIAGNOSTIC TESTS

1. Clock time check.
 2. Date check.

OPERATIONAL TESTS

1. Alcohol-free subject test result 0. 900 AC.
 2. Error recognition logic system, functioning.
Not a Successfully Completed Test Sequence printed
 3. Proper sample recognition system.
Not a Successfully Completed Test Sequence printed
Deficient sample printed.
 4. Standard Calibration Check standard 0. 100 AC. Result 0. 027 AC.

Instrument is operating properly and accurately. YES NO

COMMENTS

FUNCTIONING PROPERLY

SIGNATURE S. BROWN

INTOXILYZER 8000
Location: PRESCOTT P.D
Serial Number: 80-000337
Core Version: 8105.48
01/10/2016

21:29:45

Standard Lot#: 25514100A2
Last Changed By: M. MEDINA #356

QAS: S. HILTON #372
PRESCOTT PD


Operator: S. HILTON #372
PRESCOTT PD

Subject: SQAP,TEST
DOB: 12/19/1973
Sex: M Weight: 175

15 Minute Deprivation Period? Yes

Test	g/210L	Time
Air Blank	0.000	21:31:07
Diagnostic Test	Pass	21:31:41
Air Blank	0.000	21:32:08
0.100 Cal Check	0.095	21:32:29
Air Blank	0.000	21:32:59
Subject Test	0.000	21:33:16
Air Blank	0.000	21:33:43
Wait	ABT*	
Air Blank	0.000	21:34:17

*Sequence Aborted



Not a Successfully Completed Test Sequence

INTOXILYZER 8000
Location: PRESCOTT P.D
Serial Number: 80-000337
Core Version: 8105.48
01/10/2016

21:34:46

Standard Lot#: 25514100A2
Last Changed By: M. MEDINA #356

QAS: S. HILTON #372
PRESCOTT PD

Operator: S. HILTON #372
PRESCOTT PD

Subject: SQAP, TEST
DOB: 12/19/1973
Sex: M Weight: 175

15 Minute Deprivation Period? Yes

Test	g/210L	Time
Air Blank	IPS*	21:35:28
Air Blank	0.000	21:35:56

*Improper Sample

Air 372

Not a Successfully Completed Test Sequence

INTOXILYZER 8000

Location: PRESCOTT P.D

Serial Number: 80-000337

Core Version: 8105.48

01/10/2016

21:36:21

Standard Lot#: 25514100A2

Last Changed By: M. MEDINA #356

QAS: S. HILTON #372

PRESCOTT PD

Operator: S. HILTON #372

PRESCOTT PD

Subject: SQAP, TEST

DOB: 12/19/1973

Sex: M Weight: 175

15 Minute Deprivation Period? Yes

Test	g/210L	Time
Air Blank	RFI*	21:37:01
Air Blank	0.000	21:37:30

*RFI Detect

SMA 372

Not a Successfully Completed Test Sequence

INTOXILYZER 8000
Location: PRESCOTT P.D
Serial Number: 80-000337
Core Version: 8105.48
01/10/2016

21:37:50

Standard Lot#: 25514100A2
Last Changed By: M. MEDINA #356

QAS: S. HILTON #372
PRESCOTT PD

Operator: S. HILTON #372
PRESCOTT PD

Subject: SQAP,TEST
DOB: 12/19/1973
Sex: M Weight: 175

15 Minute Deprivation Period? Yes

Test	g/210L	Time
Air Blank	0.000	21:38:58
Diagnostic Test	Pass	21:39:32
Air Blank	0.000	21:39:59
0.100 Cal Check	0.095	21:40:20
Air Blank	0.000	21:40:50
Subject Test	0.277*	21:43:57
Air Blank	0.000	21:44:39
Air Blank	0.000	21:45:06
Subject Test	XXX**	21:45:34
Air Blank	0.000	21:46:10

*Deficient Sample

**Invalid Sample

Begin new deprivation period



Not a Successfully Completed Test Sequence

EXHIBIT G-2
THIS REPORT PREPARED PURSUANT TO DUTY IMPOSED BY A.A.C. R13-10-104 (A)

ARIZONA DEPARTMENT OF PUBLIC SAFETY

STANDARD QUALITY ASSURANCE PROCEDURES
INTOXILYZER MODEL 8000

STANDARD CALIBRATION CHECK PROCEDURE

QA SPECIALIST S. Hutton AGENCY Prescott P.D.
DATE 1/10/16 TIME 2:44
INTOXILYZER SERIAL # 80-000337 LOCATION 222 S. ARIZONA ST.

- 1. Ensure that gas tank is attached to instrument and contains a standard alcohol concentration solution 1.00 AC.
OR
Pour a standard alcohol concentration solution _____ AC, into a clean dry simulator and assemble the simulator. Ensure that a tight seal has been made. Turn on the simulator and allow temperature to reach 34°C ± 0.2°C
- 2. Intoxilyzer 8000 display reads "PUSH BUTTON TO START"
- 3. Go to the "Control Testing Menu". Select "D" for dry control test or "W" for wet control test. After selection is made press ENTER.
- 4. Air blank completed.
- 5. Calibration check completed. Test results 0.994 AC.
- 6. Air blank completed.
- 7. Remove printed record. Attach the record to the completed checklist.

SIGNATURE S. Hutton

Prescott PD CHANGE 5-18-06
7-10-000 R127

