

EXHIBIT G-2  
THIS REPORT PREPARED PURSUANT TO DUTY IMPOSED BY A.A.C. R13-10-104 (A)

ARIZONA DEPARTMENT OF PUBLIC SAFETY  
STANDARD QUALITY ASSURANCE PROCEDURES  
INTOXILYZER MODEL 8000

STANDARD CALIBRATION CHECK PROCEDURE

QA SPECIALIST T. A. [Signature] 252 AGENCY PRESIDENT PD  
DATE 05/28/16 TIME 1919  
INTOXILYZER SERIAL # 80-000387 LOCATION PRESIDENT PD

- (✓) 1. Ensure that gas tank is attached to instrument and contains a standard alcohol concentration solution. 100 AC.  
OR  
Pour a standard alcohol concentration solution AC into a clean dry simulator and assemble the simulator. Ensure that a tight seal has been made. Turn on the simulator and allow temperature to reach  $34^{\circ}\text{C} \pm 0.2^{\circ}\text{C}$ .
- (✓) 2. Intoxilyzer 8000 display reads "PUSH BUTTON TO START"
- (✓) 3. Go to the "Control Testing Menu". Select D for dry control test or "W" for wet control test. After selection is made press ENTER.
- (✓) 4. Air blank completed.
- (✓) 5. Calibration check completed. Test results 0.098 AC.
- (✓) 6. Air blank completed.
- (✓) 7. Remove printed record. Attach the record to the completed checklist.

SIGNATURE [Signature] 252

INTOXILYZER 8000

Location: PRESCOTT P.D

Serial Number: 80-000337

Core Version: 8105.48

05/28/2016

19:19:09

-----  
Standard Lot#: 22315100A1

Last Changed By: S. HILTON #372  
-----

QAS: T. ATTENBERGER #252

PRESCOTT PD  
-----

Operator: T. ATTENBERGER #252

PRESCOTT PD  
-----

Test	g/210L	Time
Air Blank	0.000	19:19:43
0.100 Cal Check	0.098	19:20:04
Air Blank	0.000	19:20:34

Preventative Maintenance  
Performed -- Timer Reset

A handwritten signature in black ink, appearing to read 'T. Attenberger' followed by the number '252'. The signature is stylized and cursive.

EXHIBIT G-4  
THIS REPORT PREPARED PURSUANT TO DUTY IMPOSED BY A.A.C. R13-10-104 (A)

ARIZONA DEPARTMENT OF PUBLIC SAFETY

STANDARD QUALITY ASSURANCE PROCEDURES  
INTOXILYZER MODEL 8000

STANDARD QUALITY ASSURANCE PROCEDURE

QA SPECIALIST T. Anderson 252 AGENCY Prescott PD

DATE 05/28/16 TIME 1922

INTOXILYZER SERIAL # 80-030337 LOCATION Prescott PD

() 1. Display Reads "PUSH BUTTON TO START"

DIAGNOSTIC TESTS

- () 1. Clock time check.
- () 2. Date check.

OPERATIONAL TESTS

- () 1. Alcohol-free subject test result 0. 000 AC.
- () 2. Error recognition logic system functioning.  
Not a Successfully Completed Test Sequence printed
- () 3. Proper sample recognition system.  
Not a Successfully Completed Test Sequence printed  
Deficient sample printed.
- () 4. Standard Calibration Check standard 0. 100 AC. Result 0. 095 AC.

Instrument is operating properly and accurately. YES  NO

COMMENTS

ALL FUNCTIONS CORRECTLY.

SIGNATURE [Signature]

INTOXILYZER 8000  
Location: PRESCOTT P.D  
Serial Number: 80-000337  
Core Version: 8105.48  
05/28/2016

19:22:11

---

Standard Lot#: 22315100A1  
Last Changed By: S. HILTON #372

---

QAS: T. ATTENBERGER #252  
PRESCOTT PD

---

Operator: T. ATTENBERGER #252  
PRESCOTT PD

---

Subject: SQAP,TEST  
DOB: 01/01/1900  
Sex: M Weight: 300

---

15 Minute Deprivation Period? Yes

---

Test	g/210L	Time
Air Blank	0.000	19:23:05
Diagnostic Test	Pass	19:23:39
Air Blank	0.000	19:24:06
0.100 Cal Check	0.098	19:24:27
Air Blank	0.000	19:24:57
Subject Test	0.000	19:25:20
Air Blank	0.000	19:25:47
Wait	ABT*	
Air Blank	0.000	19:26:23

\*Sequence Aborted



Not a Successfully Completed Test Sequence

INTOXILYZER 8000

Location: PRESCOTT P.D

Serial Number: 80-000337

Core Version: 8105.48

05/28/2016

19:26:48

---

Standard Lot#: 22315100A1

Last Changed By: S. HILTON

#372

---

QAS: T. ATTENBERGER

#252

PRESCOTT PD

---

Operator: T. ATTENBERGER

#252

PRESCOTT PD

---

Subject: SQAP,TEST

DOB: 01/01/1900

Sex: M

Weight: 300

---

15 Minute Deprivation Period? Yes

---

Test	g/210L	Time
Air Blank	IPS*	19:27:27
Air Blank	0.000	19:27:55

\*Improper Sample

A handwritten signature in black ink, appearing to be 'T. Attenberger', with a stylized flourish at the end.

**Not a Successfully Completed Test Sequence**

INTOXILYZER 8000

Location: PRESCOTT P.D

Serial Number: 80-000337

Core Version: 8105.48

05/28/2016

19:28:15

---

Standard Lot#: 22315100A1

Last Changed By: S. HILTON

#372

---

QAS: T. ATTENBERGER

#252

PRESCOTT PD

---

Operator: T. ATTENBERGER

#252

PRESCOTT PD

---

Subject: SQAP, TEST

DOB: 01/01/1900

Sex: M

Weight: 300

---

15 Minute Deprivation Period? Yes

---

Test	g/210L	Time
Air Blank	RFI*	19:28:59
Air Blank	0.000	19:29:28

\*RFI Detect

A handwritten signature in black ink, appearing to be "T. Attenberger" with the number "252" written below it.

Not a Successfully Completed Test Sequence

INTOXILYZER 8000

Location: PRESCOTT P.D

Serial Number: 80-000337

Core Version: 8105.48

05/28/2016

19:29:55

---

Standard Lot#: 22315100A1

Last Changed By: S. HILTON #372

---

QAS: T. ATTENBERGER #252

PRESCOTT PD

---

Operator: T. ATTENBERGER #252

PRESCOTT PD

---

Subject: SQAP,TEST

DOB: 01/01/1900

Sex: M

Weight: 300

---

15 Minute Deprivation Period? Yes


---

Test	g/210L	Time
Air Blank	0.000	19:30:48
Diagnostic Test	Pass	19:31:22
Air Blank	0.000	19:31:49
0.100 Cal Check	0.096	19:32:10
Air Blank	0.000	19:32:41
Subject Test	0.022*	19:35:48
Air Blank	0.000	19:36:19
Air Blank	0.000	19:36:46
Subject Test	XXX**	19:38:06
Air Blank	0.000	19:38:40

\*Deficient Sample

\*\*Invalid Sample

Begin new deprivation period



**Not a Successfully Completed Test Sequence**

EXHIBIT G-2  
THIS REPORT PREPARED PURSUANT TO DUTY IMPOSED BY A.A.C. R13-10-104 (A)

ARIZONA DEPARTMENT OF PUBLIC SAFETY  
STANDARD QUALITY ASSURANCE PROCEDURES  
INTOXILYZER MODEL 8000  
STANDARD CALIBRATION CHECK PROCEDURE

QA SPECIALIST T. Andrew Badger 252 AGENCY Prescott PD  
DATE 05/28/16 TIME 1939  
INTOXILYZER SERIAL # 80-020337 LOCATION Prescott PD

- (✓) 1. Ensure that gas tank is attached to instrument and contains a standard alcohol concentration solution 100 AC.  
OR  
Pour a standard alcohol concentration solution \_\_\_\_\_ AC, into a clean dry simulator and assemble the simulator. Ensure that a tight seal has been made. Turn on the simulator and allow temperature to reach  $34^{\circ}\text{C} \pm 0.2^{\circ}\text{C}$
- (✓) 2. Intoxilyzer 8000 display reads "PUSH BUTTON TO START"
- (✓) 3. Go to the "Control Testing Menu". Select (D) for dry control test or "W" for wet control test. After selection is made press ENTER.
- (✓) 4. Air blank completed.
- (✓) 5. Calibration check completed. Test results 0. 095 AC.
- (✓) 6. Air blank completed.
- (✓) 7. Remove printed record. Attach the record to the completed checklist.

SIGNATURE [Signature] 252

INTOXILYZER 8000

Location: PRESCOTT P.D

Serial Number: 80-000337

Core Version: 8105.48

05/28/2016

19:39:19

-----  
Standard Lot#: 22315100A1

Last Changed By: S. HILTON #372  
-----

QAS: T. ATTENBERGER #252

PRESCOTT PD  
-----

Operator: T. ATTENBERGER #252

PRESCOTT PD  
-----

Test	g/210L	Time
Air Blank	0.000	19:39:50
0.100 Cal Check	0.095	19:40:11
Air Blank	0.000	19:40:41

Preventative Maintenance  
Performed -- Timer Reset

