

EXHIBIT G-2
THIS REPORT PREPARED PURSUANT TO DUTY IMPOSED BY A.A.C. R13-10-104 (A)

ARIZONA DEPARTMENT OF PUBLIC SAFETY
STANDARD QUALITY ASSURANCE PROCEDURES
INTOXILYZER MODEL 8000

STANDARD CALIBRATION CHECK PROCEDURE

QA SPECIALIST S. HILTON AGENCY PRESOTT P.D.
DATE 10/12/17 TIME 1425
INTOXILYZER SERIAL # 80-000337 LOCATION 222 S. MARINA ST.

- Smt* 1. Ensure that gas tank is attached to instrument and contains a standard alcohol concentration solution .100 AC.
OR
Pour a standard alcohol concentration solution _____ AC, into a clean dry simulator and assemble the simulator. Ensure that a tight seal has been made. Turn on the simulator and allow temperature to reach $34^{\circ}\text{C} \pm 0.2^{\circ}\text{C}$
- Smt* 2. Intoxilyzer 8000 display reads "PUSH BUTTON TO START"
Smt 3. Go to the "Control Testing Menu". Select "D" for dry control test or "W" for wet control test. After selection is made press ENTER.
Smt 4. Air blank completed.
Smt 5. Calibration check completed. Test results 0. 099 AC.
Smt 6. Air blank completed.
Smt 7. Remove printed record. Attach the record to the completed checklist.

SIGNATURE S. Hilton

DPS Form Exh G-2 (Rev 05-01)

Prescott PID CHANGE 5-12-06

INTOXILYZER 8000
Location: PRESCOTT P.D
Serial Number: 80-000337
Core Version: 8105.48
10/12/2017 14:25:53

Standard Lot#: 27016100A1
Last Changed By: S. HILTON #372

QAS: S. HILTON #372
PRESCOTT PD

Operator: S. HILTON #372
PRESCOTT PD

Test	g/210L	Time
Air Blank	0.000	14:26:25
0.100 Cal Check	0.099	14:26:46
Air Blank	0.000	14:27:17

Preventative Maintenance
Performed -- Timer Reset



EXHIBIT G-2
THIS REPORT PREPARED PURSUANT TO DUTY IMPOSED BY A.A.C. R13-10-104 (A)

ARIZONA DEPARTMENT OF PUBLIC SAFETY
STANDARD QUALITY ASSURANCE PROCEDURES
INTOXILYZER MODEL 8000

STANDARD CALIBRATION CHECK PROCEDURE

QA SPECIALIST S. Helton AGENCY Prescott P.D.
DATE 10/28/17 TIME 1400
INTOXILYZER SERIAL # 80-000337 LOCATION 222 S. M. ARENA

- SH* 1. Ensure that gas tank is attached to instrument and contains a standard alcohol concentration solution .100 AC.
OR
Pour a standard alcohol concentration solution _____ AC, into a clean dry simulator and assemble the simulator. Ensure that a tight seal has been made. Turn on the simulator and allow temperature to reach $34^{\circ}\text{C} \pm 0.2^{\circ}\text{C}$
- SH* 2. Intoxilyzer 8000 display reads "PUSH BUTTON TO START"
- SH* 3. Go to the "Control Testing Menu". Select "D" for dry control test or "W" for wet control test. After selection is made press ENTER.
- SH* 4. Air blank completed.
- SH* 5. Calibration check completed. Test results 0. 100 AC.
- SH* 6. Air blank completed.
- SH* 7. Remove printed record. Attach the record to the completed checklist.

SIGNATURE *S. Helton*

Prescott PD CHANGE 5-18-06

INTOXILYZER 8000

Location: PRESCOTT P.D

Serial Number: 80-000337

Core Version: 8105.48

10/28/2017

14:00:26

Standard Lot#: 27016100A1

Last Changed By: S. HILTON #372

QAS: S. HILTON #372

PRESCOTT PD

Operator: S. HILTON #372

PRESCOTT PD

Test	g/210L	Time
Air Blank	0.000	14:00:58
0.100 Cal Check	0.100	14:01:19
Air Blank	0.000	14:01:49

Preventative Maintenance
Performed -- Timer Reset

SM

EXHIBIT G-4
THIS REPORT PREPARED PURSUANT TO DUTY IMPOSED BY A.A.C. R13-10-104 (A)

ARIZONA DEPARTMENT OF PUBLIC SAFETY

STANDARD QUALITY ASSURANCE PROCEDURES
INTOXILYZER MODEL 8000

STANDARD QUALITY ASSURANCE PROCEDURE

QA SPECIALIST S. Hutton AGENCY Prescott P.D.
DATE 10/28/17 TIME 1403
INTOXILYZER SERIAL # 80-000337 LOCATION 222 S. MARINA ST.

1. Display Reads "PUSH BUTTON TO START"

DIAGNOSTIC TESTS

1. Clock time check.
 2. Date check.

OPERATIONAL TESTS

1. Alcohol-free subject test result 0. 000 AC.
 2. Error recognition logic system functioning.
Not a Successfully Completed Test Sequence printed
 3. Proper sample recognition system.
Not a Successfully Completed Test Sequence printed
Deficient sample printed.
 4. Standard Calibration Check standard 0. 100 AC. Result 0. 094 AC.

Instrument is operating properly and accurately. YES NO

COMMENTS

functioning properly

SIGNATURE S. Hutton

INTOXILYZER 8000
Location: PRESCOTT P.D
Serial Number: 80-000337
Core Version: 8105.48
10/28/2017

14:03:36

Standard Lot#: 27016100A1
Last Changed By: S. HILTON #372

QAS: S. HILTON #372
PRESCOTT PD

Operator: S. HILTON #372
PRESCOTT PD

Subject: SQAP,TEST
DOB: 12/19/1973
Sex: M Weight: 175

15 Minute Deprivation Period? Yes

Test	g/210L	Time
Air Blank	0.000	14:04:42
Diagnostic Test	Pass	14:05:16
Air Blank	0.000	14:05:43
0.100 Cal Check	0.099	14:06:04
Air Blank	0.000	14:06:34
Subject Test	0.000	14:06:49
Air Blank	ABT*	14:07:16
Air Blank	0.000	14:07:45

*Sequence Aborted



Not a Successfully Completed Test Sequence

INTOXILYZER 8000

Location: PRESCOTT P.D

Serial Number: 80-000337

Core Version: 8105.48

10/28/2017

14:08:09

Standard Lot#: 27016100A1

Last Changed By: S. HILTON

#372

QAS: S. HILTON

#372

PRESCOTT PD

Operator: S. HILTON

#372

PRESCOTT PD

Subject: SQAP, TEST

DOB: 12/19/1973

Sex: M

Weight: 175

15 Minute Deprivation Period? Yes

Test	g/210L	Time
Air Blank	IPS*	14:08:46
Air Blank	0.000	14:09:14

*Improper Sample



Not a Successfully Completed Test Sequence

INTOXILYZER 8000

Location: PRESCOTT P.D

Serial Number: 80-000337

Core Version: 8105.48

10/28/2017

14:10:05

Standard Lot#: 27016100A1

Last Changed By: S. HILTON

#372

QAS: S. HILTON

#372

PRESCOTT PD

Operator: S. HILTON

#372

PRESCOTT PD

Subject: SQAP, TEST

DOB: 12/19/1973

Sex: M

Weight: 175

15 Minute Deprivation Period? Yes

Test	g/210L	Time
Air Blank	RFI*	14:10:51
Air Blank	0.000	14:11:19

*RFI Detect

SMA

Not a Successfully Completed Test Sequence

INTOXILYZER 8000
Location: PRESCOTT P.D
Serial Number: 80-000337
Core Version: 8105.48
10/28/2017

14:11:55

Standard Lot#: 27016100A1
Last Changed By: S. HILTON #372

QAS: S. HILTON #372
PRESCOTT PD

Operator: S. HILTON #372
PRESCOTT PD

Subject: SQAP, TEST
DOB: 12/19/1973
Sex: M Weight: 175

15 Minute Deprivation Period? Yes

Test	g/210L	Time
Air Blank	0.000	14:12:51
Diagnostic Test	Pass	14:13:25
Air Blank	0.000	14:13:52
0.100 Cal Check	0.099	14:14:13
Air Blank	0.000	14:14:43
Subject Test	0.264*	14:17:51
Air Blank	0.000	14:18:34
Air Blank	0.000	14:19:01
Subject Test	XXX**	14:19:26
Air Blank	0.000	14:20:12

*Deficient Sample

**Invalid Sample

Begin new deprivation period



Not a Successfully Completed Test Sequence

EXHIBIT G-2
THIS REPORT PREPARED PURSUANT TO DUTY IMPOSED BY A.A.C. R13-10-104 (A)

ARIZONA DEPARTMENT OF PUBLIC SAFETY
STANDARD QUALITY ASSURANCE PROCEDURES
INTOXILYZER MODEL 8000
STANDARD CALIBRATION CHECK PROCEDURE

QA SPECIALIST S. HAZEN AGENCY Prescott P.D.
DATE 10/28/17 TIME 1420
INTOXILYZER SERIAL # 50-000337 LOCATION 222 SIMONA ST.

- (S) 1. Ensure that gas tank is attached to instrument and contains a standard alcohol concentration solution .100 AC.
OR
Pour a standard alcohol concentration solution _____ AC, into a clean dry simulator and assemble the simulator. ~~Ensure that a tight seal has been made. Turn on the simulator and allow temperature to reach 34°C ± 0.2°C~~
- (S) 2. Intoxilyzer 8000 display reads "PUSH BUTTON TO START"
- (S) 3. Go to the "Control Testing Menu". Select "D" for dry control test or "W" for wet control test. After selection is made press ENTER.
- (S) 4. Air blank completed.
- (S) 5. Calibration check completed. Test results 0. 094 AC.
- (S) 6. Air blank completed.
- (S) 7. Remove printed record. Attach the record to the completed checklist.

SIGNATURE S. HAZEN

Prescott PD CHANGE 5-18-06
7 10 200 11:17

INTOXILYZER 8000
Location: PRESCOTT P.D
Serial Number: 80-000337
Core Version: 8105.48
10/28/2017 14:20:52

Standard Lot#: 27016100A1
Last Changed By: S. HILTON #372

QAS: S. HILTON #372
PRESCOTT PD

Operator: S. HILTON #372
PRESCOTT PD

Test	g/210L	Time
Air Blank	0.000	14:21:24
0.100 Cal Check	0.094	14:21:44
Air Blank	0.000	14:22:15

Preventative Maintenance
Performed -- Timer Reset

