

EXHIBIT G-2
THIS REPORT PREPARED PURSUANT TO DUTY IMPOSED BY A.A.C. R13-10-104 (A)

ARIZONA DEPARTMENT OF PUBLIC SAFETY
STANDARD QUALITY ASSURANCE PROCEDURES
INTOXILYZER MODEL 8000
STANDARD CALIBRATION CHECK PROCEDURE

QA SPECIALIST I.D. BARDON AGENCY PRESCOTT P.D.
DATE MAY 24, 2018 TIME 1955
INTOXILYZER SERIAL # 80-000337 LOCATION PPD

- (✓) 1. Ensure that gas tank is attached to instrument and contains a standard alcohol concentration solution .100 AC.
OR
Pour a standard alcohol concentration solution _____ AC, into a clean dry simulator and assemble the simulator. Ensure that a tight seal has been made. Turn on the simulator and allow temperature to reach $34^{\circ}\text{C} \pm 0.2^{\circ}\text{C}$
- (✓) 2. Intoxilyzer 8000 display reads "PUSH BUTTON TO START"
- (✓) 3. Go to the "Control Testing Menu". Select "D" for dry control test or "W" for wet control test. After selection is made press ENTER.
- (✓) 4. Air blank completed.
- (✓) 5. Calibration check completed. Test results 0.0986 AC.
- (✓) 6. Air blank completed.
- (✓) 7. Remove printed record. Attach the record to the completed checklist.

SIGNATURE [Signature]

DPS Form Exh G-2 (Rev 05-01)

PRESCOTT P.D. CHANGE 5-12-06
SD NAVIER #127

INTOXILYZER 8000
Location: PRESCOTT P.D
Serial Number: 80-000337
Core Version: 8105.48
05/24/2018 19:58:25

Standard Lot#: 27016100A1
Last Changed By: D. BARDON #351

QAS: D. BARDON #351
PRESCOTT PD

Operator: D. BARDON #351
PRESCOTT PD

Test	g/210L	Time
Air Blank	0.000	19:59:19
0.100 Cal Check	0.098	19:59:40
Air Blank	0.000	20:00:10

Preventative Maintenance
Performed -- Timer Reset

EXHIBIT G-4
THIS REPORT PREPARED PURSUANT TO DUTY IMPOSED BY A.A.C. R13-10-104 (A)

ARIZONA DEPARTMENT OF PUBLIC SAFETY

STANDARD QUALITY ASSURANCE PROCEDURES
INTOXILYZER MODEL 8000

STANDARD QUALITY ASSURANCE PROCEDURE

QA SPECIALIST D. BARDON AGENCY PRESCOTT P.D.
DATE MAY 24, 2018 TIME 2004
INTOXILYZER SERIAL # 80-000337 LOCATION PPD

() 1. Display Reads "PUSH BUTTON TO START"

DIAGNOSTIC TESTS

- () 1. Clock time check.
- () 2. Date check.

OPERATIONAL TESTS

- () 1. Alcohol-free subject test result 0. 000 AC.
- () 2. Error recognition logic system functioning.
Not a Successfully Completed Test Sequence printed
- () 3. Proper sample recognition system.
Not a Successfully Completed Test Sequence printed
Deficient sample printed.
- () 4. Standard Calibration Check standard 0. 100 AC. Result 0. 094 AC.

Instrument is operating properly and accurately. YES NO

COMMENTS

SIGNATURE [Signature]

PRESCOTT PD CANAL 5 18:06

INTOXILYZER 8000
Location: PRESCOTT P.D
Serial Number: 80-000337
Core Version: 8105.48
05/24/2018

20:06:45

Standard Lot#: 27016100A1
Last Changed By: D. BARDON #351

QAS: D. BARDON #351
PRESCOTT PD

Operator: D. BARDON #351
PRESCOTT PD

Subject: SQAP,TEST
DOB: 11/11/1999
Sex: M Weight: 200

15 Minute Deprivation Period? Yes

Test	g/210L	Time
Air Blank	0.000	20:07:36
Diagnostic Test	Pass	20:08:10
Air Blank	0.000	20:08:37
0.100 Cal Check	0.097	20:08:58
Air Blank	0.000	20:09:29
Subject Test	0.000	20:09:48
Air Blank	ABT*	20:09:55
Air Blank	IPS**	20:10:05

*Sequence Aborted

**Improper Sample

Not a Successfully Completed Test Sequence

INTOXILYZER 8000
Location: PRESCOTT P.D
Serial Number: 80-000337
Core Version: 8105.48
05/24/2018

20:12:59

Standard Lot#: 27016100A1
Last Changed By: D. BARDON #351

QAS: D. BARDON #351
PRESCOTT PD

Operator: D. BARDON #351
PRESCOTT PD

Subject: SQAP, TEST
DOB: 11/11/1999
Sex: M Weight: 200

15 Minute Deprivation Period? Yes

Test	g/210L	Time
Air Blank	0.000	20:13:50
Diagnostic Test	Pass	20:14:24
Air Blank	0.000	20:14:52
0.100 Cal Check	0.096	20:15:13
Air Blank	0.000	20:15:43
Subject Test	0.037*	20:18:50
Air Blank	0.000	20:19:23
Air Blank	0.000	20:19:51
Subject Test	XXX**	20:20:14
Air Blank	RFI***	20:20:21

*Deficient Sample

**Invalid Sample

Begin new deprivation period

***RFI Detect

Not a Successfully Completed Test Sequence

EXHIBIT G-2
THIS REPORT PREPARED PURSUANT TO DUTY IMPOSED BY A.A.C. R13-10-104 (A)

ARIZONA DEPARTMENT OF PUBLIC SAFETY
STANDARD QUALITY ASSURANCE PROCEDURES
INTOXILYZER MODEL 8000
STANDARD CALIBRATION CHECK PROCEDURE

QA SPECIALIST D. BARON AGENCY PRESCOTT P.D
DATE MAY 24, 2018 TIME 2024
INTOXILYZER SERIAL # 80-000337 LOCATION PPD

- (✓) 1. Ensure that gas tank is attached to instrument and contains a standard alcohol concentration solution .100 AC.
OR
Pour a standard alcohol concentration solution _____ AC, into a clean dry simulator and assemble the simulator. Ensure that a tight seal has been made. Turn on the simulator and allow temperature to reach 34°C ± 0.2°C
- (✓) 2. Intoxilyzer 8000 display reads "PUSH BUTTON TO START"
- (✓) 3. Go to the "Control Testing Menu". Select "D" for dry control test or "W" for wet control test. After selection is made press ENTER.
- (✓) 4. Air blank completed.
- () 5. Calibration check completed. Test results 0. 096 AC.
- (✓) 6. Air blank completed.
- (✓) 7. Remove printed record. Attach the record to the completed checklist.

SIGNATURE [Signature]

DPS Form Exh G-2 (Rev 05-01)

Prescott PD CHANGE 5-12-06
SD NAVTEL #127

INTOXILYZER 8000

Location: PRESCOTT P.D

Serial Number: 80-000337

Core Version: 8105.48

05/24/2018

20:23:53

Standard Lot#: 27016100A1

Last Changed By: D. BARDON #351

QAS: D. BARDON #351

PRESCOTT PD

Operator: D. BARDON #351

PRESCOTT PD

Test	g/210L	Time
Air Blank	0.000	20:24:24
0.100 Cal Check	0.096	20:24:44
Air Blank	0.000	20:25:15

Preventative Maintenance
Performed -- Timer Reset