

EXHIBIT G-2
THIS REPORT PREPARED PURSUANT TO DUTY IMPOSED BY A.A.C. R13-10-104 (A)

ARIZONA DEPARTMENT OF PUBLIC SAFETY
STANDARD QUALITY ASSURANCE PROCEDURES
INTOXILYZER MODEL 8000

STANDARD CALIBRATION CHECK PROCEDURE

QA SPECIALIST DAN RAISS AGENCY YAVAPAI CO. S.O.
DATE 4-28-08 TIME 1615
INTOXILYZER SERIAL # 80-000289 LOCATION YCSO - MOBILE

- (DL) 1. Ensure that gas tank is attached to instrument and contains a standard alcohol concentration solution 100 AC.
OR
Pour a standard alcohol concentration solution 100 AC, into a clean dry simulator and assemble the simulator. Ensure that a ~~light~~ seal has been made. Turn on the simulator and allow temperature to reach 34°C ± 0.2°C
- (DL) 2. Intoxilyzer 8000 display reads "PUSH BUTTON TO START"
- (DL) 3. Go to the "Control Testing Menu". Select "D" for dry control test or "W" for wet control test. After selection is made press ENTER.
- (DL) 4. Air blank completed.
- (DL) 5. Calibration check completed. Test results 0. 096 AC.
- (DL) 6. Air blank completed.
- (DL) 7. Remove printed record. Attach the record to the completed checklist.

SIGNATURE Daniel G. Raiss #5-31/OAS #33026

DPS Form Exh G-2 (Rev 05-01)

Daniel G. Raiss #5-31
Operator's Signature
OAS #33026

INTOXILYZER 8000
 Location: YCSO-MOBILE
 Serial Number: 80-000289
 04/28/2008 16:33:39

Standard Lot#: 6339111
 Last Changed By: D. RAISS #3693

DPS: D. RAISS #3693
 YAVAPAI CSO

Test	g/210L	Time
Air Blank	0.000	16:34:07
Cal Check	0.096	16:34:25
Air Blank	0.000	16:34:55

EXHIBIT G-4
THIS REPORT PREPARED PURSUANT TO DUTY IMPOSED BY A.A.C. R13-10-104 (A)

ARIZONA DEPARTMENT OF PUBLIC SAFETY

STANDARD QUALITY ASSURANCE PROCEDURES
INTOXILYZER MODEL 8000

STANDARD QUALITY ASSURANCE PROCEDURE

QA SPECIALIST DAN RAISS AGENCY YAVAPAI Co. S.O.
DATE 4-28-08 TIME 1615
INTOXILYZER SERIAL # 80-060289 LOCATION YCSO-MOBILE

(D) 1. Display Reads "PUSH BUTTON TO START"

DIAGNOSTIC TESTS

- (D) 1. Clock time check.
(D) 2. Date check.

OPERATIONAL TESTS

- (D) 1. Alcohol-free subject test result 0. 000 AC.
(D) 2. Error recognition logic system functioning.
Not a Successfully Completed Test Sequence printed
(D) 3. Proper sample recognition system.
Not a Successfully Completed Test Sequence printed
Deficient sample printed.
(D) 4. Standard Calibration Check standard 0. 100 AC. Result 0. 096 AC.

Instrument is operating properly and accurately. YES NO

COMMENTS

ALL OK.

SIGNATURE

Daniel Q Raimo #5-31/OAS #33026

INTOXILYZER 8000
Location: YCSD-MOBILE
Serial Number: 80-000289
04/28/2008 16:18:35

Standard Lot#: 6339011
Last Changed By: D. RAJSS #3693

ONS: D. RAJSS #3693
YAUAPP1 CS0

Operator: D. RAJSS #3693
YAUAPP1 CS0

Subject: SQRP TEST,
DOB: 01/01/1901
Sex: M Weight: 1

15 Minute Deprivation Period? Yes

Test g/210L Time

Air Blank	0.000	16:20:03
Diagnostic Test	Pass	16:20:31
Air Blank	0.000	16:20:55
0.100 Cal Check	0.097	16:21:18
Air Blank	0.000	16:21:47
Subject Test	0.000	16:22:11
Air Blank	0.000	16:22:39
Five Minute Wait	RB1*	
Air Blank	0.000	16:23:18

*Sequence Aborted

Not a Successfully
Completed
Test Sequence

INTOXILYZER 8000
Location: YCSD-MOBILE
Serial Number: 80-000289
04/28/2008 16:23:46

Standard Lot#: 6339011
Last Changed By: D. RAJSS #3693

ONS: D. RAJSS #3693
YAUAPP1 CS0

Operator: D. RAJSS #3693
YAUAPP1 CS0

Subject: SQRP TEST,
DOB: 01/01/1901
Sex: M Weight: 1

15 Minute Deprivation Period? Yes

Test g/210L Time

Air Blank	IPS*	16:24:12
Air Blank	0.000	16:24:40

*Improper Sample

Not a Successfully
Completed
Test Sequence

INTOXILYZER 8000
Location: YCSD-MOBILE
Serial Number: 80-000289
04/28/2008 16:24:47

Standard Lot#: 6339011
Last Changed By: D. RAJSS #3693

ONS: D. RAJSS #3693
YAUAPP1 CS0

Operator: D. RAJSS #3693
YAUAPP1 CS0

Subject: SQRP TEST,
DOB: 01/01/1901
Sex: M Weight: 1

15 Minute Deprivation Period? Yes

Test g/210L Time

Air Blank	RF1*	16:25:13
Air Blank	0.000	16:25:41

*RF1 Detect

Not a Successfully
Completed
Test Sequence

INTOXILYZER 8000
Location: YCSD-MOBILE
Serial Number: 80-000289
04/28/2008 16:25:47

Standard Lot#: 6339011
Last Changed By: D. RAJSS #3693

ONS: D. RAJSS #3693
YAUAPP1 CS0

Operator: D. RAJSS #3693
YAUAPP1 CS0

Subject: SQRP TEST,
DOB: 01/01/1901
Sex: M Weight: 1

15 Minute Deprivation Period? Yes

Test g/210L Time

Air Blank	0.000	16:26:28
Diagnostic Test	Pass	16:26:57
Air Blank	0.000	16:27:24
0.100 Cal Check	0.097	16:27:43
Air Blank	0.000	16:28:12
Subject Test	0.100*	16:31:21
Air Blank	0.000	16:31:52
Air Blank	0.000	16:32:19
Subject Test	XXX**	16:32:32
Air Blank	0.000	16:33:05

*Deficient Sample

**Invalid Sample

Begin new deprivation period

Not a Successfully
Completed
Test Sequence