

**COPY**

EXHIBIT G-2  
THIS REPORT PREPARED PURSUANT TO DUTY IMPOSED BY A.A.C. R13-10-104 (A)

ARIZONA DEPARTMENT OF PUBLIC SAFETY  
STANDARD QUALITY ASSURANCE PROCEDURES  
INTOXILYZER MODEL 8000  
STANDARD CALIBRATION CHECK PROCEDURE

QA SPECIALIST J. Murphy AGENCY Yavapai County Sheriff's Office  
DATE 12/23/2014 TIME 10:23:49  
INTOXILYZER SERIAL # 80-006238 LOCATION DUI Truck

- (X) 1. Ensure that gas tank is attached to instrument and contains a standard alcohol concentration solution .100 AC.  
OR  
Pour a standard alcohol concentration solution        AC, into a clean dry simulator and assemble the simulator. Ensure that a tight seal has been made. Turn on the simulator and allow temperature to reach 34°C ± 0.2°C
- (X) 2. Intoxilyzer 8000 display reads "PUSH BUTTON TO START"
- (X) 3. Go to the "Control Testing Menu". Select "D" for dry control test or "W" for wet control test. After selection is made press ENTER.
- (X) 4. Air blank completed.
- (X) 5. Calibration check completed. Test results 0. 097 AC.
- (X) 6. Air blank completed.
- (X) 7. Remove printed record. Attach the record to the completed checklist.

SIGNATURE J. Murphy NS6/5597

INTOXILYZER 8000  
Location: YAVAPAI COUNTY SO  
Serial Number: 80-006228  
Core Version: 8105.48  
12/23/2014 10:23:49

-----  
Standard Lot#: 12  
Last Changed By: 1. 1 #1  
-----

QAS: J. MURRAY5597 #N56  
YAVAPAI SHERIFFS  
-----

Operator: J. MURRAY5597 #N56  
YAVAPAI SHERIFFS  
-----

Test	g/210L	Time
Air Blank	0.000	10:24:28
0.100 Cal Check	0.097	10:24:49
Air Blank	0.000	10:25:18

Preventative Maintenance  
Performed -- Timer Reset

EXHIBIT G-4  
THIS REPORT PREPARED PURSUANT TO DUTY IMPOSED BY A.A.C. R13-10-104 (A)

ARIZONA DEPARTMENT OF PUBLIC SAFETY  
STANDARD QUALITY ASSURANCE PROCEDURES  
INTOXILYZER MODEL 8000

STANDARD QUALITY ASSURANCE PROCEDURE

QA SPECIALIST J. MURPHY AGENCY Yavapai County Sheriff's Office  
DATE 12/23/2014 TIME 10:27:33  
INTOXILYZER SERIAL # 80-006228 LOCATION DUI Truck

( ) 1. Display Reads "PUSH BUTTON TO START"

DIAGNOSTIC TESTS

- 1. Clock time check.
- 2. Date check.

OPERATIONAL TESTS

- 1. Alcohol-free subject test result 0. 000 AC.
- 2. Error recognition logic system functioning.  
Not a Successfully Completed Test Sequence printed
- 3. Proper sample recognition system.  
Not a Successfully Completed Test Sequence printed  
Deficient sample printed.
- 4. Standard Calibration Check standard 0. 100 AC. Result 0. 096 AC.

Instrument is operating properly and accurately. YES  NO

COMMENTS

Ø13141ØØA3 -1000 PSI  
exp. - 2/1/2016

SIGNATURE J. Murphy NS6/5597

INTOXILYZER 8000  
Location: YAVAPAI COUNTY SO  
Serial Number: 80-006228  
Core Version: 8105.48  
12/23/2014 10:27:33

---

Standard Lot#: 12  
Last Changed By: 1. 1 #1

---

QAS: J. MURRAY5597 #N56  
YAVAPAI SHERIFFS

---

Operator: J. MURRAY5597 #N56  
YAVAPAI SHERIFFS

---

Subject: SQAP,TEST  
DOB: 10/20/1960  
Sex: F Weight: 115

---

15 Minute Deprivation Period? Yes

---

Test	g/210L	Time
Air Blank	0.000	10:29:00
Diagnostic Test	Pass	10:29:35
Air Blank	0.000	10:30:02
0.100 Cal Check	0.098	10:30:23
Air Blank	0.000	10:30:52
Subject Test	0.000	10:31:22
Air Blank	0.000	10:31:49
Wait	ABT*	
Air Blank	RFI**	10:32:06

\*Sequence Aborted

\*\*RFI Detect

**Not a Successfully Completed Test Sequence**

INTOXILYZER 8000  
Location: YAVAPAI COUNTY SO  
Serial Number: 80-006228  
Core Version: 8105.48  
12/23/2014 10:32:11

---

Standard Lot#: 12  
Last Changed By: 1. 1 #1

---

QAS: J. MURRAY5597 #N56  
YAVAPAI SHERIFFS

---

Operator: J. MURRAY5597 #N56  
YAVAPAI SHERIFFS

---

Subject: SQAP,TEST  
DOB: 10/20/1960  
Sex: F Weight: 115

---

15 Minute Deprivation Period? Yes

---

Test	g/210L	Time
Air Blank	IPS*	10:33:16
Air Blank	0.000	10:33:44

\*Improper Sample

**Not a Successfully Completed Test Sequence**

INTOXILYZER 8000  
Location: YAVAPAI COUNTY SO  
Serial Number: 80-006228  
Core Version: 8105.48  
12/23/2014 10:34:20

---

Standard Lot#: 12  
Last Changed By: 1. 1 #1

---

QAS: J. MURRAY5597 #N56  
YAVAPAI SHERIFFS

---

Operator: J. MURRAY5597 #N56  
YAVAPAI SHERIFFS

---

Subject: SQAP,TEST  
DOB: 10/20/1960  
Sex: F Weight: 115

---

15 Minute Deprivation Period? Yes

---

Test	g/210L	Time
Air Blank	0.000	10:35:13
Diagnostic Test	Pass	10:35:48
Air Blank	0.000	10:36:16
0.100 Cal Check	0.097	10:36:36
Air Blank	0.000	10:37:06
Subject Test	0.426*	10:40:13
Air Blank	0.000	10:40:46
Air Blank	0.000	10:41:13
Subject Test	XXX**	10:41:42
Air Blank	0.000	10:42:16

\*Deficient Sample  
\*\*Invalid Sample  
Begin new deprivation period

**Not a Successfully Completed Test Sequence**

EXHIBIT G-2  
THIS REPORT PREPARED PURSUANT TO DUTY IMPOSED BY A.A.C. R13-10-104 (A)

ARIZONA DEPARTMENT OF PUBLIC SAFETY  
STANDARD QUALITY ASSURANCE PROCEDURES  
INTOXILYZER MODEL 8000  
STANDARD CALIBRATION CHECK PROCEDURE

QA SPECIALIST J. Murray AGENCY Yavapai County Sheriff's Office  
DATE 12/23/2014 TIME 10:42:42  
INTOXILYZER SERIAL # 80-006228 LOCATION DUI Truck

1. Ensure that gas tank is attached to instrument and contains a standard alcohol concentration solution .100 AC.  
OR  
Pour a standard alcohol concentration solution        AC, into a clean dry simulator and assemble the simulator. Ensure that a tight seal has been made. Turn on the simulator and allow temperature to reach  $34^{\circ}\text{C} \pm 0.2^{\circ}\text{C}$
2. Intoxilyzer 8000 display reads "PUSH BUTTON TO START"
3. Go to the "Control Testing Menu". Select "D" for dry control test or "W" for wet control test. After selection is made press ENTER.
4. Air blank completed.
5. Calibration check completed. Test results 0. 096 AC.
6. Air blank completed.
7. Remove printed record. Attach the record to the completed checklist.

SIGNATURE J. Murray NS6/5597

INTOXILYZER 8000  
Location: YAVAPAI COUNTY SO  
Serial Number: 80-006228  
Core Version: 8105.48  
12/23/2014 10:42:42

-----  
Standard Lot#: 12  
Last Changed By: 1. 1 #1  
-----

QAS: J. MURRAY5597 #N56  
YAVAPAI SHERIFFS  
-----

Operator: J. MURRAY5597 #N56  
YAVAPAI SHERIFFS  
-----

Test	g/210L	Time
Air Blank	0.000	10:43:19
0.100 Cal Check	0.096	10:43:40
Air Blank	0.000	10:44:09

Preventative Maintenance  
Performed -- Timer Reset