

COPY

EXHIBIT G-2
THIS REPORT PREPARED PURSUANT TO DUTY IMPOSED BY A.A.C. R13-10-104 (A)

ARIZONA DEPARTMENT OF PUBLIC SAFETY
STANDARD QUALITY ASSURANCE PROCEDURES
INTOXILYZER MODEL 8000
STANDARD CALIBRATION CHECK PROCEDURE

QA SPECIALIST Scott E Bickley Sr. AGENCY Yavapai County SD
DATE 5-26-16 TIME 07:03:46
INTOXILYZER SERIAL # 80-006228 LOCATION DUI Van

- () 1. Ensure that gas tank is attached to instrument and contains a standard alcohol concentration solution 0.100 AC.
OR
Pour a standard alcohol concentration solution — AC, into a clean dry simulator and assemble the simulator. Ensure that a tight seal has been made. Turn on the simulator and allow temperature to reach $34^{\circ}\text{C} \pm 0.2^{\circ}\text{C}$
- () 2. Intoxilyzer 8000 display reads "PUSH BUTTON TO START"
- () 3. Go to the "Control Testing Menu". Select "D" for dry control test or "W" for wet control test. After selection is made press ENTER.
- () 4. Air blank completed.
- () 5. Calibration check completed. Test results 0.099 AC.
- () 6. Air blank completed.
- () 7. Remove printed record. Attach the record to the completed checklist.

SIGNATURE Scott E Bickley Sr.

DPS Form Exh G-2 (Rev 05-01)

Serial # 28515100A1
Exp 11-5-2017
Prrs 10700 psi

INTOXILYZER 8000

Location: YAVAPAI COUNTY SO

Serial Number: 80-006228

Core Version: 8105.48

05/26/2016

07:00:50

Standard Lot#: 01314100A3

Last Changed By: T. TEST #0000

QAS: S. BLAKLEY SR #L6861

YAVAPAI CSO

Operator: S. BLAKLEY SR #L6861

YAVAPAI CSO

Test	g/210L	Time
Air Blank	0.000	07:01:24
0.100 Cal Check	0.100	07:01:45
Air Blank	0.000	07:02:15

Preventative Maintenance
Performed -- Timer Reset

INTOXILYZER 8000

Location: YAVAPAI COUNTY SO

Serial Number: 80-006228

Core Version: 8105.48

05/26/2016

07:03:14

Standard Lot#: 01314100A3

Last Changed By: T. TEST #0000

QAS: S. BLAKLEY SR #L6861

YAVAPAI CSO

Operator: S. BLAKLEY SR #L6861

YAVAPAI CSO

Test	g/210L	Time
Air Blank	0.000	07:03:46
0.100 Cal Check	0.099	07:04:07
Air Blank	0.000	07:04:36

Preventative Maintenance
Performed -- Timer Reset