

EXHIBIT G-2
THIS REPORT PREPARED PURSUANT TO DUTY IMPOSED BY A.A.C. R13-10-104 (A)

ARIZONA DEPARTMENT OF PUBLIC SAFETY
STANDARD QUALITY ASSURANCE PROCEDURES
INTOXILYZER MODEL 8000
STANDARD CALIBRATION CHECK PROCEDURE

QA SPECIALIST M. MEDINA AGENCY PRESCOTT
DATE 9/23/14 TIME 0907
INTOXILYZER SERIAL # 80-006230 LOCATION PRESCOTT P.D.

- () 1. Ensure that gas tank is attached to instrument and contains a standard alcohol concentration solution 100 AC.
OR
Pour a standard alcohol concentration solution _____ AC, into a clean dry simulator and assemble the simulator. Ensure that a tight seal has been made. Turn on the simulator and allow temperature to reach 34°C ± 0.2°C
- () 2. Intoxilyzer 8000 display reads "PUSH BUTTON TO START" CONTACT QAS 31 DAYMAIT.
- () 3. Go to the "Control Testing Menu". Select "D" for dry control test or "W" for wet control test. After selection is made press ENTER.
- () 4. Air blank completed.
- () 5. Calibration check completed. Test results 0. 099 AC.
- () 6. Air blank completed.
- () 7. Remove printed record. Attach the record to the completed checklist.

SIGNATURE [Signature]

PRESCOTT PD CHANGE 5-12-06
2 11 251 R127

INTOXILYZER 8000

Location: PRESCOTT PD

Serial Number: 80-006230

Core Version: 8105.48

09/23/2014

08:07:24

Standard Lot#: 1234567890

Last Changed By: A. A #A

QAS: M. MEDINA #356

PRESCOTT PD

Operator: M. MEDINA #356

PRESCOTT PD

Test	g/210L	Time
Air Blank	0.000	08:07:59
0.100 Cal Check	0.099	08:08:19
Air Blank	0.000	08:08:49

Preventative Maintenance
Performed -- Timer Reset

EXHIBIT G-4
THIS REPORT PREPARED PURSUANT TO DUTY IMPOSED BY A.A.C. R13-10-104 (A)

ARIZONA DEPARTMENT OF PUBLIC SAFETY

STANDARD QUALITY ASSURANCE PROCEDURES
INTOXILYZER MODEL 8000

STANDARD QUALITY ASSURANCE PROCEDURE

QA SPECIALIST M. MEDINA AGENCY PRESCOTT

DATE 9/23/14 TIME 0811

INTOXILYZER SERIAL # 80-006230 LOCATION PRESCOTT P.D.

() 1. Display Reads "PUSH BUTTON TO START"

DIAGNOSTIC TESTS

() 1. Clock time check.

() 2. Date check.

OPERATIONAL TESTS

() 1. Alcohol-free subject test result 0. 000 AC.

() 2. Error recognition logic system functioning.

Not a Successfully Completed Test Sequence printed

() 3. Proper sample recognition system.

Not a Successfully Completed Test Sequence printed

Deficient sample printed.

() 4. Standard Calibration Check standard 0. 100 AC. Result 0. 098 AC.

Instrument is operating properly and accurately. YES NO

COMMENTS

SIGNATURE 

9-23-14
M

INTOXILYZER 8000
Location: PRESCOTT PD
Serial Number: 80-006230
Core Version: 8105.48
09/23/2014 08:11:29

Standard Lot#: 1234567890
Last Changed By: A. A #A

QAS: M. MEDINA #356
PRESCOTT PD

Operator: M. MEDINA #356
PRESCOTT PD

Subject: TEST, SQAP
DOB: 02/02/2000
Sex: M Weight: 180

15 Minute Deprivation Period? Yes

Test	g/210L	Time
Air Blank	ABT*	08:12:12
Air Blank	RFI**	08:12:33

*Sequence Aborted
**RFI Detect

Not a Successfully Completed Test Sequence

M. MEDINA
9/23/14

INTOXILYZER 8000
Location: PRESCOTT PD
Serial Number: 80-006230
Core Version: 8105.48
09/23/2014 08:12:56

Standard Lot#: 1234567890
Last Changed By: A. A #A

QAS: M. MEDINA #356
PRESCOTT PD

Operator: M. MEDINA #356
PRESCOTT PD

Subject: TEST, SQAP
DOB: 02/02/2000
Sex: M Weight: 180

15 Minute Deprivation Period? Yes

Test	g/210L	Time
Air Blank	0.000	08:13:57
Diagnostic Test	Pass	08:14:32
Air Blank	0.000	08:14:59
0.100 Cal Check	0.098	08:15:20
Air Blank	0.000	08:15:49
Subject Test	0.000	08:16:08
Air Blank	0.000	08:16:35
Wait	ABT*	
Air Blank	RFI**	08:16:58

*Sequence Aborted

**RFI Detect

Not a Successfully Completed Test Sequence

M. MEDINA
9/23/14

INTOXILYZER 8000
Location: PRESCOTT PD
Serial Number: 80-006230
Core Version: 8105.48
09/23/2014 08:17:27

Standard Lot#: 1234567890
Last Changed By: A. A #A

QAS: M. MEDINA #356
PRESCOTT PD

Operator: M. MEDINA #356
PRESCOTT PD

Subject: TEST, SQAP
DOB: 02/02/2000
Sex: M Weight: 180

15 Minute Deprivation Period? Yes

Test	g/210L	Time
Air Blank	IPS*	08:18:20
Air Blank	0.000	08:18:49

*Improper Sample

Not a Successfully Completed Test Sequence

M. MEDINA
9/23/14

INTOXILYZER 8000
Location: PRESCOTT PD
Serial Number: 80-006230
Core Version: 8105.48
09/23/2014 08:19:12

Standard Lot#: 1234567890
Last Changed By: A. A #A

QAS: M. MEDINA #356
PRESCOTT PD

Operator: M. MEDINA #356
PRESCOTT PD

Subject: TEST, SQAP
DOB: 02/02/2000
Sex: M Weight: 180

15 Minute Deprivation Period? Yes

Test	g/210L	Time
Air Blank	0.000	08:20:05
Diagnostic Test	Pass	08:20:40
Air Blank	0.000	08:21:08
0.100 Cal Check	0.099	08:21:28
Air Blank	0.000	08:21:58
Subject Test	0.022*	08:25:05
Air Blank	0.000	08:25:35
Air Blank	0.000	08:26:02
Subject Test	XXX**	08:26:27
Air Blank	0.000	08:26:59

*Deficient Sample

**Invalid Sample

Begin new deprivation period

Not a Successfully Completed Test Sequence

EXHIBIT G-2
THIS REPORT PREPARED PURSUANT TO DUTY IMPOSED BY A.A.C. R13-10-104 (A)

ARIZONA DEPARTMENT OF PUBLIC SAFETY
STANDARD QUALITY ASSURANCE PROCEDURES
INTOXILYZER MODEL 8000
STANDARD CALIBRATION CHECK PROCEDURE

QA SPECIALIST K. MEDINA AGENCY PRESCOTT

DATE 9/23/14 TIME _____

INTOXILYZER SERIAL # 80-006230 LOCATION PRESCOTT P.D.

- () 1. Ensure that gas tank is attached to instrument and contains a standard alcohol concentration solution 100 AC.
OR
Pour a standard alcohol concentration solution _____ AC, into a clean dry simulator and assemble the simulator. Ensure that a tight seal has been made. Turn on the simulator and allow temperature to reach $34^{\circ}\text{C} \pm 0.2^{\circ}\text{C}$.
- () 2. Intoxilyzer 8000 display reads "PUSH BUTTON TO START"
- () 3. Go to the "Control Testing Menu". Select "D" for dry control test or "W" for wet control test. After selection is made press ENTER.
- () 4. Air blank completed.
- () 5. Calibration check completed. Test results 0. 098 AC.
- () 6. Air blank completed.
- () 7. Remove printed record. Attach the record to the completed checklist.

SIGNATURE 

PRESCOTT P.D CHANGE 5-18-06
NO. 177

M. MEDINA
9/23/14

INTOXILYZER 8000
Location: PRESCOTT PD
Serial Number: 80-006230
Core Version: 8105.48
09/23/2014 08:27:42

Standard Lot#: 1234567890
Last Changed By: A. A #A

QAS: M. MEDINA #356
PRESCOTT PD

Operator: M. MEDINA #356
PRESCOTT PD

Test	g/210L	Time
Air Blank	0.000	08:28:13
0.100 Cal Check	0.098	08:28:33
Air Blank	0.000	08:29:03

Preventative Maintenance
Performed -- Timer Reset