

SMA
12/12/14

EXHIBIT G-2
THIS REPORT PREPARED PURSUANT TO DUTY IMPOSED BY A.A.C. R13-10-104 (A)

ARIZONA DEPARTMENT OF PUBLIC SAFETY
STANDARD QUALITY ASSURANCE PROCEDURES
INTOXILYZER MODEL 8000

STANDARD CALIBRATION CHECK PROCEDURE

QA SPECIALIST S. Helton AGENCY Prescott P.D.
DATE 12/12/2014 TIME 1217
INTOXILYZER SERIAL # 80-006230 LOCATION P.P.D.

- 1. Ensure that gas tank is attached to instrument and contains a standard alcohol concentration solution .100 AC.
OR
~~Pour a standard alcohol concentration solution _____ AC, insert in dry simulator and assemble the simulator. Ensure that a tight seal has been made. Turn on the simulator and allow temperature to reach 34°C ± 0.2°C~~
- 2. Intoxilyzer 8000 display reads "PUSH BUTTON TO START"
- 3. Go to the "Control Testing Menu". Select "D" for dry control test or "W" for wet control test. After selection is made press ENTER.
- 4. Air blank completed.
- 5. Calibration check completed. Test results 0. 199 AC.
- 6. Air blank completed.
- 7. Remove printed record. Attach the record to the completed checklist.

SIGNATURE S. Helton 372

SMA 372
12/12/14

INTOXILYZER 8000
Location: PRESCOTT PD
Serial Number: 80-006230
Core Version: 8105.48
12/12/2014 12:23:00

Standard Lot#: 1234567890
Last Changed By: A. A #A

QAS: S. HILTON #372
PRESCOTT PD

Operator: S. HILTON #372
PRESCOTT PD

Test	g/210L	Time
Air Blank	0.000	12:23:41
0.100 Cal Check	0.099	12:24:02
Air Blank	0.000	12:24:31

Preventative Maintenance
Performed -- Timer Reset

SMA 372
12/12/14

EXHIBIT G-4
THIS REPORT PREPARED PURSUANT TO DUTY IMPOSED BY A.A.C. R13-10-104 (A)

ARIZONA DEPARTMENT OF PUBLIC SAFETY

STANDARD QUALITY ASSURANCE PROCEDURES
INTOXILYZER MODEL 8000

STANDARD QUALITY ASSURANCE PROCEDURE

QA SPECIALIST S. Hutton AGENCY Prescott P.D.
DATE 12/12/14 TIME 1235
INTOXILYZER SERIAL # 80-006230 LOCATION P.P.D.

~~SMA~~ 1. Display Reads "PUSH BUTTON TO START"

DIAGNOSTIC TESTS

- ~~SMA~~ 1. Clock time check.
- ~~SAT~~ 2. Date check.

OPERATIONAL TESTS

- ~~SMA~~ 1. Alcohol-free subject test result 0. 000 AC.
- ~~SMA~~ 2. Error recognition logic system functioning.
Not a Successfully Completed Test Sequence printed
- ~~SMA~~ 3. Proper sample recognition system.
Not a Successfully Completed Test Sequence printed
Deficient sample printed.
- ~~SMA~~ 4. Standard Calibration Check standard 0. 100 AC. Result 0. 099 AC.

Instrument is operating properly and accurately. YES NO

COMMENTS OPERATION PROPERLY

SIGNATURE S. Hutton

Prescott PD CALLER 12:06

SMA
12/10/14

INTOXILYZER 8000
Location: PRESCOTT PD
Serial Number: 80-006230
Core Version: 8105.48
12/12/2014 12:28:41

Standard Lot#: 1234567890
Last Changed By: A. A #A

QAS: S. HILTON #372
PRESCOTT PD

Operator: S. HILTON #372
PRESCOTT PD

Subject: SQAP,TEST
DOB: 01/01/2014
Sex: M Weight: 175

15 Minute Deprivation Period? Yes

Test	g/210L	Time
Air Blank	0.000	12:30:14
Diagnostic Test	Pass	12:30:47
Air Blank	0.000	12:31:15
0.100 Cal Check	0.099	12:31:36
Air Blank	0.000	12:32:05
Subject Test	0.000	12:32:28
Air Blank	0.000	12:32:56
Wait	ABT*	
Air Blank	0.000	12:33:30

*Sequence Aborted

Not a Successfully Completed Test Sequence

SWA
12/12/14

INTOXILYZER 8000
Location: PRESCOTT PD
Serial Number: 80-006230
Core Version: 8105.48
12/12/2014 12:34:13

Standard Lot#: 1234567890
Last Changed By: A. A #A

QAS: S. HILTON #372
PRESCOTT PD

Operator: S. HILTON #372
PRESCOTT PD

Subject: SQAP, TEST
DOB: 01/01/2014
Sex: M Weight: 175

15 Minute Deprivation Period? Yes

Test	g/210L	Time
Air Blank	IPS*	12:35:00
Air Blank	0.000	12:35:28

*Improper Sample

Not a Successfully Completed Test Sequence

SMA
12/12/14

INTOXILYZER 8000
Location: PRESCOTT PD
Serial Number: 80-006230
Core Version: 8105.48
12/12/2014 12:36:40

Standard Lot#: 1234567890
Last Changed By: A. A #A

QAS: S. HILTON #372
PRESCOTT PD

Operator: S. HILTON #372
PRESCOTT PD

Subject: SQAP, TEST
DOB: 01/01/2014
Sex: M Weight: 175

15 Minute Deprivation Period? Yes

Test	g/210L	Time
Air Blank	RFI*	12:37:37
Air Blank	0.000	12:38:05

*RFI Detect

Not a Successfully Completed Test Sequence

SMA
12/12/14

INTOXILYZER 8000
Location: PRESCOTT PD
Serial Number: 80-006230
Core Version: 8105.48
12/12/2014 12:42:56

Standard Lot#: 1234567890
Last Changed By: A. A #A

QAS: S. HILTON #372
PRESCOTT PD

Operator: S. HILTON #372
PRESCOTT PD

Subject: SQAP, TEST
DOB: 01/01/2014
Sex: M Weight: 175

15 Minute Deprivation Period? Yes

Test	g/210L	Time
Air Blank	0.000	12:44:14
Diagnostic Test	Pass	12:44:48
Air Blank	0.000	12:45:15
0.100 Cal Check	0.098	12:45:36
Air Blank	0.000	12:46:05
Subject Test	0.118*	12:49:13
Air Blank	0.000	12:49:46
Air Blank	0.000	12:50:14
Subject Test	XXX**	12:50:38
Air Blank	0.000	12:51:13

*Deficient Sample
**Invalid Sample
Begin new deprivation period

Not a Successfully Completed Test Sequence

sent 12/12/14

EXHIBIT G-2
THIS REPORT PREPARED PURSUANT TO DUTY IMPOSED BY A.A.C. R13-10-104 (A)
ARIZONA DEPARTMENT OF PUBLIC SAFETY
STANDARD QUALITY ASSURANCE PROCEDURES
INTOXILYZER MODEL 8000
STANDARD CALIBRATION CHECK PROCEDURE

QA SPECIALIST S. Helt AGENCY Prescott P.P.
DATE 12/12/14 TIME 1253
INTOXILYZER SERIAL # 80-006230 LOCATION P.P.D.

- (H) 1. Ensure that gas tank is attached to instrument and contains a standard alcohol concentration solution .100 AC.
OR
Pour a standard alcohol concentration solution _____ AC, into a clean dry simulator and assemble the simulator. ~~Ensure that a tight seal has been made. Turn on the simulator and allow temperature to reach 34°C ± 0.2°C~~
- (H) 2. Intoxilyzer 8000 display reads "PUSH BUTTON TO START"
- (H) 3. Go to the "Control Testing Menu". Select "D" for dry control test or "W" for wet control test. After selection is made press ENTER.
- (H) 4. Air blank completed.
- (H) 5. Calibration check completed. Test results 0.097 AC.
- (H) 6. Air blank completed.
- (H) 7. Remove printed record. Attach the record to the completed checklist.

SIGNATURE S. Helt 372

Prescott PD CHANGE 5-12-06

SMA
12/12/14

INTOXILYZER 8000
Location: PRESCOTT PD
Serial Number: 80-006230
Core Version: 8105.48
12/12/2014 12:54:57

Standard Lot#: 1234567890
Last Changed By: A. A #A

QAS: S. HILTON #372
PRESCOTT PD

Operator: S. HILTON #372
PRESCOTT PD

Test	g/210L	Time
Air Blank	0.000	12:55:31
0.100 Cal Check	0.097	12:55:52
Air Blank	0.000	12:56:21

Preventative Maintenance
Performed -- Timer Reset