

EXHIBIT G-2
THIS REPORT PREPARED PURSUANT TO DUTY IMPOSED BY A.A.C. R13-10-104 (A)

ARIZONA DEPARTMENT OF PUBLIC SAFETY
STANDARD QUALITY ASSURANCE PROCEDURES
INTOXILYZER MODEL 8000

STANDARD CALIBRATION CHECK PROCEDURE

QA SPECIALIST S. HILTON AGENCY Prescott P.D.
DATE 6/13/15 TIME 1045
INTOXILYZER SERIAL # 80-006230 LOCATION DUT Room

- SH* 1. Ensure that gas tank is attached to instrument and contains a standard alcohol concentration solution .100 AC.
OR
Pour a standard alcohol concentration solution .100 AC, into a clean dry simulator and assemble the simulator. Ensure that a tight seal has been made. Turn on the simulator and allow temperature to reach 32°C ± 0.2°C
- SH* 2. Intoxilyzer 8000 display reads "PUSH BUTTON TO START"
- SH* 3. Go to the "Control Testing Menu". Select "D" for dry control test or "W" for wet control test. After selection is made press ENTER.
- SH* 4. Air blank completed.
- SH* 5. Calibration check completed. Test results 0. 097 AC.
- SH* 6. Air blank completed.
- SH* 7. Remove printed record. Attach the record to the completed checklist.

SIGNATURE *S. Hilton*

INTOXILYZER 8000

Location: PRESCOTT PD

Serial Number: 80-006230

Core Version: 8105.48

06/13/2015

10:45:08

Standard Lot#: 1234567890

Last Changed By: A. A #A

QAS: S. HILTON #372

PRESCOTT PD

Operator: S. HILTON #372

PRESCOTT PD

Test	g/210L	Time
Air Blank	0.000	10:45:39
0.100 Cal Check	0.097	10:46:00
Air Blank	0.000	10:46:30

Preventative Maintenance
Performed -- Timer Reset

SMA 372

EXHIBIT G-4
THIS REPORT PREPARED PURSUANT TO DUTY IMPOSED BY A.A.C. R13-10-104 (A)

ARIZONA DEPARTMENT OF PUBLIC SAFETY
STANDARD QUALITY ASSURANCE PROCEDURES
INTOXILYZER MODEL 8000
STANDARD QUALITY ASSURANCE PROCEDURE

QA SPECIALIST S. Hutton AGENCY PRESCOTT P.P.
DATE 6/13/15 TIME 1049
INTOXILYZER SERIAL # 80-006230 LOCATION DUI ROOM

Set 1. Display Reads "PUSH BUTTON TO START"

DIAGNOSTIC TESTS

- Set 1. Clock time check.
- Set 2. Date check.

OPERATIONAL TESTS

- Set 1. Alcohol-free subject test result 0. 100 AC.
- Set 2. Error recognition logic system; functioning.
Not a Successfully Completed Test Sequence printed
- Set 3. Proper sample recognition system.
Not a Successfully Completed Test Sequence printed
Deficient sample printed.
- Set 4. Standard Calibration Check standard 0. 100 AC. Result 0. .098 AC.

Instrument is operating properly and accurately. YES NO

COMMENTS

- FUNCTIONING PROPERLY
- PLACED IN SERVICE

SIGNATURE S. Hutton

PRESCOTT PD CAROLE S 12:06

INTOXILYZER 8000
Location: PRESCOTT PD
Serial Number: 80-006230
Core Version: 8105.48
06/13/2015

10:49:11

Standard Lot#: 1234567890
Last Changed By: A. A

#A

QAS: S. HILTON
PRESCOTT PD

#372

Operator: S. HILTON
PRESCOTT PD

#372


Subject: SQAP, TEST
DOB: 12/19/1973
Sex: M

Weight: 175

15 Minute Deprivation Period? Yes

Test	g/210L	Time
Air Blank	0.000	10:50:22
Diagnostic Test	Pass	10:50:57
Air Blank	0.000	10:51:25
0.100 Cal Check	0.097	10:51:45
Air Blank	0.000	10:52:15
Subject Test	0.000	10:52:38
Air Blank	0.000	10:53:05
Wait	ABT*	
Air Blank	0.000	10:53:40

*Sequence Aborted



Not a Successfully Completed Test Sequence

INTOXILYZER 8000
Location: PRESCOTT PD
Serial Number: 80-006230
Core Version: 8105.48
06/13/2015 10:54:11

Standard Lot#: 1234567890
Last Changed By: A. A #A

QAS: S. HILTON #372
PRESCOTT PD

Operator: S. HILTON #372
PRESCOTT PD

Subject: SQAP, TEST
DOB: 12/19/1973
Sex: M Weight: 175

15 Minute Deprivation Period? Yes

Test	g/210L	Time
Air Blank	IPS*	10:55:02
Air Blank	0.000	10:55:30

*Improper Sample

SIM 32

Not a Successfully Completed Test Sequence

INTOXILYZER 8000
Location: PRESCOTT PD
Serial Number: 80-006230
Core Version: 8105.48
06/13/2015 10:55:59

Standard Lot#: 1234567890
Last Changed By: A. A #A

QAS: S. HILTON #372
PRESCOTT PD

Operator: S. HILTON #372
PRESCOTT PD

Subject: SQAP, TEST
DOB: 12/19/1973
Sex: M Weight: 175

15 Minute Deprivation Period? Yes

Test	g/210L	Time
Air Blank	RFI*	10:56:41
Air Blank	0.000	10:57:10

*RFI Detect

SM 3M

Not a Successfully Completed Test Sequence

INTOXILYZER 8000
Location: PRESCOTT PD
Serial Number: 80-006230
Core Version: 8105.48
06/13/2015

10:57:38

Standard Lot#: 1234567890
Last Changed By: A. A #A

QAS: S. HILTON #372
PRESCOTT PD

Operator: S. HILTON #372
PRESCOTT PD

Subject: SQAP,TEST
DOB: 12/19/1973
Sex: M Weight: 175


15 Minute Deprivation Period? Yes

Test	g/210L	Time
Air Blank	0.000	10:58:33
Diagnostic Test	Pass	10:59:08
Air Blank	0.000	10:59:35
0.100 Cal Check	0.098	10:59:56
Air Blank	0.000	11:00:25
Subject Test	0.171*	11:03:33
Air Blank	0.000	11:04:07
Air Blank	0.000	11:04:35
Subject Test	XXX**	11:05:00
Air Blank	0.000	11:05:36

*Deficient Sample

**Invalid Sample

Begin new deprivation period



Not a Successfully Completed Test Sequence

EXHIBIT G-2
THIS REPORT PREPARED PURSUANT TO DUTY IMPOSED BY A.A.C. R13-10-104 (A)

ARIZONA DEPARTMENT OF PUBLIC SAFETY
STANDARD QUALITY ASSURANCE PROCEDURES
INTOXILYZER MODEL 8000
STANDARD CALIBRATION CHECK PROCEDURE

QA SPECIALIST S. Hinton AGENCY Prescott P.D.
DATE 6/13/15 TIME 1106
INTOXILYZER SERIAL # 80-006230 LOCATION D.E. Room

- (S) 1. Ensure that gas tank is attached to instrument and contains a standard alcohol concentration solution. 100 AC.
OR
Pour a standard alcohol concentration solution AC. into a clean dry simulator and assemble the simulator. Ensure that a tight seal has been made. Turn on the simulator and allow temperature to reach $34^{\circ}\text{C} \pm 0.2^{\circ}\text{C}$
- (S) 2. Intoxilyzer 8000 display reads "PUSH BUTTON TO START"
(S) 3. Go to the "Control Testing Menu". Select "D" for dry control test or "W" for wet control test. After selection is made press ENTER.
(S) 4. Air blank completed.
(S) 5. Calibration check completed. Test results 0. 098 AC.
(S) 6. Air blank completed.
(S) 7. Remove printed record. Attach the record to the completed checklist.

SIGNATURE S. Hinton

DPS Form Exh G-2 (Rev 05-01)

INTOXILYZER 8000

Location: PRESCOTT PD

Serial Number: 80-006230

Core Version: 8105.48

06/13/2015

11:06:49

Standard Lot#: 1234567890

Last Changed By: A. A #A

QAS: S. HILTON #372


PRESCOTT PD

Operator: S. HILTON #372

PRESCOTT PD

Test	g/210L	Time
Air Blank	0.000	11:07:21
0.100 Cal Check	0.098	11:07:42
Air Blank	0.000	11:08:11

Preventative Maintenance
Performed -- Timer Reset

A handwritten signature in black ink, appearing to be 'SA' with a long horizontal stroke extending to the right.