


EXHIBIT G-2
THIS REPORT PREPARED PURSUANT TO DUTY IMPOSED BY A.A.C. R13-10-104 (A)

ARIZONA DEPARTMENT OF PUBLIC SAFETY
STANDARD QUALITY ASSURANCE PROCEDURES
INTOXILYZER MODEL 8000
STANDARD CALIBRATION CHECK PROCEDURE

QA SPECIALIST D. BARDON AGENCY PRESCOTT P.D.
DATE 10.11.18 TIME 1344
INTOXILYZER SERIAL # 80-006230 LOCATION PRESCOTT P.D.

1. Ensure that gas tank is attached to instrument and contains a standard alcohol concentration solution 0.100 AC.
OR
Pour a standard alcohol concentration solution _____ AC, into a clean dry simulator and assemble the simulator. Ensure that a tight seal has been made. Turn on the simulator and allow temperature to reach $34^{\circ}\text{C} \pm 0.2^{\circ}\text{C}$
2. Intoxilyzer 8000 display reads "PUSH BUTTON TO START"
3. Go to the "Control Testing Menu". Select "D" for dry control test or "W" for wet control test. After selection is made press ENTER.
4. Air blank completed.
5. Calibration check completed. Test results 0. 092 AC.
6. Air blank completed.
7. Remove printed record. Attach the record to the completed checklist.

SIGNATURE 

INTOXILYZER 8000
Location: PRESCOTT PD
Serial Number: 80-006230
Core Version: 8105.48
10/11/2018 13:44:25

Standard Lot#: 27016100A1
Last Changed By: S. HILTON #372

QAS: D. BARDON #351
PRESCOTT PD

Operator: D. BARDON #351
PRESCOTT PD

Test g/210L Time

Air Blank 0.000 13:44:56
0.100 Cal Check 0.092 13:45:17
Air Blank 0.000 13:45:46

Preventative Maintenance
Performed -- Timer Reset

EXHIBIT G-4
THIS REPORT PREPARED PURSUANT TO DUTY IMPOSED BY A.A.C. R13-10-104 (A)

ARIZONA DEPARTMENT OF PUBLIC SAFETY
STANDARD QUALITY ASSURANCE PROCEDURES
INTOXILYZER MODEL 8000
STANDARD QUALITY ASSURANCE PROCEDURE

QA SPECIALIST D. BARDON AGENCY PRESCOTT P.D.
DATE 10.11.18 TIME 1319
INTOXILYZER SERIAL # 80-006230 LOCATION PRESCOTT P.D.

() 1. Display Reads "PUSH BUTTON TO START"

DIAGNOSTIC TESTS


- () 1. Clock time check.
- () 2. Date check.

OPERATIONAL TESTS

- () 1. Alcohol-free subject test result 0. 000 AC.
- () 2. Error recognition logic system functioning.
Not a Successfully Completed Test Sequence printed
- () 3. Proper sample recognition system.
Not a Successfully Completed Test Sequence printed
Deficient sample printed.
- () 4. Standard Calibration Check standard 0. 100 AC. Result 0. 093 AC.

Instrument is operating properly and accurately. YES NO

COMMENTS

SIGNATURE  #351

INTOXILYZER 8000
Location: PRESCOTT PD
Serial Number: 80-006230
Core Version: 8105.48
10/11/2018

13:49:32

Standard Lot#: 27016100A1
Last Changed By: S. HILTON #372

QAS: D. BARDON #351
PRESCOTT PD

Operator: D. BARDON #351
PRESCOTT PD

Subject: SQAP, TEST
DOB: 09/09/1999
Sex: M Weight: 210

15 Minute Deprivation Period? Yes

Test	g/210L	Time
Air Blank	0.000	13:50:32
Diagnostic Test	Pass	13:51:07
Air Blank	0.000	13:51:34
0.100 Cal Check	0.092	13:51:55
Air Blank	0.000	13:52:24
Subject Test	0.000	13:52:51
Air Blank	0.000	13:53:18
Wait	ABT*	
Air Blank	0.000	13:53:52

*Sequence Aborted

Not a Successfully Completed Test Sequence

INTOXILYZER 8000
Location: PRESCOTT PD
Serial Number: 80-006230
Core Version: 8105.48
10/11/2018

13:54:27

Standard Lot#: 27016100A1
Last Changed By: S. HILTON #372

QAS: D. BARDON #351
PRESCOTT PD

Operator: D. BARDON #351
PRESCOTT PD

Subject: SQAP,TEST
DOB: 09/09/1999
Sex: M Weight: 210

15 Minute Deprivation Period? Yes

Test	g/210L	Time
Air Blank	IPS*	13:55:05
Air Blank	0.000	13:55:34

*Improper Sample

Not a Successfully Completed Test Sequence

INTOXILYZER 8000
Location: PRESCOTT PD
Serial Number: 80-006230
Core Version: 8105.48
10/11/2018

13:56:01

Standard Lot#: 27016100A1
Last Changed By: S. HILTON #372

QAS: D. BARDON #351
PRESCOTT PD

Operator: D. BARDON #351
PRESCOTT PD

Subject: SQAP, TEST
DOB: 09/09/1999
Sex: M Weight: 210

15 Minute Deprivation Period? Yes

Test	g/210L	Time
Air Blank	RFI*	13:56:45
Air Blank	0.000	13:57:13

*RFI Detect

Not a Successfully Completed Test Sequence

INTOXILYZER 8000
Location: PRESCOTT PD
Serial Number: 80-006230
Core Version: 8105.48
10/11/2018

13:57:33

Standard Lot#: 27016100A1
Last Changed By: S. HILTON #372

QAS: D. BARDON #351
PRESCOTT PD

Operator: D. BARDON #351
PRESCOTT PD

Subject: SQAP,TEST
DOB: 09/09/1999
Sex: M Weight: 210

15 Minute Deprivation Period? Yes

Test	g/210L	Time
Air Blank	0.000	13:58:30
Diagnostic Test	Pass	13:59:05
Air Blank	0.000	13:59:32
0.100 Cal Check	0.093	13:59:53
Air Blank	0.000	14:00:22
Subject Test	0.028*	14:03:30
Air Blank	0.000	14:04:00
Air Blank	0.000	14:04:28
Subject Test	XXX**	14:04:54
Air Blank	0.000	14:05:26

*Deficient Sample

**Invalid Sample

Begin new deprivation period

Not a Successfully Completed Test Sequence

EXHIBIT G-2
THIS REPORT PREPARED PURSUANT TO DUTY IMPOSED BY A.A.C. R13-10-104 (A)

ARIZONA DEPARTMENT OF PUBLIC SAFETY
STANDARD QUALITY ASSURANCE PROCEDURES
INTOXILYZER MODEL 8000
STANDARD CALIBRATION CHECK PROCEDURE

QA SPECIALIST D. BARDON AGENCY PRESCOTT P.D.

DATE 10.11.18 TIME 1406

INTOXILYZER SERIAL # 80-006230 LOCATION PRESCOTT P.D.

1. Ensure that gas tank is attached to instrument and contains a standard alcohol concentration solution 100 AC.
OR
Pour a standard alcohol concentration solution _____ AC, into a clean dry simulator and assemble the simulator. Ensure that a tight seal has been made. Turn on the simulator and allow temperature to reach $34^{\circ}\text{C} \pm 0.2^{\circ}\text{C}$.
2. Intoxilyzer 8000 display reads "PUSH BUTTON TO START"
3. Go to the "Control Testing Menu". Select "D" for dry control test or "W" for wet control test. After selection is made press ENTER.
4. Air blank completed.
5. Calibration check completed. Test results 0. 093 AC.
6. Air blank completed.
7. Remove printed record. Attach the record to the completed checklist.

SIGNATURE [Signature]

INTOXILYZER 8000
Location: PRESCOTT PD
Serial Number: 80-006230
Core Version: 8105.48
10/11/2018

14:06:57

Standard Lot#: 27016100A1
Last Changed By: S. HILTON #372

QAS: D. BARDON #351
PRESCOTT PD

Operator: D. BARDON #351
PRESCOTT PD

Test	g/210L	Time
Air Blank	0.000	14:07:28
0.100 Cal Check	0.093	14:07:49
Air Blank	0.000	14:08:19

Preventative Maintenance
Performed -- Timer Reset

EXHIBIT G-2
THIS REPORT PREPARED PURSUANT TO DUTY IMPOSED BY A.A.C. R13-10-104 (A)

ARIZONA DEPARTMENT OF PUBLIC SAFETY
STANDARD QUALITY ASSURANCE PROCEDURES
INTOXILYZER MODEL 8000
STANDARD CALIBRATION CHECK PROCEDURE

QA SPECIALIST D. BARDON AGENCY PRESCOTT POLICE DEPT
DATE 10/18/18 TIME 1203
INTOXILYZER SERIAL # 60-006230 LOCATION PRESCOTT P.D

- (✓) 1. Ensure that gas tank is attached to instrument and contains a standard alcohol concentration solution .100 AC.
OR
Pour a standard alcohol concentration solution _____ AC, into a clean dry simulator and assemble the simulator. Ensure that a tight seal has been made. Turn on the simulator and allow temperature to reach $34^{\circ}\text{C} \pm 0.2^{\circ}\text{C}$
- (✓) 2. Intoxilyzer 8000 display reads "PUSH BUTTON TO START"
- (✓) 3. Go to the "Control Testing Menu". Select "D" for dry control test or "W" for wet control test. After selection is made press ENTER.
- (✓) 4. Air blank completed.
- (✓) 5. Calibration check completed. Test results 0. 095 AC.
- (✓) 6. Air blank completed.
- (✓) 7. Remove printed record. Attach the record to the completed checklist.

SIGNATURE [Signature] 351

INTOXILYZER 8000
Location: PRESCOTT PD
Serial Number: 80-006230
Core Version: 8105.48
10/18/2018

12:03:17

Standard Lot#: 22817100A3
Last Changed By: D. BARDON #351

QAS: D. BARDON #351
PRESCOTT PD

Operator: D. BARDON #351
PRESCOTT PD

Test	g/210L	Time
Air Blank	0.000	12:03:48
0.100 Cal Check	0.095	12:04:08
Air Blank	0.000	12:04:38

Preventative Maintenance
Performed -- Timer Reset