

EXHIBIT G-2
THIS REPORT PREPARED PURSUANT TO DUTY IMPOSED BY A.A.C. R-13-10-104 (A)

ARIZONA DEPARTMENT OF PUBLIC SAFETY
STANDARD QUALITY ASSURANCE PROCEDURES
INTOXILYZER MODEL 8000

STANDARD CALIBRATION CHECK PROCEDURE

QA SPECIALIST W. Lopez 30 AGENCY SEDONA P.D.

DATE 05/02/14 TIME 0410

INTOXILYZER SERIAL # 80-000924 LOCATION SEDONA P.D.

(bt) 1. Ensure that gas tank is attached to instrument and contains a standard alcohol concentration solution . 100 AC.

Pour a standard alcohol concentration solution . OR N/A AC, into a clean dry simulator and assemble the simulator. Ensure that a tight seal has been made. Turn on the simulator and allow temperature to reach 34°C ± 0.2°C

- (bt) 2. Intoxilyzer 8000 display reads "PUSH BUTTON TO START"
- (bt) 3. Go to the "Control Testing Menu". Select "D" for dry control test or "W" for wet control test. After selection is made press ENTER.
- (bt) 4. Air blank complete.
- (bt) 5. Calibration check complete. Test results 0. 096 AC.
- (bt) 6. Air blank complete.
- (bt) 7. Remove printed record. Attach the record to the complete checklist.

SIGNATURE W. Lopez 30

DPS Form Exb G-2 (Rev 05-01)

PASS

I hereby certify that the above and forgoing is a true and correct copy of the record of periodic maintenance for the intoxilyzer 8000 referenced therein maintained by the Sedona Police Department pursuant to the requirements of the Arizona Department of Public Safety.

Dated this 2 day of MAY, 20 14.

W. Lopez 30
Quality Assurance Specialist

INTOXILYZER 8000
 Location: SEDONA
 Serial Number: 80-000924
 Core Version: 8105.48
 05/02/2014 04:10:04

 Standard Lot#: 01513100A1
 Last Changed By: K. HUDSPETH #474

QAS: W. LOPEZ #30
 SEDONA PD

Operator: W. LOPEZ #30
 SEDONA PD

Test	g/210L	Time
Air Blank	0.000	04:10:37
0.100 Cal Check	0.096	04:10:58
Air Blank	0.000	04:11:28

Preventative Maintenance
 Performed -- Timer Reset

I hereby certify that the above and forgoing is a true and correct copy of the record of periodic maintenance for the intoxilyzer 8000 referenced therein maintained by the Sedona Police Department pursuant to the requirements of the Arizona Department of Public Safety.

Dated this 2 day of MAY, 2014.

Will 30
 Quality Assurance Specialist

EXHIBIT G-2

THIS REPORT PREPARED PURSUANT TO DUTY IMPOSED BY A.A.C. R-13-10-104 (A)

ARIZONA DEPARTMENT OF PUBLIC SAFETY

STANDARD QUALITY ASSURANCE PROCEDURES

INTOXILYZER MODEL 8000

STANDARD CALIBRATION CHECK PROCEDURE

QA SPECIALIST W. Lopez #30 AGENCY SEONA P.O.

DATE 05/02/14 TIME 0353

INTOXILYZER SERIAL # 80-000924 LOCATION SEONA P.O.

() 1. Ensure that gas tank is attached to instrument and contains a standard alcohol concentration solution . 100 AC.

Pour a standard alcohol concentration solution . ^{OR} N/A AC, into a clean dry simulator and assemble the simulator. Ensure that a tight seal has been made. Turn on the simulator and allow temperature to reach $34^{\circ}\text{C} \pm 0.2^{\circ}\text{C}$

() 2. Intoxilyzer 8000 display reads "PUSH BUTTON TO START"

() 3. Go to the "Control Testing Menu". Select "D" for dry control test or "W" for wet control test. After selection is made press ENTER.

() 4. Air blank complete.

() 5. Calibration check complete. Test results 0. _____ AC.

() 6. Air blank complete.

() 7. Remove printed record. Attach the record to the complete checklist.

SIGNATURE _____

DPS Form Exb G-2 (Rev 05-01)

DSP TEST FAILED TWICE
ANALYTICAL TEST FAILED TWICE

TURND OFF INTOX. 0405 HOURS AND RE-STARTED.

BOTH TEST PASSED