

BEFORE SQAP

EXHIBIT G-2
THIS REPORT PREPARED PURSUANT TO DUTY IMPOSED BY A.A.C. R-13-10-104 (A)

ARIZONA DEPARTMENT OF PUBLIC SAFETY

STANDARD QUALITY ASSURANCE PROCEDURES
INTOXILYZER MODEL 8000

STANDARD CALIBRATION CHECK PROCEDURE

QA SPECIALIST KEVIN HUDSPETH #474 AGENCY SEDONA P.D.

DATE 12/28/2014 TIME 0908

INTOXILYZER SERIAL # 80-000924 LOCATION SEDONA P.D.

(KH) 1. Ensure that gas tank is attached to instrument and contains a standard alcohol concentration solution . 100 AC.

OR
Pour a standard alcohol concentration solution . N/A AC, into a clean dry simulator and assemble the simulator. Ensure that a tight seal has been made. Turn on the simulator and allow temperature to reach 34°C ± 0.2°C

(KH) 2. Intoxilyzer 8000 display reads "PUSH BUTTON TO START"

(KH) 3. Go to the "Control Testing Menu". Select "D" for dry control test or "W" for wet control test. After selection is made press ENTER.


(KH) 4. Air blank complete.

(KH) 5. Calibration check complete. Test results 0. 097 AC.

(KH) 6. Air blank complete.

(KH) 7. Remove printed record. Attach the record to the complete checklist.


SIGNATURE

 474 (PASSED)

DPS Form Exb G-2 (Rev 05-01)

I hereby certify that the above and foregoing is a true and correct copy of the record of periodic maintenance for the Intoxilyzer 8000 referenced therein maintained by the Sedona Police Department pursuant to the requirements of the Arizona Department of Public Safety.

Date this 28TH day of DECEMBER 20 14

 474
Quality Assurance Specialist

INTOXILYZER 8000
Location: SEDONA
Serial Number: 80-000924
Core Version: 8105.48
12/28/2014

09:08:20

Standard Lot#: 01513100A1
Last Changed By: K. HUDSPETH #474

QAS: K. HUDSPETH #474
SEDONA PD

Operator: K. HUDSPETH #474
SEDONA PD

Test	g/210L	Time
Air Blank	0.000	09:08:52
0.100 Cal Check	0.097	09:09:13
Air Blank	0.000	09:09:42

[Signature] 474 PASSED

Preventative Maintenance
Performed -- Timer Reset

I hereby certify that the above and foregoing is a true and correct copy of the record of periodic maintenance for the intoxilyzer 8000 referenced therein maintained by the Sedona Police Department pursuant to the requirements of the Arizona Department of Public Safety.

Date this 28TH day of DECEMBER 2014

[Signature] 474
Quality Assurance Specialist

SQAP

EXHIBIT G-4
THIS REPORT PREPARED PURSUANT TO DUTY IMPOSED BY A.A.C. R13-10-104 (A)

ARIZONA DEPARTMENT OF PUBLIC SAFETY
STANDARD QUALITY ASSURANCE PROCEDURES
INTOXILYZER MODEL 8000

STANDARD QUALITY ASSURANCE PROCEDURE

QA SPECIALIST KEVIN HUDSPETH #474 AGENCY SEDONA P.D.
DATE 12-28-2014 TIME 0913
INTOXILYZER SERIAL # 80-000924 LOCATION SEDONA P.D.

(KH) 1. Display Reads "PUSH BUTTON TO START"

DIAGNOSTIC TESTS


- (KH) 1. Clock time check.
- (KH) 2. Date check.

OPERATIONAL TESTS

- (KH) 1. Alcohol-free subject test result 0. 000 AC.
- (KH) 2. Error recognition logic system functioning.
Not a Successfully Completed Test Sequence printed
- (KH) 3. Proper sample recognition system.
Not a Successfully Completed Test Sequence printed
Deficient sample printed.
- (KH) 4. Standard Calibration Check standard 0. 100 AC. Result 0. 097 AC.


Instrument is operating properly and accurately. YES NO

COMMENTS

SIGNATURE  474

DPS Form Exh G-4 (Rev 05-01)

I hereby certify that the above and foregoing is a true and correct copy of the record of periodic maintenance for the intoxilyzer 8000 referenced therein maintained by the Sedona Police Department pursuant to the requirements of the Arizona Department of Public Safety.

Date this 28th day of DECEMBER 20 14
 474
Quality Assurance Specialist

SCAP

INTOXILYZER 8000
Location: SEDONA
Serial Number: 80-000924
Core Version: 8105.48
12/28/2014

09:19:51

Standard Lot#: 01513100A1
Last Changed By: K. HUDSPETH #474

QAS: K. HUDSPETH #474
SEDONA PD

Operator: K. HUDSPETH #474
SEDONA PD

Subject: SQAP, TEST
DOB: 04/13/1967
Sex: M Weight: 250

15 Minute Deprivation Period? Yes

Test	g/210L	Time
Air Blank	IPS*	09:20:42
Air Blank	0.000	09:21:10

*Improper Sample

[Signature] 474

(PASSED)

Blow AT WREWG-TIME ✓

I hereby certify that the above and foregoing is a true and correct copy of the record of periodic maintenance for the intoxilyzer 8000 referenced therein maintained by the Sedona Police Department pursuant to the requirements of the Arizona Department of Public Safety.

Date this 28TH day of DECEMBER 20 14

[Signature] 474
Quality Assurance Specialist

Not a Successfully Completed Test Sequence

SQAP

INTOXILYZER 8000
Location: SEDONA
Serial Number: 80-000924
Core Version: 8105.48
12/28/2014

09:22:58

Standard Lot#: 01513100A1
Last Changed By: K. HUDSPETH #474

QAS: K. HUDSPETH #474
SEDONA PD

Operator: K. HUDSPETH #474
SEDONA PD

Subject: SQAP, TEST
DOB: 04/13/1967
Sex: M Weight: 250

15 Minute Deprivation Period? Yes

Test	g/210L	Time
Air Blank	RFI*	09:23:45
Air Blank	0.000	09:24:14

*RFI Detect

K. Hudspeth 474

(PASSED)

RFI DETECT ✓

I hereby certify that the above and foregoing is a true and correct copy of the record of periodic maintenance for the Intoxilyzer 8000 referenced therein maintained by the Sedona Police Department pursuant to the requirements of the Arizona Department of Public Safety.

Date this 28th day of DECEMBER 2014

[Signature] 474
Quality Assurance Specialist

Not a Successfully Completed Test Sequence

SQAP

INTOXILYZER 8000
Location: SEDONA
Serial Number: 80-000924
Core Version: 8105.48
12/28/2014

09:36:31

Standard Lot#: 01513100A1
Last Changed By: K. HUDSPETH #474

QAS: K. HUDSPETH #474
SEDONA PD

Operator: K. HUDSPETH #474
SEDONA PD

Subject: SQAP,TEST
DOB: 04/13/1967
Sex: M Weight: 250

15 Minute Deprivation Period? Yes

Test	g/210L	Time
Air Blank	0.000	09:37:46
Diagnostic Test	Pass	09:38:21
Air Blank	0.000	09:38:48
0.100 Cal Check	0.097	09:39:09
Air Blank	0.000	09:39:39
Subject Test	0.078*	09:42:46
Air Blank	0.000	09:43:18
Air Blank	0.000	09:43:46
Subject Test	XXX**	09:44:14
Air Blank	0.000	09:44:48

*Deficient Sample

**Invalid Sample

Begin new deprivation period

I hereby certify that the above and foregoing is a true and correct copy of the record of periodic maintenance for the intoxilyzer 8000 referenced therein maintained by the Sedona Police Department pursuant to the requirements of the Arizona Department of Public Safety.

[Signature] 474

DEFICIENT SAMPLE ✓
INVALID SAMPLE ✓

Date this 28TH day of DECEMBER 2014
[Signature] 474

(PASSED)

Quality Assurance Specialist

Not a Successfully Completed Test Sequence

AFTER SQAP

EXHIBIT G-2

THIS REPORT PREPARED PURSUANT TO DUTY IMPOSED BY A.A.C. R-13-10-104 (A)

ARIZONA DEPARTMENT OF PUBLIC SAFETY

STANDARD QUALITY ASSURANCE PROCEDURES
INTOXILYZER MODEL 8000

STANDARD CALIBRATION CHECK PROCEDURE

QA SPECIALIST KEVIN HUDSPETH #474 AGENCY SEDONA P.D.

DATE 12-28-2014 TIME 0947

INTOXILYZER SERIAL # 80-000924 LOCATION SEDONA P.D.

(Kt) 1. Ensure that gas tank is attached to instrument and contains a standard alcohol concentration solution . 100 AC.

Pour a standard alcohol concentration solution . ^{OR} N/A AC, into a clean dry simulator and assemble the simulator. Ensure that a tight seal has been made. Turn on the simulator and allow temperature to reach 34°C ± 0.2°C

(Kt) 2. Intoxilyzer 8000 display reads "PUSH BUTTON TO START"

(Kt) 3. Go to the "Control Testing Menu". Select "D" for dry control test or "W" for wet control test. After selection is made press ENTER.


(Kt) 4. Air blank complete.

(Kt) 5. Calibration check complete. Test results 0. 097 AC.

(Kt) 6. Air blank complete.

(Kt) 7. Remove printed record. Attach the record to the complete checklist.

SIGNATURE

 474

DPS Form Exb G-2 (Rev 05-01)

I hereby certify that the above and foregoing is a true and correct copy of the record of periodic maintenance for the intoxilyzer 8000 referenced therein maintained by the Sedona Police Department pursuant to the requirements of the Arizona Department of Public Safety.

Date this 28TH day of DECEMBER 20 14

 474
Quality Assurance Specialist

AFTER SQAP

INTOXILYZER 8000
Location: SEDONA
Serial Number: 80-000924
Core Version: 8105.48
12/28/2014

09:47:11

Standard Lot#: 01513100A1
Last Changed By: K. HUDSPETH #474

QAS: K. HUDSPETH #474
SEDONA PD

Operator: K. HUDSPETH #474
SEDONA PD

Test	g/210L	Time
Air Blank	0.000	09:47:43
0.100 Cal Check	0.097	09:48:04
Air Blank	0.000	09:48:33

K. HUDSPETH 474

Preventative Maintenance
Performed -- Timer Reset

I hereby certify that the above and foregoing is a true and correct copy of the record of periodic maintenance for the Intoxilyzer 8000 referenced therein maintained by the Sedona Police Department pursuant to the requirements of the Arizona Department of Public Safety.

Date this 28TH day of DECEMBER 20 14

[Signature] 474
Quality Assurance Specialist