

BEFORE SQAP

EXHIBIT G-2  
THIS REPORT PREPARED PURSUANT TO DUTY IMPOSED BY A.A.C. R-13-10-104 (A)

ARIZONA DEPARTMENT OF PUBLIC SAFETY  
STANDARD QUALITY ASSURANCE PROCEDURES  
INTOXILYZER MODEL 8000

STANDARD CALIBRATION CHECK PROCEDURE

QA SPECIALIST W. Lopez 30 AGENCY SEDONA P.D.  
DATE 03-06-15 TIME 1810  
INTOXILYZER SERIAL # 80-000924 LOCATION S.P.D.

- 1. Ensure that gas tank is attached to instrument and contains a standard alcohol concentration solution . 100 AC.  
 Pour a standard alcohol concentration solution . N/A AC, into a clean dry simulator and assemble the simulator. Ensure that a tight seal has been made. Turn on the simulator and allow temperature to reach  $34^{\circ}\text{C} \pm 0.2^{\circ}\text{C}$
- (✓) Intoxilyzer 8000 display reads "PUSH BUTTON TO START"
- (✓) Go to the "Control Testing Menu". Select "D" for dry control test or "W" for wet control test. After selection is made press ENTER.
- (✓) Air blank complete.
- (✓) Calibration check complete. Test results 0.099 AC.
- (✓) Air blank complete.
- (✓) Remove printed record. Attach the record to the complete checklist.

SIGNATURE [Signature] 30

DPS Form Exb G-2 (Rev 05-01)  
CHANGE NEW TANK # 14014100A1 (CMI) INC.  
Exp. 7-5-16 CYL # 35  
OLD TANK # 01513100A1 (CMI) INC.  
Exp. 3-1-15 CYL # 015

I hereby certify that the above and forgoing is a true and correct copy of the record of periodic maintenance for the intoxilyzer 8000 referenced therein maintained by the Sedona Police Department pursuant to the requirements of the Arizona Department of Public Safety.

Dated this 06 day of MARCH, 2015.  
[Signature] 30  
Quality Assurance Specialist

INTOXILYZER 8000  
Location: SEDONA  
Serial Number: 80-000924  
Core Version: 8105.48  
03/06/2015 18:17:49

-----  
Standard Lot#: 14014100A1  
Last Changed By: W. LOPEZ #30  
-----

QAS: W. LOPEZ #30 *W/L 30*  
SEDONA PD  
-----

Operator: W. LOPEZ #30  
SEDONA PD  
-----

Test	g/210L	Time
Air Blank	0.000	18:18:18
0.100 Cal Check	0.099	18:18:39
Air Blank	0.000	18:19:09

Preventative Maintenance  
Performed -- Timer Reset

*BEFORE SQAP*

EXHIBIT G-4  
THIS REPORT PREPARED PURSUANT TO DUTY IMPOSED BY A.A.C. R13-10-104 (A)

ARIZONA DEPARTMENT OF PUBLIC SAFETY  
STANDARD QUALITY ASSURANCE PROCEDURES  
INTOXILYZER MODEL 8000

STANDARD QUALITY ASSURANCE PROCEDURE

QA SPECIALIST W. Lopez 30 AGENCY SEDONA P.O.  
DATE 03-06-15 TIME 1824<sup>HR</sup>  
INTOXILYZER SERIAL # 80-000924 LOCATION SEDONA P.O.

(  ) 1. Display Reads "PUSH BUTTON TO START"

DIAGNOSTIC TESTS

(  ) 1. Clock time check.  
(  ) 2. Date check.

OPERATIONAL TESTS

(  ) 1. Alcohol-free subject test result 0. 000 AC.  
(  ) 2. Error recognition logic system functioning.  
Not a Successfully Completed Test Sequence printed  
(  ) 3. Proper sample recognition system.  
Not a Successfully Completed Test Sequence printed  
Deficient sample printed.  
(  ) 4. Standard Calibration Check standard 0. 100 AC. Result 0. 098 AC.

Instrument is operating properly and accurately. YES  NO

COMMENTS CHANGE OUT TANK # 01513100A1 TO TANK # 14014100A1

SIGNATURE W. Lopez 30

INTOXILYZER 8000  
 Location: SEDONA  
 Serial Number: 80-000924  
 Core Version: 8105.48  
 03/06/2015 18:22:12

Standard Lot#: 14014100A1  
 Last Changed By: W. LOPEZ #30

QAS: W. LOPEZ *W* #30  
 SEDONA PD

Operator: W. LOPEZ #30  
 SEDONA PD

Subject: SQAP, TEST  
 DOB: 06/06/1965  
 Sex: M Weight: 210

15 Minute Deprivation Period? Yes

Test	g/210L	Time
Air Blank	0.000	18:23:44
Diagnostic Test	Pass	18:24:19
Air Blank	0.000	18:24:46
0.100 Cal Check	0.099	18:25:07
Air Blank	0.000	18:25:37
Subject Test	0.000	18:26:05
Air Blank	ABT*	18:26:17
Air Blank	0.000	18:26:45

\*Sequence Aborted

*ALCOHOL FREE + PUSH START TEST BUTTON*

**Not a Successfully Completed Test Sequence**



INTOXILYZER 8000  
 Location: SEDONA  
 Serial Number: 80-000924  
 Core Version: 8105.48  
 03/06/2015 18:28:27

Standard Lot#: 14014100A1  
 Last Changed By: W. LOPEZ #30

QAS: W. LOPEZ *WLP* #30  
 SEDONA PD

Operator: W. LOPEZ #30  
 SEDONA PD

Subject: SQAP, TEST  
 DOB: 06/06/1965  
 Sex: M Weight: 210

15 Minute Deprivation Period? Yes

Test	g/210L	Time
Air Blank	IPS*	18:29:31
Air Blank	0.000	18:30:00

\*Improper Sample

*BLOW INTO BREATH TUBE  
 DURING AIRBLANK.*

**Not a Successfully Completed Test Sequence**

INTOXILYZER 8000  
 Location: SEDONA  
 Serial Number: 80-000924  
 Core Version: 8105.48  
 03/06/2015 18:31:34

Standard Lot#: 14014100A1  
 Last Changed By: W. LOPEZ #30

QAS: W. LOPEZ #30  
 SEDONA PD

Operator: W. LOPEZ #30  
 SEDONA PD

Subject: SQAP,TEST  
 DOB: 06/06/1965  
 Sex: M Weight: 210

15 Minute Deprivation Period? Yes

Test	g/210L	Time
Air Blank	RFI*	18:32:36
Air Blank	0.000	18:33:04

\*RFI Detect  
*Key Read by face plate*

**Not a Successfully Completed Test Sequence**

INTOXILYZER 8000  
Location: SEDONA  
Serial Number: 80-000924  
Core Version: 8105.48  
03/06/2015 18:34:26

Standard Lot#: 14014100A1  
Last Changed By: W. LOPEZ #30

QAS: W. LOPEZ #30  
SEDONA PD

Operator: W. LOPEZ #30  
SEDONA PD

Subject: SQAP,TEST  
DOB: 06/06/1965  
Sex: M Weight: 210

15 Minute Deprivation Period? Yes

Test	g/210L	Time
Air Blank	0.000	18:35:58
Diagnostic Test	Pass	18:36:33
Air Blank	0.000	18:37:00
0.100 Cal Check	0.098	18:37:21
Air Blank	0.000	18:37:51
Subject Test	0.058*	18:40:58
Air Blank	0.000	18:41:30
Air Blank	0.000	18:41:57
Subject Test	XXX**	18:42:13
Air Blank	0.000	18:42:47

\*Deficient Sample

\*\*Invalid Sample

Begin new deprivation period

*O2O NOT Blow Long ENOUGH - HAD MOUTH ALCOHOL PRESENT*

**Not a Successfully Completed Test Sequence**

AFTER SQAP

**EXHIBIT G-2**  
**THIS REPORT PREPARED PURSUANT TO DUTY IMPOSED BY A.A.C. R-13-10-104 (A)**

**ARIZONA DEPARTMENT OF PUBLIC SAFETY**

**STANDARD QUALITY ASSURANCE PROCEDURES**  
**INTOXILYZER MODEL 8000**

**STANDARD CALIBRATION CHECK PROCEDURE**

QA SPECIALIST W. Lopez 30 AGENCY SEONA P.O.

DATE 03-06-15 TIME 1844

INTOXILYZER SERIAL # 80-000924 LOCATION SEONA P.O.

( bl ) 1. Ensure that gas tank is attached to instrument and contains a standard alcohol concentration solution . 100 AC.

Pour a standard alcohol concentration solution . <sup>OR</sup> N/A AC, into a clean dry simulator and assemble the simulator. Ensure that a tight seal has been made. Turn on the simulator and allow temperature to reach 34°C ± 0.2°C

( bl ) Intoxilyzer 8000 display reads "PUSH BUTTON TO START"


( bl ) 3. Go to the "Control Testing Menu". Select "D" for dry control test or "W" for wet control test. After selection is made press ENTER.

( bl ) 4. Air blank complete.

( bl ) Calibration check complete. Test results 0. 098 AC.

( bl ) Air blank complete.

( bl ) Remove printed record. Attach the record to the complete checklist.

SIGNATURE  30

NOTE SOAP

INTOXILYZER 8000  
Location: SEDONA  
Serial Number: 80-000924  
Core Version: 8105.48  
03/06/2015 18:44:13

-----  
Standard Lot#: 14014100A1  
Last Changed By: W. LOPEZ #30

-----  
QAS: W. LOPEZ #30 *WLP*  
SEDONA PD

-----  
Operator: W. LOPEZ #30  
SEDONA PD

-----  
Test g/210L Time  
-----  
Air Blank 0.000 18:44:44  
0.100 Cal Check 0.098 18:45:05  
Air Blank 0.000 18:45:35

Preventative Maintenance  
Performed -- Timer Reset