

BEFORE  
SQAP

1 of 2

EXHIBIT G-2

THIS REPORT PREPARED PURSUANT TO DUTY IMPOSED BY A.A.C. R-13-10-104 (A)

ARIZONA DEPARTMENT OF PUBLIC SAFETY

STANDARD QUALITY ASSURANCE PROCEDURES  
INTOXILYZER MODEL 8000

STANDARD CALIBRATION CHECK PROCEDURE

QA SPECIALIST W. Lopez 30 AGENCY SEDONA P.D.

DATE 08-08-15 TIME 1836

INTOXILYZER SERIAL # 80-000924 LOCATION SEDONA P.D.

(). Ensure that gas tank is attached to instrument and contains a standard alcohol concentration solution .100 AC.

Pour a standard alcohol concentration solution N/A <sup>OR</sup> AC, into a clean dry simulator and assemble the simulator. Ensure that a tight seal has been made. Turn on the simulator and allow temperature to reach  $34\text{ }^{\circ}\text{C} \pm 0.2\text{ }^{\circ}\text{C}$

- (). Intoxilyzer 8000 display reads "PUSH BUTTON TO START"
- (). Go to the "Control Testing Menu". Select "D" for dry control test or "W" for wet control test. After selection is made press ENTER.
- (). Air blank complete.
- (). Calibration check complete. Test results 0.098 AC.
- (). Air blank complete.
- (). Remove printed record. Attach the record to the complete checklist.

SIGNATURE [Signature] 30

DPS Form Exh G-2 (Rev 05-01)

PASS

I hereby certify that the above and forgoing is a true and correct copy of the record of periodic maintenance for the intoxilyzer 8000 referenced therein maintained by the Sedona Police Department pursuant to the requirements of the Arizona Department of Public Safety.

Dated this 08 day of AUG., 2015.

[Signature] 30  
Quality Assurance Specialist

BEFORE  
SQAP

2 of 2

INTOXILYZER 8000  
Location: SEDONA  
Serial Number: 80-000924  
Core Version: 8105.48  
08/08/2015 18:36:30

-----  
Standard Lot#: 14014100A1  
Last Changed By: W. LOPEZ #30  
-----

QAS: W. LOPEZ #30 *WLP*  
SEDONA PD  
-----

Operator: W. LOPEZ #30  
SEDONA PD  
-----

Test	g/210L	Time
Air Blank	0.000	18:37:14
0.100 Cal Check	0.098	18:37:35
Air Blank	0.000	18:38:04

Preventative Maintenance  
Performed -- Timer Reset

EXHIBIT G-4

THIS REPORT PREPARED PURSUANT TO DUTY IMPOSED BY A.A.C. R-13-10-104 (A)

ARIZONA DEPARTMENT OF PUBLIC SAFETY

STANDARD QUALITY ASSURANCE PROCEDURES  
INTOXILYZER MODEL 8000

STANDARD QUALITY ASSURANCE PROCEDURE

QA SPECIALIST W. Lopez 30 AGENCY SEDONA P.O.

DATE 08-08-15 TIME 1844

INTOXILYZER SERIAL # 80-000924 LOCATION SEDONA P.O.

1. Display Reads "PUSH BUTTON TO START"

DIAGNOSTIC TESTS

- 1. Clock time check.
- 2. Date check.

OPERATIONAL TESTS

- 1. Alcohol-free subject test result 0. 000 AC.
- 2. Error recognition logic system functioning.  
Not a Successfully Completed Test Sequence printed.
- 3. Proper sample recognition system.  
Not a Successfully Completed Test Sequence printed.  
Deficient sample printed.
- 4. Standard Calibration Check standard 0. 100 AC. Result 0. 097 AC.

Instrument is operating properly and accurately. YES  NO

COMMENTS

ALL TESTS PASSED.

SIGNATURE [Signature] 30

DPS Form Exh G-4 (Rev 05-01)

I hereby certify that the above and forgoing is a true and correct copy of the record of periodic maintenance for the intoxilyzer 8000 referenced therein maintained by the Sedona Police Department pursuant to the requirements of the Arizona Department of Public Safety.

Dated this 08 day of AUG, 2015.

[Signature] 30  
Quality Assurance Specialist

INTOXILYZER 8000  
 Location: SEDONA  
 Serial Number: 80-000924  
 Core Version: 8105.48  
 08/08/2015 18:44:07

Standard Lot#: 14014100A1  
 Last Changed By: W. LOPEZ #30

QAS: W. LOPEZ #30  
 SEDONA PD

Operator: W. LOPEZ #30  
 SEDONA PD

Subject: SQAP, TEST  
 DOB: 05/05/1965  
 Sex: M Weight: 250

15 Minute Deprivation Period? Yes

Test	g/210L	Time
Air Blank	0.000	18:45:29
Diagnostic Test	Pass	18:46:04
Air Blank	0.000	18:46:31
0.100 Cal Check	0.098	18:46:52
Air Blank	0.000	18:47:22
Subject Test	0.000	18:47:41
Air Blank	0.000	18:48:08
Wait	ABT*	
Air Blank	0.000	18:48:46

\*Sequence Aborted

1. ALCOHOL FREE

2. START TEST BUTTON PUSH DURING WAIT PERIOD.

**Not a Successfully Completed Test Sequence**

INTOXILYZER 8000  
Location: SEDONA  
Serial Number: 80-000924  
Core Version: 8105.48  
08/08/2015 18:51:08

Standard Lot#: 14014100A1  
Last Changed By: W. LOPEZ #30

QAS: W. LOPEZ #30  
SEDONA PD

Operator: W. LOPEZ #30  
SEDONA PD

Subject: SQAP, TEST  
DOB: 05/05/1965  
Sex: M Weight: 250

15 Minute Deprivation Period? Yes

Test	g/210L	Time
Air Blank	IPS*	18:52:25
Air Blank	0.000	18:52:53

\*Improper Sample

1. BLOW INTO BREATH TUBE DURING FIRST AIR BLANK.

**Not a Successfully Completed Test Sequence**

INTOXILYZER 8000  
Location: SEDONA  
Serial Number: 80-000924  
Core Version: 8105.48  
08/08/2015 18:55:01

---

Standard Lot#: 14014100A1  
Last Changed By: W. LOPEZ #30

---

QAS: W. LOPEZ #30  
SEDONA PD

---

Operator: W. LOPEZ #30  
SEDONA PD

---

Subject: SQAP,TEST  
DOB: 05/05/1965  
Sex: M Weight: 250

---

15 Minute Deprivation Period? Yes

---

Test	g/210L	Time
Air Blank	RFI*	18:55:46
Air Blank	0.000	18:56:14

\*RFI Detect

1. RADIO KEYED DURING FIRST AIR BLANK

**Not a Successfully Completed Test Sequence**

INTOXILYZER 8000  
 Location: SEDONA  
 Serial Number: 80-000924  
 Core Version: 8105.48  
 08/08/2015 18:57:55

Standard Lot#: 14014100A1  
 Last Changed By: W. LOPEZ #30

QAS: W. LOPEZ *W* #30  
 SEDONA PD

Operator: W. LOPEZ #30  
 SEDONA PD

Subject: SQAP, TEST  
 DOB: 05/05/1965  
 Sex: M Weight: 250

15 Minute Deprivation Period? Yes

Test	g/210L	Time
Air Blank	0.000	18:58:57
Diagnostic Test	Pass	18:59:32
Air Blank	0.000	18:59:59
0.100 Cal Check	0.097	19:00:20
Air Blank	0.000	19:00:50
Subject Test	0.177*	19:03:57
Air Blank	0.000	19:04:30
Air Blank	0.000	19:04:57
Subject Test	XXX**	19:05:34
Air Blank	0.000	19:06:08

\*Deficient Sample

\*\*Invalid Sample

Begin new deprivation period

1. STOPPED BLOWER AFTER 1 SECOND.
2. BLEW WITH SAMPLE IN MOUTH CAUSING IT TO GO UP THEN DRIP AT A RATE OF .006%.

**Not a Successfully Completed Test Sequence**

AGTER-  
SCAP 1062

**EXHIBIT G-2**  
**THIS REPORT PREPARED PURSUANT TO DUTY IMPOSED BY A.A.C. R-13-10-104 (A)**

**ARIZONA DEPARTMENT OF PUBLIC SAFETY**

**STANDARD QUALITY ASSURANCE PROCEDURES**  
**INTOXILYZER MODEL 8000**

**STANDARD CALIBRATION CHECK PROCEDURE**

QA SPECIALIST W. Lopez AGENCY Sedona P.D.

DATE 08-08-15 TIME 1909

INTOXILYZER SERIAL # 80-000924 LOCATION Sedona P.O.

(  ) 1. Ensure that gas tank is attached to instrument and contains a standard alcohol concentration solution . 100 AC.

Pour a standard alcohol concentration solution . <sup>OR</sup> N/A AC, into a clean dry simulator and assemble the simulator. Ensure that a tight seal has been made. Turn on the simulator and allow temperature to reach 34 °C± 0.2 °C

- (  ) 2. Intoxilyzer 8000 display reads "PUSH BUTTON TO START"
- (  ) 3. Go to the "Control Testing Menu". Select "D" for dry control test or "W" for wet control test. After selection is made press ENTER.
- (  ) 4. Air blank complete.
- (  ) 5. Calibration check complete. Test results 0. 097 AC.
- (  ) 6. Air blank complete.
- (  ) 7. Remove printed record. Attach the record to the complete checklist.

SIGNATURE [Signature]

DPS Form Exh G-2 (Rev 05-01) PASS

I hereby certify that the above and forgoing is a true and correct copy of the record of periodic maintenance for the intoxilyzer 8000 referenced therein maintained by the Sedona Police Department pursuant to the requirements of the Arizona Department of Public Safety.

Dated this 08 day of Aug, 20 15. [Signature]

AFTER  
SQAP 2 of 2

INTOXILYZER 8000

Location: SEDONA

Serial Number: 80-000924

Core Version: 8105.48

08/08/2015

19:09:34

-----  
Standard Lot#: 14014100A1

Last Changed By: W. LOPEZ #30

-----  
QAS: W. LOPEZ #30 *WLP*

SEDONA PD

-----  
Operator: W. LOPEZ #30

SEDONA PD

-----

Test	g/210L	Time
Air Blank	0.000	19:10:05
0.100 Cal Check	0.097	19:10:26
Air Blank	0.000	19:10:55

-----

Preventative Maintenance  
Performed -- Timer Reset