

EXHIBIT G-2
THIS REPORT PREPARED PURSUANT TO DUTY IMPOSED BY A.A.C. R-13-10-104 (A)

ARIZONA DEPARTMENT OF PUBLIC SAFETY
STANDARD QUALITY ASSURANCE PROCEDURES
INTOXILYZER MODEL 8000

STANDARD CALIBRATION CHECK PROCEDURE

QA SPECIALIST KEVIN HUDSPETH #474 AGENCY SEDONA POLICE DEPARTMENT

DATE 04/26/2018 TIME 0052

INTOXILYZER SERIAL # 80-000924 LOCATION SEDONA P.D.

(✓) 1. Ensure that gas tank is attached to instrument and contains a standard alcohol concentration solution .100 AC.
OR
Pour a standard alcohol concentration solution N/A AC, into a clean dry simulator and assemble the simulator. Ensure that a tight seal has been made. Turn on the simulator and allow temperature to reach 34°C± 0.2°C

(✓) 2. Intoxilyzer 8000 display reads "PUSH BUTTON TO START"


(✓) 3. Go to the "Control Testing Menu". Select "D" for dry control test or "W" for wet control test. After selection is made press ENTER.

(✓) 4. Air blank complete.

() 5. Calibration check complete. Test results 0. 099 AC.

(✓) 6. Air blank complete.

(✓) 7. Remove printed record. Attach the record to the complete checklist.

SIGNATURE  #474 PASSED

QAS 30-DAY
2 of 2

INTOXILYZER 8000
Location: SEDONA
Serial Number: 80-000924
Core Version: 8105.48
04/26/2018 00:52:38

Standard Lot#: 02617100A3
Last Changed By: K. HUDSPETH #474

QAS: K. HUDSPETH #474
SEDONA PD

Operator: K. HUDSPETH #474
SEDONA PD

Test g/210L Time

Air Blank 0.000 00:53:10
0.100 Cal Check 0.099 00:53:31
Air Blank 0.000 00:54:00

Preventative Maintenance
Performed -- Timer Reset