

EXHIBIT G-2
THIS REPORT PREPARED PURSUANT TO DUTY IMPOSED BY A.A.C. R13-10-104 (A)

ARIZONA DEPARTMENT OF PUBLIC SAFETY
STANDARD QUALITY ASSURANCE PROCEDURES
INTOXILYZER MODEL 8000
STANDARD CALIBRATION CHECK PROCEDURE

QA SPECIALIST RICHARD J. BYRNE III AGENCY AZ DPS
DATE 02/04/08 TIME 1651
INTOXILYZER SERIAL # 80-000929 DPS LOCATION WICKENBURG P.D.

- (98) 1. Ensure that gas tank is attached to instrument and contains a standard alcohol concentration solution 0.100 AC.
OR
Pour a standard alcohol concentration solution N/A AC, into a clean dry simulator and assemble the simulator. Ensure that a tight seal has been made. Turn on the simulator and allow temperature to reach $34^{\circ}\text{C} \pm 0.2^{\circ}\text{C}$
- (99) 2. Intoxilyzer 8000 display reads "PUSH BUTTON TO START"
- (144) 3. Go to the "Control Testing Menu". Select "D" for dry control test or "W" for wet control test. After selection is made press ENTER.
- (298) 4. Air blank completed.
- (643) 5. Calibration check completed. Test results 0.100 AC.
- (248) 6. Air blank completed.
- (448) 7. Remove printed record. Attach the record to the completed checklist.

SIGNATURE Richard J. Byrne III # 6386

DPS Form Exh G-2 (Rev 05-01)

STANDARD CHANGED
LOT # 618601I Exp 07/07/08

INTOXILYZER 8000

Location: WICKENBURG PD

Serial Number: 80-000929

02/04/2008

16:51:54

Standard Lot#: 618601I

Last Changed By: R. BYRNE, III #6386

QAS: R. BYRNE, III #6386

AZDPS

Test	g/210L	Time
Air Blank	0.000	16:52:22
Cal Check	0.100	16:52:40
Air Blank	0.000	16:53:10

Richard A. Byrne III #6386

Operator's Signature