

EXHIBIT G-2
THIS REPORT PREPARED PURSUANT TO DUTY IMPOSED BY A.A.C. R13-10-104 (A)

ARIZONA DEPARTMENT OF PUBLIC SAFETY
STANDARD QUALITY ASSURANCE PROCEDURES
INTOXILYZER MODEL 8000
STANDARD CALIBRATION CHECK PROCEDURE

QA SPECIALIST RICHARD J. BYRNE III AGENCY AZ DPS

DATE 06/20/08 TIME 2231

INTOXILYZER SERIAL # 80-000929 LOCATION WICKENBURG P.D.

- (RB) 1. Ensure that gas tank is attached to instrument and contains a standard alcohol concentration solution 100 AC.
OR
Pour a standard alcohol concentration solution N/A AC, into a clean dry simulator and assemble the simulator. Ensure that a tight seal has been made. Turn on the simulator and allow temperature to reach $34^{\circ}\text{C} \pm 0.2^{\circ}\text{C}$
- (RB) 2. Intoxilyzer 8000 display reads "PUSH BUTTON TO START"
- (RB) 3. Go to the "Control Testing Menu". Select "D" for dry control test or "W" for wet control test. After selection is made press ENTER.
- (RB) 4. Air blank completed.
- (RB) 5. Calibration check completed. Test results 0. 100 AC.
- (RB) 6. Air blank completed.
- (RB) 7. Remove printed record. Attach the record to the completed checklist.

SIGNATURE Richard J. Byrne III 6386

INTOXILYZER 8000

Location: WICKENBURG PD

Serial Number: 80-000929

06/02/2008

22:31:16

Standard Lot#: 618601I

Last Changed By: R. BYRNE, III #6386

QAS: R. BYRNE, III #6386

AZDPS

Test	g/210L	Time
Air Blank	0.000	22:31:44
Cal Check	0.100	22:32:02
Air Blank	0.000	22:32:32

 6386

Operator's Signature

EXHIBIT G-4
THIS REPORT PREPARED PURSUANT TO DUTY IMPOSED BY A.A.C. R13-10-104 (A)

ARIZONA DEPARTMENT OF PUBLIC SAFETY

STANDARD QUALITY ASSURANCE PROCEDURES
INTOXILYZER MODEL 8000

STANDARD QUALITY ASSURANCE PROCEDURE

QA SPECIALIST RICHARD J. BYRNE III AGENCY AZ DPS

DATE 06/02/08 TIME 2233

INTOXILYZER SERIAL # 80-020929 DPS LOCATION WICKENBURG P.D.

(RJB) 1. Display Reads "PUSH BUTTON TO START"

DIAGNOSTIC TESTS

(RJB) 1. Clock time check.

(RJB) 2. Date check.

OPERATIONAL TESTS

(RJB) 1. Alcohol-free subject test result 0.000 AC.

(RJB) 2. Error recognition logic system functioning.
Not a Successfully Completed Test Sequence printed

(RJB) 3. Proper sample recognition system.
Not a Successfully Completed Test Sequence printed
Deficient sample printed.

(RJB) 4. Standard Calibration Check standard 0.100 AC. Result 0.099 AC.

Instrument is operating properly and accurately. YES RJB NO _____

COMMENTS

LOT # 618601I EXP 07/07/08

SIGNATURE Richard J. Byrne III 6386

INTOXILYZER 8000
Location: WICKENBURG PD
Serial Number : 80-000929
06/02/2008

22:33:14

Standard Lot#: 618601I
Last Changed By: R. BYRNE, III #6386

QAS: R. BYRNE, III #6386
AZDPS

Operator: R. BYRNE, III #6386
AZDPS

Subject: SQAP,TEST

DOB: 08/08/1968

Sex: M

Weight: 185

15 Minute Deprivation Period? Yes

Test	g/210L	Time
Air Blank	0.000	22:34:34
Diagnostic Test	Pass	22:35:01
Air Blank	0.000	22:35:29
0.100 Cal Check	0.100	22:35:48
Air Blank	0.000	22:36:17
Subject Test	0.000	22:36:42
Air Blank	ABT*	22:36:49
Air Blank	0.000	22:37:18

*Sequence Aborted

Not a Successfully Completed Test Sequence

INTOXILYZER 8000
Location: WICKENBURG PD
Serial Number : 80-000929
06/02/2008

22:37:50

Standard Lot#: 618601I
Last Changed By: R. BYRNE, III #6386

QAS: R. BYRNE, III #6386
AZDPS

Operator: R. BYRNE, III #6386
AZDPS

Subject: SQAP,TEST

DOB: 08/08/1968

Sex: M Weight: 185

15 Minute Deprivation Period? Yes

Test	g/210L	Time
Air Blank	IPS*	22:38:48
Air Blank	RFI**	22:38:58

*Improper Sample

**RFI Detect

Not a Successfully Completed Test Sequence

INTOXILYZER 8000
Location: WICKENBURG PD
Serial Number : 80-000929
06/02/2008

22:39:22

Standard Lot#: 618601I
Last Changed By: R. BYRNE, III #6386

QAS: R. BYRNE, III #6386
AZDPS

Operator: R. BYRNE III #6386
AZDPS

Subject: SQAP,TEST

DOB: 08/08/1968

Sex: M

Weight: 185

15 Minute Deprivation Period? Yes

Test	g/210L	Time
Air Blank	0.000	22:41:03
Diagnostic Test	Pass	22:41:31
Air Blank	0.000	22:41:58
0.100 Cal Check	0.100	22:42:17
Air Blank	0.000	22:42:46
Subject Test	0.082*	22:45:55
Air Blank	0.000	22:46:26
Air Blank	0.000	22:46:53
Subject Test	XXX**	22:47:25
Air Blank	0.000	22:47:58

*Deficient Sample

**Invalid Sample

Begin new deprivation period

Not a Successfully Completed Test Sequence

EXHIBIT G-2
THIS REPORT PREPARED PURSUANT TO DUTY IMPOSED BY A.A.C. R13-10-104 (A)

ARIZONA DEPARTMENT OF PUBLIC SAFETY
STANDARD QUALITY ASSURANCE PROCEDURES
INTOXILYZER MODEL 8000
STANDARD CALIBRATION CHECK PROCEDURE

QA SPECIALIST RICHARD J. BYRNE III AGENCY AZ DPS

DATE 06/20/08 TIME 2249

INTOXILYZER SERIAL # 80-000929 DPS LOCATION WICKENBURG P.D.

- (SAB) 1. Ensure that gas tank is attached to instrument and contains a standard alcohol concentration solution .100 AC.
OR
Pour a standard alcohol concentration solution N/A AC, into a clean dry simulator and assemble the simulator. Ensure that a tight seal has been made. Turn on the simulator and allow temperature to reach $34^{\circ}\text{C} \pm 0.2^{\circ}\text{C}$
- (SAB) 2. Intoxilyzer 8000 display reads "PUSH BUTTON TO START"
- (SAB) 3. Go to the "Control Testing Menu". Select "D" for dry control test or "W" for wet control test. After selection is made press ENTER.
- (SAB) 4. Air blank completed.
- (SAB) 5. Calibration check completed. Test results 0.099 AC.
- (SAB) 6. Air blank completed.
- (SAB) 7. Remove printed record. Attach the record to the completed checklist.

SIGNATURE Richard J. Byrne III 6386

INTOXILYZER 8000

Location: WICKENBURG PD

Serial Number: 80-000929

06/02/2008

22:49:09

Standard Lot#: 6186011

Last Changed By: R. BYRNE, III #6386

QAS: R. BYRNE, III #6386

AZDPS

Test	g/210L	Time
Air Blank	0.000	22:49:36
Cal Check	0.099	22:49:55
Air Blank	0.000	22:50:25


Operator's Signature