

EXHIBIT G-2
THIS REPORT PREPARED PURSUANT TO DUTY IMPOSED BY A.A.C. R13-10-104 (A)

ARIZONA DEPARTMENT OF PUBLIC SAFETY

STANDARD QUALITY ASSURANCE PROCEDURES
INTOXILYZER MODEL 8000

STANDARD CALIBRATION CHECK PROCEDURE

QA SPECIALIST RICHARD J. BYRNE III 12386 AGENCY AZ DPS

DATE 07/02/08 TIME 1519

INTOXILYZER SERIAL # 80-000929 DPS LOCATION WICKENBURG P.D.

- (R/B) 1. Ensure that gas tank is attached to instrument and contains a standard alcohol concentration solution .100 AC.
OR
Pour a standard alcohol concentration solution N/A AC, into a clean dry simulator and assemble the simulator. Ensure that a tight seal has been made. Turn on the simulator and allow temperature to reach $34^{\circ}\text{C} \pm 0.2^{\circ}\text{C}$
- (R/B) 2. Intoxilyzer 8000 display reads "PUSH BUTTON TO START"
- (R/B) 3. Go to the "Control Testing Menu". Select "D" for dry control test or "W" for wet control test. After selection is made press ENTER.
- (R/B) 4. Air blank completed.
- (R/B) 5. Calibration check completed. Test results 0.100 AC.
- (R/B) 6. Air blank completed.
- (R/B) 7. Remove printed record. Attach the record to the completed checklist.

SIGNATURE Richard J. Byrne III # 12386

DPS Form Exh G-2 (Rev 05-01)

STANDARD CHANGED 07/02/08
LOT # 735501 G EXP 12/28/09


INTOXILYZER 8000
Location: WICKENBURG PD
Serial Number: 80-000929
07/02/2008

15:19:49

Standard Lot#: 735501G
Last Changed By: R. BYRNE, III #6386

QAS: R. BYRNE, III #6386
AZDPS

Test	g/210L	Time
Air Blank	0.000	15:20:16
Cal Check	0.100	15:20:35
Air Blank	0.000	15:21:04



#6386

Operator's Signature