

EXHIBIT G-2
THIS REPORT PREPARED PURSUANT TO DUTY IMPOSED BY A.A.C. R13-10-104 (A)

ARIZONA DEPARTMENT OF PUBLIC SAFETY
STANDARD QUALITY ASSURANCE PROCEDURES
INTOXILYZER MODEL 8000
STANDARD CALIBRATION CHECK PROCEDURE

QA SPECIALIST RICHARD J. BYRNE III AGENCY AZ DPS

DATE 12/07/08 TIME 2224

INTOXILYZER SERIAL # 80-000929 LOCATION WICKENBURG P.D.

- (RJB) 1. Ensure that gas tank is attached to instrument and contains a standard alcohol concentration solution .100 AC.
OR
Pour a standard alcohol concentration solution N/A AC, into a clean dry simulator and assemble the simulator. Ensure that a tight seal has been made. Turn on the simulator and allow temperature to reach $34^{\circ}\text{C} \pm 0.2^{\circ}\text{C}$
- (RJB) 2. Intoxilyzer 8000 display reads "PUSH BUTTON TO START"
- (RJB) 3. Go to the "Control Testing Menu". Select "D" for dry control test or "W" for wet control test. After selection is made press ENTER.
- (RJB) 4. Air blank completed.
- (RJB) 5. Calibration check completed. Test results 0. 100 AC.
- (RJB) 6. Air blank completed.
- (RJB) 7. Remove printed record. Attach the record to the completed checklist.

SIGNATURE Richard J. Byrne III # 6386

INTOXILYZER 8000

Location: WICKENBURG PD

Serial Number: 80-000929

12/07/2008

22:24:02

Standard Lot#: 735501G

Last Changed By: R. BYRNE, III #6386

QAS: R. BYRNE, III #6386

AZDPS

Test	g/210L	Time
Air Blank	0.000	22:24:30
Cal Check	0.100	22:24:48
Air Blank	0.000	22:25:18

Richard J. Byrne III #6386

Operator's Signature

EXHIBIT G-4
THIS REPORT PREPARED PURSUANT TO DUTY IMPOSED BY A.A.C. R13-10-104 (A)

ARIZONA DEPARTMENT OF PUBLIC SAFETY
STANDARD QUALITY ASSURANCE PROCEDURES
INTOXILYZER MODEL 8000

STANDARD QUALITY ASSURANCE PROCEDURE

QA SPECIALIST RICHARD J. BYRNE III AGENCY AZ DPS
DATE 12/07/08 TIME 2228
INTOXILYZER SERIAL # 80-000929 DPS LOCATION WICKENBURG P.D.

~~(P/B)~~ 1. Display Reads "PUSH BUTTON TO START"

DIAGNOSTIC TESTS

- ~~(P/B)~~ 1. Clock time check.
~~(P/B)~~ 2. Date check.

OPERATIONAL TESTS

- ~~(P/B)~~ 1. Alcohol-free subject test result 0. 000 AC.
~~(P/B)~~ 2. Error recognition logic system functioning.
Not a Successfully Completed Test Sequence printed
~~(P/B)~~ 3. Proper sample recognition system.
Not a Successfully Completed Test Sequence printed
Deficient sample printed.
~~(P/B)~~ 4. Standard Calibration Check standard 0. 100 AC. Result 0. 099 AC.

Instrument is operating properly and accurately. YES _____ NO _____

COMMENTS

LOT # 735501 A, EXP 12/28/09

SIGNATURE Richard J. Byrne III #6386

INTOXILYZER 8000

Location: WICKENBURG PD

Serial Number : 80-000929

12/07/2008

22:28:04

Standard Lot#: 735501G

Last Changed By: R. BYRNE, III #6386

QAS: R. BYRNE, III #6386

AZDPS

Operator: R. BYRNE, III #6386

AZDPS

Subject: SQAP,TEST

DOB: 08/08/1978

Sex: M

Weight: 185

15 Minute Deprivation Period? Yes

Test	g/210L	Time
Air Blank	0.000	22:29:27
Diagnostic Test	Pass	22:29:55
Air Blank	0.000	22:30:22
0.100 Cal Check	0.100	22:30:41
Air Blank	0.000	22:31:10
Subject Test	0.000	22:31:33
Air Blank	ABT*	22:31:42
Air Blank	0.000	22:32:11

*Sequence Aborted

Not a Successfully Completed Test Sequence

INTOXILYZER 8000

Location: WICKENBURG PD

Serial Number : 80-000929

12/07/2008

22:33:00

Standard Lot#: 735501G

Last Changed By: R. BYRNE, III #6386

QAS: R. BYRNE, III #6386

AZDPS

Operator: R. BYRNE, III #6386

AZDPS

Subject: SQAP,TEST

DOB: 08/08/1978

Sex: M

Weight: 185

15 Minute Deprivation Period? Yes

Test	g/210L	Time
Air Blank	IPS*	22:33:52
Air Blank	RFI**	22:34:00

*Improper Sample

**RFI Detect

Not a Successfully Completed Test Sequence

INTOXILYZER 8000
Location: WICKENBURG PD
Serial Number : 80-000929
12/07/2008

22:34:29

Standard Lot#: 735501G
Last Changed By: R. BYRNE, III #6386

QAS: R. BYRNE, III #6386
AZDPS

Operator: R. BYRNE, III #6386
AZDPS

Subject: SQAP,TEST

DOB: 08/08/1978

Sex: M

Weight: 185

15 Minute Deprivation Period? Yes

Test	g/210L	Time
Air Blank	0.000	22:35:29
Diagnostic Test	Pass	22:35:56
Air Blank	0.000	22:36:24
0.100 Cal Check	0.100	22:36:43
Air Blank	0.000	22:37:12
Subject Test	0.253*	22:40:20
Air Blank	0.000	22:40:53
Air Blank	0.000	22:41:20
Subject Test	XXX**	22:42:29
Air Blank	0.000	22:43:03

*Deficient Sample

**Invalid Sample

Begin new deprivation period

Not a Successfully Completed Test Sequence

EXHIBIT G-2
THIS REPORT PREPARED PURSUANT TO DUTY IMPOSED BY A.A.C. R13-10-104 (A)

ARIZONA DEPARTMENT OF PUBLIC SAFETY
STANDARD QUALITY ASSURANCE PROCEDURES
INTOXILYZER MODEL 8000
STANDARD CALIBRATION CHECK PROCEDURE

QA SPECIALIST RICHARD J. BYRNE III AGENCY AZ DPS

DATE 12/07/08 TIME 2244

INTOXILYZER SERIAL # 80000929 DPS LOCATION WICKENBURG P.D.

- (144) 1. Ensure that gas tank is attached to instrument and contains a standard alcohol concentration solution .100 AC.
OR
Pour a standard alcohol concentration solution N/A AC, into a clean dry simulator and assemble the simulator. Ensure that a tight seal has been made. Turn on the simulator and allow temperature to reach $34^{\circ}\text{C} \pm 0.2^{\circ}\text{C}$
- (144) 2. Intoxilyzer 8000 display reads "PUSH BUTTON TO START"
- (146) 3. Go to the "Control Testing Menu". Select "D" for dry control test or "W" for wet control test. After selection is made press ENTER.
- (144) 4. Air blank completed.
- (145) 5. Calibration check completed. Test results 0.099 AC.
- (145) 6. Air blank completed.
- (146) 7. Remove printed record. Attach the record to the completed checklist.

SIGNATURE Richard J. Byrne III #6396

INTOXILYZER 8000

Location: WICKENBURG PD

Serial Number: 80-000929

12/07/2008

22:44:58

Standard Lot#: 735501G

Last Changed By: R. BYRNE, III #6386

QAS: R. BYRNE, III #6386

AZDPS

Test	g/210L	Time
Air Blank	0.000	22:45:25
Cal Check	0.099	22:45:44
Air Blank	0.000	22:46:13

Richard M. Byrne III #6386

Operator's Signature