

EXHIBIT G-2  
THIS REPORT PREPARED PURSUANT TO DUTY IMPOSED BY A.A.C. R13-10-104 (A)

ARIZONA DEPARTMENT OF PUBLIC SAFETY  
STANDARD QUALITY ASSURANCE PROCEDURES  
INTOXILYZER MODEL 8000  
STANDARD CALIBRATION CHECK PROCEDURE

QA SPECIALIST RICHARD J BYRNE III AGENCY AZ DPS

DATE 07/09/09 TIME 1936

INTOXILYZER SERIAL # 80-006929 DPS LOCATION WICKENBURG P.D.

- (AB) 1. Ensure that gas tank is attached to instrument and contains a standard alcohol concentration solution .100 AC.  
OR  
Pour a standard alcohol concentration solution N/A AC, into a clean dry simulator and assemble the simulator. Ensure that a tight seal has been made. Turn on the simulator and allow temperature to reach  $34^{\circ}\text{C} \pm 0.2^{\circ}\text{C}$
- (AB) 2. Intoxilyzer 8000 display reads "PUSH BUTTON TO START"
- (AB) 3. Go to the "Control Testing Menu". Select "D" for dry control test or "W" for wet control test. After selection is made press ENTER.
- (AB) 4. Air blank completed.
- (AB) 5. Calibration check completed. Test results 0.094 AC.
- (AB) 6. Air blank completed.
- (AB) 7. Remove printed record. Attach the record to the completed checklist.

SIGNATURE Richard J. Byrne III #6386



EXHIBIT G-4  
THIS REPORT PREPARED PURSUANT TO DUTY IMPOSED BY A.A.C. R13-10-104 (A)

ARIZONA DEPARTMENT OF PUBLIC SAFETY  
STANDARD QUALITY ASSURANCE PROCEDURES  
INTOXILYZER MODEL 8000  
STANDARD QUALITY ASSURANCE PROCEDURE

QA SPECIALIST RICHARD J. BYRNE III AGENCY AZ DPS  
DATE 02/09/09 TIME 1937  
INTOXILYZER SERIAL # 80-000929 DPS LOCATION WICKENBURG P.D

~~(PAS)~~ 1. Display Reads "PUSH BUTTON TO START"

DIAGNOSTIC TESTS

- ~~(PAS)~~ 1. Clock time check.
- ~~(PAS)~~ 2. Date check.

OPERATIONAL TESTS

- ~~(PAS)~~ 1. Alcohol-free subject test result 0.000 AC.
- ~~(PAS)~~ 2. Error recognition logic system functioning.  
Not a Successfully Completed Test Sequence printed
- ~~(PAS)~~ 3. Proper sample recognition system.  
Not a Successfully Completed Test Sequence printed  
Deficient sample printed.
- ~~(PAS)~~ 4. Standard Calibration Check standard 0.100 AC. Result 0.099 AC.

Instrument is operating properly and accurately. YES PAS NO \_\_\_\_\_

COMMENTS  
LOT # 7355016 EXP 12/28/09  
\_\_\_\_\_  
\_\_\_\_\_  
\_\_\_\_\_

SIGNATURE Richard J. Byrne III #6386

INTOXILYZER 8000

Location: WICKENBURG PD

Serial Number : 80-000929

07/09/2009

19:37:40

Standard Lot#: 735501G

Last Changed By: R. BYRNE, III #6386

OAS: R. BYRNE, III #6386

AZDPS

Operator: R. BYRNE, III #6386

AZDPS

Subject: SOAP,TEST

DOB: 08/08/1978

Sex: M

Weight: 185

15 Minute Deprivation Period? Yes

Test	g/210L	Time
Air Blank	0.000	19:38:54
Diagnostic Test	Pass	19:39:22
Air Blank	0.000	19:39:49
0.100 Cal Check	0.099	19:40:08
Air Blank	0.000	19:40:37
Subject Test	0.000	19:41:02
Air Blank	0.000	19:41:29
Five Minute Wait	ABT*	
Air Blank	0.000	19:42:09

\*Sequence Aborted

**Not a Successfully Completed Test Sequence**

INTOXILYZER 8000

Location: WICKENBURG PD

Serial Number : 80-000929

07/09/2009

19:42:34

Standard Lot#: 735501G

Last Changed By: R. BYRNE, III #6386

OAS: R. BYRNE, III #6386

A3DPS

Operator: R. BYRNE, III #6386

AZDPS

Subject: SQAP,TEST

DOB: 08/08/1978

Sex: M

Weight: 185

15 Minute Deprivation Period? Yes

Test	g/210L	Time
Air Blank	IPS*	19:43:17
Air Blank	RFI**	19:43:24

\*Improper Sample

\*\*RFI Detect

**Not a Successfully Completed Test Sequence**

INTOXILYZER 8000

Location: WICKENBURG PD

Serial Number : 80-000929

07/09/2009

19:44:47

Standard Lot#: 735501G

Last Changed By: R. BYRNE, III #6386

OAS: R. BYRNE, III #6386

AZDPS

Operator: R. BYRNE, III #6386

AZDPS

Subject: SOAP,TEST

DOB: 08/08/1978

Sex: M

Weight: 185

15 Minute Deprivation Period? Yes

Test	g/210L	Time
Air Blank	0.000	19:45:43
Diagnostic Test	Pass	19:46:16
Air Blank	0.000	19:46:43
0.100 Cal Check	0.100	19:47:02
Air Blank	0.000	19:47:31
Subject Test	0.156*	19:50:40
Air Blank	0.000	19:51:11
Air Blank	0.000	19:51:39
Subject Test	XXX**	19:52:15
Air Blank	0.000	19:52:49

\*Deficient Sample

\*\*Invalid Sample

Begin new deprivation period

Not a Successfully Completed Test Sequence

EXHIBIT G-2  
THIS REPORT PREPARED PURSUANT TO DUTY IMPOSED BY A.A.C. R13-10-104 (A)

ARIZONA DEPARTMENT OF PUBLIC SAFETY  
STANDARD QUALITY ASSURANCE PROCEDURES  
INTOXILYZER MODEL 8000  
STANDARD CALIBRATION CHECK PROCEDURE

QA SPECIALIST RICHARD J. BYRNE AGENCY AZ DPS  
DATE 07/09/09 TIME 1953  
INTOXILYZER SERIAL # 80-000928 LOCATION WICKENBURG P.D.

- (PAB) 1. Ensure that gas tank is attached to instrument and contains a standard alcohol concentration solution 1100 AC.  
OR  
Pour a standard alcohol concentration solution N/A AC, into a clean dry simulator and assemble the simulator. Ensure that a tight seal has been made. Turn on the simulator and allow temperature to reach  $34^{\circ}\text{C} \pm 0.2^{\circ}\text{C}$
- (PAB) 2. Intoxilyzer 8000 display reads "PUSH BUTTON TO START"
- (PAB) 3. Go to the "Control Testing Menu". Select "D" for dry control test or "W" for wet control test. After selection is made press ENTER.
- (PAB) 4. Air blank completed.
- (PAB) 5. Calibration check completed. Test results 0.09% AC.
- (PAB) 6. Air blank completed.
- (PAB) 7. Remove printed record. Attach the record to the completed checklist.

SIGNATURE Richard J. Byrne #6386

INTOXILYZER 8000

Location: WICKENBURG PD

Serial Number: 80-000929

07/09/2009

19:53:40

-----  
Standard Lot#: 7355016

Last Changed By: R. BYRNE, III #6386

-----  
OAS: R. BYRNE, III #6386

AZDPS  
-----

Test	g/210L	Time
Air Blank	0.000	19:54:08
Cal Check	0.099	19:54:27
Air Blank	0.000	19:54:56

*R. Byrne III* #6386

-----  
Operator's Signature