

EXHIBIT G-2  
THIS REPORT PREPARED PURSUANT TO DUTY IMPOSED BY A.A.C. R13-10-104 (A)

ARIZONA DEPARTMENT OF PUBLIC SAFETY  
STANDARD QUALITY ASSURANCE PROCEDURES  
INTOXILYZER MODEL 8000  
STANDARD CALIBRATION CHECK PROCEDURE

QA SPECIALIST RICHARD J. BYRNE III AGENCY AZ DPS  
DATE 12/09/09 TIME 0049  
INTOXILYZER SERIAL # 80-000929 DPS LOCATION WICKENBURG P.D.

- (AS) 1. Ensure that gas tank is attached to instrument and contains a standard alcohol concentration solution \_\_\_\_\_ AC.  
OR  
Pour a standard alcohol concentration solution N/A AC, into a clean dry simulator and assemble the simulator. Ensure that a tight seal has been made. Turn on the simulator and allow temperature to reach  $34^{\circ}\text{C} \pm 0.2^{\circ}\text{C}$
- (AS) 2. Intoxilyzer 8000 display reads "PUSH BUTTON TO START"
- (AS) 3. Go to the "Control Testing Menu". Select "D" for dry control test or "W" for wet control test. After selection is made press ENTER.
- (AS) 4. Air blank completed.
- (AS) 5. Calibration check completed. Test results 0. 100 AC.
- (AS) 6. Air blank completed.
- (AS) 7. Remove printed record. Attach the record to the completed checklist.

SIGNATURE Richard J. Byrne III #386

INTOXILYZER 8000

Location: WICKENBURG PD

Serial Number: 80-000929

12/09/2009

00:49:01

Standard Lot#: 735501G

Last Changed By: R. BYRNE, III #6386

QAS: R. BYRNE, III #6386

AZEPS

Test	g/210L	Time
Air Blank	0.000	00:49:29
Cal Check	0.100	00:49:47
Air Blank	0.000	00:50:17

*Richard Byrne III*

#6386

Operator's Signature