

EXHIBIT G-2  
THIS REPORT PREPARED PURSUANT TO DUTY IMPOSED BY A.A.C. R13-10-104 (A)

ARIZONA DEPARTMENT OF PUBLIC SAFETY  
STANDARD QUALITY ASSURANCE PROCEDURES  
INTOXILYZER MODEL 8000  
STANDARD CALIBRATION CHECK PROCEDURE

QA SPECIALIST RICHARD J. BYRNE III AGENCY AZ DPS  
DATE 12/28/09 TIME 2342  
INTOXILYZER SERIAL # 80-000929 DPS LOCATION WICKEN Bldg

- (49) 1. Ensure that gas tank is attached to instrument and contains a standard alcohol concentration solution 100 AC.  
OR  
Pour a standard alcohol concentration solution NA AC, into a clean dry simulator and assemble the simulator. Ensure that a tight seal has been made. Turn on the simulator and allow temperature to reach  $34^{\circ}\text{C} \pm 0.2^{\circ}\text{C}$
- (95) 2. Intoxilyzer 8000 display reads "PUSH BUTTON TO START"
- (95) 3. Go to the "Control Testing Menu". Select "D" for dry control test or "W" for wet control test. After selection is made press ENTER.
- (95) 4. Air blank completed.
- (95) 5. Calibration check completed. Test results 0.100 AC.
- (95) 6. Air blank completed.
- (95) 7. Remove printed record. Attach the record to the completed checklist.

SIGNATURE Richard J. Byrne III #6384

DPS Form Exh G-2 (Rev 05-01)

STANDARD CHANGE  
LOT # 590854  
EXP 03/05/2011

INTOXILYZER 8000

Location: WICKENBURG PD

Serial Number: 80-000929

12/28/2009

23:42:11

-----  
Standard Lot#: 590854

Last Changed By: R. BYRNE, III #6386

-----  
QAS: R. BYRNE, III #6386

AZDPS

-----  
Test                      g/210L                      Time  
-----  
Air Blank                      0.000                      23:42:38  
Cal Check                      0.100                      23:42:57  
Air Blank                      0.000                      23:43:27

*Richard Byrne III #6386*  
-----  
Operator's Signature