

EXHIBIT G-2  
THIS REPORT PREPARED PURSUANT TO DUTY IMPOSED BY A.A.C. R13-10-104 (A)

ARIZONA DEPARTMENT OF PUBLIC SAFETY  
STANDARD QUALITY ASSURANCE PROCEDURES  
INTOXILYZER MODEL 8000  
STANDARD CALIBRATION CHECK PROCEDURE

QA SPECIALIST RICHARD J. BYRNE III AGENCY AZ DPS

DATE 02/10/10 TIME 1343

INTOXILYZER SERIAL # 80-000929 DPS LOCATION WICKENBURG P. D.

- RB* 1. Ensure that gas tank is attached to instrument and contains a standard alcohol concentration solution .100 AC.  
OR  
Pour a standard alcohol concentration solution N/A AC, into a clean dry simulator and assemble the simulator. Ensure that a tight seal has been made. Turn on the simulator and allow temperature to reach  $34^{\circ}\text{C} \pm 0.2^{\circ}\text{C}$
- RB* 2. Intoxilyzer 8000 display reads "PUSH BUTTON TO START"
- RB* 3. Go to the "Control Testing Menu". Select "D" for dry control test or "W" for wet control test. After selection is made press ENTER.
- RB* 4. Air blank completed.
- RB* 5. Calibration check completed. Test results 0.100 AC.
- RB* 6. Air blank completed.
- RB* 7. Remove printed record. Attach the record to the completed checklist.

SIGNATURE *Richard J. Byrne III* # 6386

INTOXILYZER 8000

Location: WICKENBURG PD

Serial Number: 80-000929

02/10/2010

13:43:13

Standard Lot#: 590954

Last Changed By: R. BYRNE, TJI 46286

QAS: R. BYRNE, TJI 46286

AZDPS

Test	g. 210L	Time
Air Blank	0.000	13:44:03
Cal Check	0.100	13:44:57
Air Blank	0.000	13:44:59



# 6386

Operator's Signature

EXHIBIT G-4  
THIS REPORT PREPARED PURSUANT TO DUTY IMPOSED BY A.A.C. R13-10-104 (A)

ARIZONA DEPARTMENT OF PUBLIC SAFETY  
STANDARD QUALITY ASSURANCE PROCEDURES  
INTOXILYZER MODEL 8000  
STANDARD QUALITY ASSURANCE PROCEDURE

QA SPECIALIST RICHARD J. BYRNE III AGENCY A2 DPS  
DATE 02/10/10 TIME 1348  
INTOXILYZER SERIAL # 80-000929 DPS LOCATION WICKENBURG P. D.

~~(RJB)~~ 1. Display Reads "PUSH BUTTON TO START"

DIAGNOSTIC TESTS

- ~~(RJB)~~ 1. Clock time check.  
~~(RJB)~~ 2. Date check.

OPERATIONAL TESTS

- ~~(RJB)~~ 1. Alcohol-free subject test result 0.000 AC.  
~~(RJB)~~ 2. Error recognition logic system functioning.  
Not a Successfully Completed Test Sequence printed  
~~(RJB)~~ 3. Proper sample recognition system.  
Not a Successfully Completed Test Sequence printed  
Deficient sample printed.  
~~(RJB)~~ 4. Standard Calibration Check standard 0.100 AC. Result 0.100 AC.

Instrument is operating properly and accurately. YES RJB NO \_\_\_\_\_

COMMENTS

LOT # 590854 EXP: 03/05/11

SIGNATURE

Richard J. Byrne III #6386

INTOXILYZER 8000

Location: WICKENBURG PD

Serial Number : 80-000929

02/10/2010

13:48:13

Standard Lot#: 590854

Last Changed By: R. BYRNE, III #6386

QAS: R. BYRNE, III #6386

AZDPS

Operator: R. BYRNE, III #6386

AZDPS

Subject: SOAP, TEST

DOB: 08/08/1968

Sex: M

Weight: 185

15 Minute Deprivation Period? Yes

Test	g/210L	Time
Air Blank	0.000	13:49:10
Diagnostic Test	Pass	13:49:37
Air Blank	0.000	13:50:00
0.100 Cal Check	0.099	13:50:21
Air Blank	0.000	13:50:51
Subject Test	0.000	13:51:19
Air Blank	ABT*	13:51:19
Air Blank	0.000	13:51:46

\*Sequence Aborted

Not a Successfully Completed Test Sequence

INTOXILYZER 8000

Location: WICKENBURG PD

Serial Number : 80-000929

02/10/2010

13:52:29

Standard Lot#: 590854

Last Changed By: R. BYRNE, III #638

OAS: R. BYRNE, III #638

AZDES

Operator: R. BYRNE, III #638

AZDES

Subject: SOAP,TEST

DOB: 08/08/1968

Sex: M

Weight: 185

15 Minute Deprivation Period? Yes

Test	g/210L	Time
Air Blank	IPS*	13:53:08
Air Blank	RFI**	13:53:15

\*Improper Sample

\*\*RFI Detect

Not a Successfully Completed Test Sequence

INTOXILYZER 8000

Location: WICKENBURG PD

Serial Number : 80-000929

02/10/2010

13:53:44

Standard Lot#: 590854

Last Changed By: R. BYRNE, III #6386

QAS: R. BYRNE, III #6386

AEDPS

Operator: R. BYRNE, III #6386

AEDPS

Subject: SOAP,TEST

DOB: 08/08/1968

Sex: M

Weight: 185

15 Minute Deprivation Period? Yes

Test	g/210L	Time
Air Blank	0.000	13:54:43
Diagnostic Test	Pass	13:55:09
Air Blank	0.000	13:55:36
0.100 Cal Check	0.100	13:55:53
Air Blank	0.000	13:56:20
Subject Test	XXX*	13:56:47
Air Blank	0.000	13:57:14

\*Invalid Sample

Begin new deprivation period

Not a Successfully Completed Test Sequence

ENOXILYZER 8000

Location: WICKENBURG PD

Serial Number : 80-000929

02/10/2010

13:57:26

Standard Lot#: 590854

Last Changed By: R. BYRNE, III #635

QAS: R. BYRNE, III #635

AZDPS

Operator: R. BYRNE, III #635

AZDPS

Subject: SOAP,TEST

DOB: 08/08/1968

Sex: M

Weight: 185

15 Minute Deprivation Period? Yes

Test	g/210L	Time
Air Blank	0.000	13:58:20
Diagnostic Test	Pass	13:58:37
Air Blank	0.000	13:59:13
0.100 Cal Check	0.100	13:59:44
Air Blank	0.000	14:00:13
Subject Test	0.141*	14:03:20
Air Blank	0.000	14:03:52
Air Blank	0.000	14:04:21
Subject Test	XXX**	14:04:45
Air Blank	0.000	14:05:18

\*Deficient Sample

\*\*Invalid Sample

Begin new deprivation period

Not a Successfully Completed Test Sequence

EXHIBIT G-2  
THIS REPORT PREPARED PURSUANT TO DUTY IMPOSED BY A.A.C. R13-10-104 (A)

ARIZONA DEPARTMENT OF PUBLIC SAFETY  
STANDARD QUALITY ASSURANCE PROCEDURES  
INTOXILYZER MODEL 8000  
STANDARD CALIBRATION CHECK PROCEDURE

QA SPECIALIST RICHARD J. BYRNE III AGENCY AZ DPS

DATE 02/10/10 TIME 1405

INTOXILYZER SERIAL # 80-000929 DPS LOCATION WICKENBURG

- 1. Ensure that gas tank is attached to instrument and contains a standard alcohol concentration solution .100 AC.  
OR  
Pour a standard alcohol concentration solution N/A AC, into a clean dry simulator and assemble the simulator. Ensure that a tight seal has been made. Turn on the simulator and allow temperature to reach  $34^{\circ}\text{C} \pm 0.2^{\circ}\text{C}$
- 2. Intoxilyzer 8000 display reads "PUSH BUTTON TO START"
- 3. Go to the "Control Testing Menu". Select "D" for dry control test or "W" for wet control test. After selection is made press ENTER.
- 4. Air blank completed.
- 5. Calibration check completed. Test results 0.100 AC.
- 6. Air blank completed.
- 7. Remove printed record. Attach the record to the completed checklist.

SIGNATURE *Richard J. Byrne III* # 6386

ESTABLISHED 1870  
LAWSON & SONS, INC.  
Serial Number 11 000000  
01 107 000

STANDARD LOCKS 500004

Last Changed By: W. BYRNE, 117 2000

000000 000000 000000

0.00

Item	Quantity	Price
Air Blank	0.000	0.000
Air Check	0.100	0.000
Air Blank	0.000	0.000

*Richard Byrne* #6386  
Operator Signature