

EXHIBIT G-2  
THIS REPORT PREPARED PURSUANT TO DUTY IMPOSED BY A.A.C. R13-10-104 (A)

ARIZONA DEPARTMENT OF PUBLIC SAFETY  
STANDARD QUALITY ASSURANCE PROCEDURES  
INTOXILYZER MODEL 8000  
STANDARD CALIBRATION CHECK PROCEDURE

QA SPECIALIST RICHARD J. BYRNE III AGENCY AZ DPS

DATE 07/24/10 TIME 1151

INTOXILYZER SERIAL # 80-000929 DPS LOCATION WICKHAMBURG PD

- 1A** 1. Ensure that gas tank is attached to instrument and contains a standard alcohol concentration solution .100 AC.  
OR  
Pour a standard alcohol concentration solution N/A AC, into a clean dry simulator and assemble the simulator. Ensure that a tight seal has been made. Turn on the simulator and allow temperature to reach  $34^{\circ}\text{C} \pm 0.2^{\circ}\text{C}$
- 1B** 2. Intoxilyzer 8000 display reads "PUSH BUTTON TO START"
- 2B** 3. Go to the "Control Testing Menu". Select "D" for dry control test or "W" for wet control test. After selection is made press ENTER.
- 3B** 4. Air blank completed.
- 4B** 5. Calibration check completed. Test results 0.099 AC.
- 5B** 6. Air blank completed.
- 6B** 7. Remove printed record. Attach the record to the completed checklist.

SIGNATURE Richard J. Byrne III #6386

INTOXILYZER 8000

Location: WICKENBURG PD

Serial Number: 80-000929

Core Version: 8105.48

07/24/2010

11:51:31

-----  
Standard Lot#: 590854

Last Changed By: R. BYRNE, III #6386

-----  
QAS: R. BYRNE, III #6386

AZDPS

-----  
Operator: R. BYRNE, III #6386

AZDPS

-----  
Test                      g/210L                      Time  
-----  
Air Blank                      0.000                      11:52:03  
0.100 Cal Check                0.099                      11:52:24  
Air Blank                      0.000                      11:52:53

Preventative Maintenance  
Performed -- Timer Reset