

EXHIBIT G-2
THIS REPORT PREPARED PURSUANT TO DUTY IMPOSED BY A.A.C. R13-10-104 (A)

ARIZONA DEPARTMENT OF PUBLIC SAFETY
STANDARD QUALITY ASSURANCE PROCEDURES
INTOXILYZER MODEL 8000

STANDARD CALIBRATION CHECK PROCEDURE


QA SPECIALIST RICHARD J. BYRNE III AGENCY AZ DPS

DATE 09/23/10 TIME 0839

INTOXILYZER SERIAL # 80-000929 LOCATION WICKENBURG P.D.

- RJB* 1. Ensure that gas tank is attached to instrument and contains a standard alcohol concentration solution .100 AC.
OR
Pour a standard alcohol concentration solution N/A AC, into a clean dry simulator and assemble the simulator. Ensure that a tight seal has been made. Turn on the simulator and allow temperature to reach $34^{\circ}\text{C} \pm 0.2^{\circ}\text{C}$
- RJB* 2. Intoxilyzer 8000 display reads "PUSH BUTTON TO START"
- RJB* 3. Go to the "Control Testing Menu". Select "D" for dry control test or "W" for wet control test. After selection is made press ENTER.
- RJB* 4. Air blank completed.
- RJB* 5. Calibration check completed. Test results 0.100 AC.
- RJB* 6. Air blank completed.
- RJB* 7. Remove printed record. Attach the record to the completed checklist.

SIGNATURE Richard J. Byrne III #6386



INTOXILYZER 8000
 Location: WICKENBURG PD
 Serial Number: 80-000929
 Core Version: 8105.48
 09/23/2010 08:39:06

 Standard Lot#: 590854
 Last Changed By: R. BYRNE, III #6386

 QAS: R. BYRNE, III #6386
 AZDPS

 Operator: R. BYRNE, III #6386
 AZDPS

Test	g/210L	Time
Air Blank	0.000	08:39:37
0.100 Cal Check	0.100	08:39:58
Air Blank	0.000	08:40:28

Preventative Maintenance
 Performed -- Timer Reset