

EXHIBIT G-2
THIS REPORT PREPARED PURSUANT TO DUTY IMPOSED BY A.A.C. R13-10-104 (A)

ARIZONA DEPARTMENT OF PUBLIC SAFETY
STANDARD QUALITY ASSURANCE PROCEDURES
INTOXILYZER MODEL 8000
STANDARD CALIBRATION CHECK PROCEDURE

QA SPECIALIST RICHARD J. BYRNE III AGENCY AZ DPS
DATE 02/21/11 TIME 0115
INTOXILYZER SERIAL # 80-000529 DPS LOCATION WICKENBURG PD.

- (H.B.) 1. Ensure that gas tank is attached to instrument and contains a standard alcohol concentration solution - 100 AC.
OR
Pour a standard alcohol concentration solution N/A AC, into a clean dry simulator and assemble the simulator. Ensure that a tight seal has been made. Turn on the simulator and allow temperature to reach $34^{\circ}\text{C} \pm 0.2^{\circ}\text{C}$
- (P.S.) 2. Intoxilyzer 8000 display reads "PUSH BUTTON TO START"
- (H.B.) 3. Go to the "Control Testing Menu". Select "D" for dry control test or "W" for wet control test. After selection is made press ENTER.
- (P.S.) 4. Air blank completed.
- (H.B.) 5. Calibration check completed. Test results 0.100 AC.
- (P.S.) 6. Air blank completed.
- (H.B.) 7. Remove printed record. Attach the record to the completed checklist.

SIGNATURE Richard J. Byrne III #6386

INNOVITY25P 8000
Location: WICKEDDPRG PD
Serial Number: 80-000929
Core Version: 3105.43
02/21/2011

01:15:16

Standard Lot#: 590854
Last Changed By: R. BYRNE, III #6386

QAS: P. BYRNE, III #6386
AZDPS

Operator: R. BYRNE, III #6386
AZDPS

| Test | q/210L | Time |
|-----------------|--------|----------|
| Air Blank | 0.000 | 01:15:50 |
| 0.100 Cal Check | 0.100 | 01:16:10 |
| Air Blank | 0.000 | 01:16:40 |

Preventative Maintenance
Performed -- Timer Reset

EXHIBIT G-4
THIS REPORT PREPARED PURSUANT TO DUTY IMPOSED BY A.A.C. R13-10-104 (A)

ARIZONA DEPARTMENT OF PUBLIC SAFETY
STANDARD QUALITY ASSURANCE PROCEDURES
INTOXILYZER MODEL 8000
STANDARD QUALITY ASSURANCE PROCEDURE

QA SPECIALIST RICHARD J. BYRNE III AGENCY AZ DPS
DATE 02/21/11 TIME 0118
INTOXILYZER SERIAL # B000929 DPS LOCATION WICKENBURG P.D.

AB 1. Display Reads "PUSH BUTTON TO START"

DIAGNOSTIC TESTS

- AB 1. Clock time check.
AB 2. Date check.

OPERATIONAL TESTS

- AB 1. Alcohol-free subject test result 0. 000 AC.
AB 2. Error recognition logic system functioning.
Not a Successfully Completed Test Sequence printed
AB 3. Proper sample recognition system.
Not a Successfully Completed Test Sequence printed
Deficient sample printed.
AB 4. Standard Calibration Check standard 0. 100 AC. Result 0. 099 AC.

Instrument is operating properly and accurately. YES AB NO

COMMENTS

LOT # 590854 EXP 03/05/11

SIGNATURE Richard J. Byrne III #6386

INTOXILYZER 8000
Location: WICKENBURG PD
Serial Number: 80-000929
Core Version: 8105.48
02/21/2011

01:18:56

Standard Lot#: 590854
Last Changed By: R. BYRNE, III #6386

QAS: R. BYRNE, III #6386
AZDPS

Operator: R. BYRNE, III #6386
AZDPS

Subject: SQAP,TEST
DOB: 08/08/1968
Sex: M Weight: 185

15 Minute Deprivation Period? Yes

| Test | g/210L | Time |
|-----------------|--------|----------|
| Air Blank | 0.000 | 01:19:55 |
| Diagnostic Test | Pass | 01:20:30 |
| Air Blank | 0.000 | 01:20:57 |
| 0.100 Cal Check | 0.100 | 01:21:18 |
| Air Blank | 0.000 | 01:21:48 |
| Subject Test | 0.000 | 01:22:10 |
| Air Blank | ABT* | 01:22:16 |
| Air Blank | 0.000 | 01:22:44 |

*Sequence Aborted

Not a Successfully Completed Test Sequence

INTOXILYZER 8000

Location: WICKENBURG PD

Serial Number: 80-000929

Core Version: 8105.48

02/21/2011

01:22:49

Standard Lot#: 590854

Last Changed By: R. BYRNE, III #6386

QAS: R. BYRNE, III #6386

AZDPS

Operator: R. BYRNE, III #6386

AZDPS

Subject: SQAP,TEST

DOB: 08/08/1968

Sex: M

Weight: 185

15 Minute Deprivation Period? Yes

| Test | g/210L | Time |
|-----------|--------|----------|
| Air Blank | IPS* | 01:23:30 |
| Air Blank | RFI** | 01:23:37 |

*Improper Sample

**RFI Detect

Not a Successfully Completed Test Sequence

INTOXILYZER 8000

Location: WICKENBURG PD

Serial Number: 80-000929

Core Version: 8105.48

02/21/2011

01:32:46

Standard Lot#: 590854

Last Changed By: R. BYRNE, III #6386

OAS: R. BYRNE, III #6386

AZDPS

Operator: R. BYRNE, III #6386

AZDPS

Subject: SQAP,TEST

DOB: 08/08/1968

Sex: M

Weight: 185

15 Minute Deprivation Period? Yes

| Test | g/210L | Time |
|-----------------|--------|----------|
| Air Blank | 0.000 | 01:33:40 |
| Diagnostic Test | Pass | 01:34:15 |
| Air Blank | 0.000 | 01:34:43 |
| 0.100 Cal Check | 0.100 | 01:35:03 |
| Air Blank | 0.000 | 01:35:33 |
| Subject Test | 0.393* | 01:38:40 |
| Air Blank | 0.000 | 01:39:13 |
| Air Blank | 0.000 | 01:39:41 |
| Subject Test | XXX** | 01:39:50 |
| Air Blank | 0.000 | 01:40:33 |

*Deficient Sample

**Invalid Sample

Begin new deprivation period

Not a Successfully Completed Test Sequence

EXHIBIT G-2
THIS REPORT PREPARED PURSUANT TO DUTY IMPOSED BY A.A.C. R13-10-104 (A)

ARIZONA DEPARTMENT OF PUBLIC SAFETY
STANDARD QUALITY ASSURANCE PROCEDURES
INTOXILYZER MODEL 8000
STANDARD CALIBRATION CHECK PROCEDURE

QA SPECIALIST RICHARD J. BYRNE AGENCY AZ DPS
DATE 02/21/11 TIME 0140
INTOXILYZER SERIAL # 80-000929 DPS LOCATION WICKENBURG P.D.

- (RJB) 1. Ensure that gas tank is attached to instrument and contains a standard alcohol concentration solution 0.100 AC.
OR
Pour a standard alcohol concentration solution N/A AC, into a clean dry simulator and assemble the simulator. Ensure that a tight seal has been made. Turn on the simulator and allow temperature to reach $34^{\circ}\text{C} \pm 0.2^{\circ}\text{C}$
- (RJB) 2. Intoxilyzer 8000 display reads "PUSH BUTTON TO START"
- (RJB) 3. Go to the "Control Testing Menu". Select "D" for dry control test or "W" for wet control test. After selection is made press ENTER.
- (RJB) 4. Air blank completed.
- (RJB) 5. Calibration check completed. Test results 0.095 AC.
- (RJB) 6. Air blank completed.
- (RJB) 7. Remove printed record. Attach the record to the completed checklist.

SIGNATURE Richard J. Byrne III #6386

EPH0X1057ER 8200
Local Job: A113048080 89
Serial Number: 80-000927
Core Version: 6105.49
02/21/2011 01:49:50

Standard Job#: 590574
Last Changed By: W. BYRNE, FI 46286
Date: K. BYRNE, FI 46286
AZDPC
Operator: P. BYRNE, FI 46286
ADDS

| Test | 9/2/06 | Time |
|-----------------|--------|----------|
| Air Blank | 0.600 | 01:11:21 |
| 0.100 Cal Check | 0.699 | 01:11:32 |
| Air Blank | 0.000 | 01:12:10 |

Preventative Maintenance
Performed -- Timer Reset