

EXHIBIT G-2  
THIS REPORT PREPARED PURSUANT TO DUTY IMPOSED BY A.A.C. R13-10-104 (A)

ARIZONA DEPARTMENT OF PUBLIC SAFETY  
STANDARD QUALITY ASSURANCE PROCEDURES  
INTOXILYZER MODEL 8000  
STANDARD CALIBRATION CHECK PROCEDURE

QA SPECIALIST RICHARD J. BYRNE III AGENCY AZ DPS

DATE 07/05/11 TIME 2346

INTOXILYZER SERIAL # 80-000929 DPS LOCATION WICKENBURG P.D.

- RJB* 1. Ensure that gas tank is attached to instrument and contains a standard alcohol concentration solution 100 AC.  
OR  
Pour a standard alcohol concentration solution N/A AC, into a clean dry simulator and assemble the simulator. Ensure that a tight seal has been made. Turn on the simulator and allow temperature to reach  $34^{\circ}\text{C} \pm 0.2^{\circ}\text{C}$
- RJB* 2. Intoxilyzer 8000 display reads "PUSH BUTTON TO START"
- RJB* 3. Go to the "Control Testing Menu". Select "D" for dry control test or "W" for wet control test. After selection is made press ENTER.
- RJB* 4. Air blank completed.
- RJB* 5. Calibration check completed. Test results 0.098 AC.
- RJB* 6. Air blank completed.
- RJB* 7. Remove printed record. Attach the record to the completed checklist.

SIGNATURE Richard J. Byrne III #6386

INTOXILYZER 8000

Location: WICKENBURG PD

Serial Number: 80-000929

Core Version: 8105.48

07/05/2011

23:46:56

-----  
Standard Lot#: 25310100A1

Last Changed By: R. BYRNE, III #6386  
-----

QAS: R. BYRNE, III #6386

AZDPS  
-----

Operator: R. BYRNE, III #6386

AZDPS  
-----

Test	g/210L	Time
Air Blank	0.000	23:47:28
0.100 Cal Check	0.098	23:47:49
Air Blank	0.000	23:48:18

Preventative Maintenance  
Performed -- Timer Reset