

EXHIBIT G-2
THIS REPORT PREPARED PURSUANT TO DUTY IMPOSED BY A.A.C. R13-10-104 (A)

ARIZONA DEPARTMENT OF PUBLIC SAFETY
STANDARD QUALITY ASSURANCE PROCEDURES
INTOXILYZER MODEL 8000
STANDARD CALIBRATION CHECK PROCEDURE

QA SPECIALIST RICHARD S. BYRNE III AGENCY AZ DPS

DATE 12/24/11 TIME 2342

INTOXILYZER SERIAL # 80-000924 D5 LOCATION WICKERBURG P.D.

- (AB) 1. Ensure that gas tank is attached to instrument and contains a standard alcohol concentration solution .100 AC.
OR
Pour a standard alcohol concentration solution N/A AC, into a clean dry simulator and assemble the simulator. Ensure that a tight seal has been made. Turn on the simulator and allow temperature to reach $34^{\circ}\text{C} \pm 0.2^{\circ}\text{C}$
- (AB) 2. Intoxilyzer 8000 display reads "PUSH BUTTON TO START"
- (AB) 3. Go to the "Control Testing Menu". Select "D" for dry control test or "W" for wet control test. After selection is made press ENTER.
- (AB) 4. Air blank completed.
- (AB) 5. Calibration check completed. Test results 0.098 AC.
- (AB) 6. Air blank completed.
- (AB) 7. Remove printed record. Attach the record to the completed checklist.

SIGNATURE Richard Byrne III # 16386

↑↑↑↑↑

INNOVATIZER 8000
Location: WICKENBURG PD;
Serial Number: 80-000929
Core Version: 6105.48
12/04/2011 23:42:15

Standard Job#: 25310100A1
Last Changed By: R. BYRNE, LII #6386

OAS: R. BYRNE LII #6386
AZ DPS

Operator: R. BYRNE LII #6386
AZ DPS

Test	g/210L	Time
Air Blank	0.000	23:42:47
0.100 Cal Check	0.098	23:43:08
Air Blank	0.000	23:43:38

Preventative Maintenance
Performed -- Timer Reset

RECEIVED

DEC 5 2011

WICKENBURG POLICE DEPARTMENT